

# AMSafe® Pre-Filled Normal Saline Flush Syringe

## 0.9% Sodium Chloride Injection, USP

### Features

- Provides all convenient size syringes & volumes
- DEHP-Free, PVC-Free, Latex Free
- Sterile, non-pyrogenic and non-toxic
- Offers sterile fluid path (dust cover) and sterile field configurations
- Reduces rebound
- Eliminates needles and vials
- Eliminates potentials of cross-contamination
- Complies with all PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared



30 syringes/inner box

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

Item Code	Description	Qty/Unit
<b>IVF1203</b>	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF1205</b>	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF1210</b>	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF2020</b>	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case

### Sterile Field Flush Syringes (Sterile Fluid Pathway and Sterile Syringe)

Item Code	Description	Qty/Unit
<b>IVF1203SS</b>	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF1205SS</b>	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF1210SS</b>	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case
<b>IVF2020SS</b>	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case

## AMShield® Disinfecting Cap

The AMShield® Disinfecting Cap for needleless connectors is an alcohol-containing cap that twists onto ports for disinfection and protection.

### Features

- Fits needle-free valves
- Stays securely in place
- Contains 70% Isopropyl Alcohol
- Helps reduce CLABSI (central line-associated bloodstream infections)
- Available individually packaged or as IV pole multi-strips (10caps/strip) to fit an institution's requirements
- Bright orange color is easy to identify when installed on the IV line for visual verification of compliance monitoring
- Sterile, DEHP-free, latex-free, single-use-only



### AMShield® Disinfecting Cap

Item Code	Description	Quantity/Unit	
<b>DCF2000</b>	AMShield® Disinfecting Cap, individually packaged	200 Caps/Box	10 Boxes/Case
<b>DCF2500S</b>	AMShield® Disinfecting Cap, 10 Caps mounted on a strip	25 Strips/Box	10 Boxes/Case

# AMSure® New Products

All products are DEHP-Free, PVC-Free, Latex-Free, Preservative-Free, and Non-Pyrogenic

## AMSure® Sterile Water Prefilled Humidifier, USP

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- Prefilled humidifiers
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



## AMSure® Sterile Water for Inhalation, USP

- 500 mL and 750 mL, and 1000 mL neutralizer bottles Sterile Water, USP
- 12 bottles per case
- Sterile water USP for Inhalation Therapy
- Used in neutralizers and humidifiers
- 3 mL and 5 mL vials Sterile Water, USP
- Single-use vials with graduation
- Easy open twist-off for drop dispensing
- Used for inhalation and suctioning



## AMSure® Sodium Chloride Inhalation Solution (0.9%), USP

- Available in three sizes: 3 mL, 5 mL and 15 mL
- Embossed single-use vials with inverted milliliter graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



## AMSure® Device and Wound Irrigation

- Available in 100 mL Sterile Water, USP and also Sodium Chloride Solution 0.9%, USP bottle
- Easy open twist screw top container
- Used for device and wound irrigation

## Order Information

Item Code	Description	Qty/Unit
<b>AS0352</b>	350 mL Sterile Water for Inhalation Prefilled Humidifier	20/Case
<b>AS0552</b>	550 mL Sterile Water for Inhalation Prefilled Humidifier	12/Case
<b>AS1565</b>	500 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1765</b>	750 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1065</b>	1000 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1066</b>	1000 mL Sterile Water with Screw Cap for Inhalation, USP	12/Case
<b>ASR0035</b>	3 mL Sterile Water for Inhalation, USP	1000/Case
<b>ASR0055</b>	5 mL Sterile Water for Inhalation, USP	1000/Case
<b>ASR0039</b>	3 mL Sodium Chloride Inhalation Solution (0.9%), USP	1000/Case
<b>ASR0059</b>	5 mL Sodium Chloride Inhalation Solution (0.9%), USP	1000/Case
<b>ASR0159</b>	15 mL Sodium Chloride Inhalation Solution (0.9%), USP	144/Case
<b>AS399</b>	100 mL Sterile Water for Irrigation, USP	48/Case
<b>AS398</b>	100 mL Sodium Chloride Solution 0.9%, USP	48/Case



**Amsino International, Inc.**

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • info@amsino.com

To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSafe® and AMSure® are registered trademarks of Amsino International, Inc. © Copyright 2021, Amsino International, Inc. All rights reserved.