

# AMSafe® Pre-Filled Saline Flush Syringe

## 0.9% Sodium Chloride Injection, USP

### Features

- DEHP-free, PVC-free, not manufactured with natural rubber latex
- Sterile, non-pyrogenic
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared



## 12 mL and 20 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

Item Code	Description	Qty/Unit
<b>IVF1203</b>	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1205</b>	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1210</b>	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF2020</b>	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

### Sterile Field Ready Flush Syringes

Item Code	Description	Qty/Unit
<b>IVF1203SS</b>	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1205SS</b>	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1210SS</b>	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF2020SS</b>	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)



## 10 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

Item Code	Description	Qty/Unit
<b>IVF1003</b>	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1005</b>	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1010</b>	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
<b>IVF1010T</b>	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1010TM</b>	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

# AMSure® Water & Saline Medical Device Products

## AMSure® Sterile Water Prefilled Humidifier

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



## AMSure® Sterile Water for Inhalation

- 500 mL and 750 mL, and 1000 mL Sterile Water, USP
- Sterile water USP for Inhalation Therapy
- Used in nebulizers and humidifiers



## AMSure® Sterile Water and Sodium Chloride Inhalation Solution (0.9%)

- Available in three sizes: 3 mL, 5 mL and 15 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



## AMSure® Device and Wound Irrigation

- Available in 100 mL Sterile Water, USP and Sodium Chloride Solution 0.9%, USP
- Easy open twist screw top container
- Used for device and wound irrigation



## AMSure® Hanging Spikeable Bottle

- Sterile Water for Inhalation Therapy
- For use with active humidification systems
- Use with vented spikes only

## Order Information

Item Code	Description	Qty/Unit
<b>AS0352</b>	350 mL Sterile Water for Inhalation Prefilled Humidifier	20/Case
<b>AS0552</b>	550 mL Sterile Water for Inhalation Prefilled Humidifier	12/Case
<b>AS1565</b>	500 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1566</b>	500 mL Sterile Water Inhalation w/non-sterile cap	12/Case
<b>AS1765</b>	750 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1065</b>	1000 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1066</b>	1000 mL Sterile Water Inhalation w/non-sterile cap	12/Case
<b>ASR0035</b>	3 mL Sterile Water for Inhalation, USP	20/Case
<b>ASR0055</b>	5 mL Sterile Water for Inhalation, USP	12/Case
<b>ASR0039</b>	3 mL Sodium Chloride Inhalation Solution (0.9%), USP	1000/Case
<b>ASR0059</b>	5 mL Sodium Chloride Inhalation Solution (0.9%), USP	1000/Case
<b>ASR0159</b>	15 mL Sodium Chloride Inhalation Solution (0.9%), USP	144/Case
<b>AS399</b>	100 mL Sterile Water for Irrigation, USP	48/Case
<b>AS398</b>	100 mL Sodium Chloride Solution 0.9%, USP	48/Case
<b>AS1000</b>	1000 mL, Sterile Water for inhalation, Hanging Spikeable Bottle	12/Case



### AMSINO International, Inc.

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • info@amsino.com

To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSafe® and AMSure® are registered trademarks of AMSINO International, Inc. © Copyright 2022, AMSINO International, Inc. All rights reserved.