

AMSafe-3[®]

Designed for any I.V. environment

AMSafe-3[®] is a 3-in-1 multi-drip chamber I.V. set that allows providers to select a **10, 15 or 60 drop setting**. AMSafe-3[®] is Latex-free, DEHP-free and available with or without Pre-Attached Extension Sets.

- Drop setting can be adjusted during medication delivery
- Designed for efficiency where inventory reduction and cost containment are essential
- Selectable drop settings eliminate the need to carry different I.V. sets



AMSafe-3[®] is **NOW** available with two types of needle-free connectors

Both types of needle-free connectors have clear housings for visible flush confirmation. Their smooth access surface allows for effective disinfection and cleaning.








AMSafe[®] Needle-Free Connector






MicroClave[®] Clear Neutral Displacement, Needle-Free Connector

AMSafe-3® with AMSafe® Needle-Free Connectors

Item Code	Description	Length	Priming Volume (Approx.)	Qty/Unit
D38301	 I.V. Admin. Set, Roller Clamp, Pre-Pierced and Needle-Free Y-sites, and Rotating Male Luer Lock	83" (211 cm)	28 mL	50/Case
D38302	 I.V. Admin. Set, Roller Clamp, (2) Pre-Pierced and (2) Needle-Free Y-sites, Pre-Attached Extension Set, Slide Clamp, and Rotating Male Luer Lock	91" (231 cm)	30 mL	50/Case
D38303	 I.V. Admin. Set, Roller Clamp, Backcheck Valve, (1) Pre-Pierced Y site, (1) Needle-Free Y-site, Pre-Attached Ext Set with Pinch Clamp and Rotating Male Luer	98" (249 cm)	32 mL	50/Case
D38304	 I.V. Admin. Set, Roller Clamp, Backcheck Valve, (2) Pre-Pierced Y-sites, Pre-Attached Ext Set, Pinch Clamp, (1) Needle-Free Y-site and Rotating Male Luer	103" (262 cm)	34 mL	50/Case
D38305	 I.V. Admin. Set, Roller Clamp, (2) Pre-Pierced and (2) Needle-Free Y-sites, Backcheck Valve, Pre-Attached Extension Set, Slide Clamp, and Rotating Male Luer Lock	91" (231 cm)	30 mL	50/Case

AMSafe-3® with MicroClave® Clear Needle-Free Connectors

Item Code	Description	Length	Priming Volume (Approx.)	Qty/Unit
MR8301	 I.V. Admin. Set, Roller Clamp, Pre-Pierced and MicroClave Clear Y-sites, and Rotating Male Luer	83" (211 cm)	28 mL	50/Case
MR8302	 I.V. Admin. Set, Roller Clamp, (2) Pre-Pierced and (2) MicroClave Clear Y-sites, Pre-Attached Extension Set, Slide Clamp, and Rotating Male Luer	91" (231 cm)	30 mL	50/Case
MR8305	 I.V. Admin. Set, Roller Clamp, (2) Pre-Pierced and (2) MicroClave Clear Y-Sites, Backcheck Valve, Pre-Attached Extension Set, Slide Clamp, and Rotating Male Luer	91" (231 cm)	30 mL	50/Case



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768
800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. MicroClave is a registered trademark of ICU Medical. AMSafe-3 and AMSafe are registered trademarks of Amsino International, Inc. © Copyright 2016, Amsino International, Inc. All rights reserved.