


# AMSafe-3<sup>®</sup>

## I.V. Solutions for an Ever-changing Environment

AMSafe-3<sup>®</sup> is a 3-in-1 multi-drip chamber I.V. set that allows for a selectable **10, 15 or 60 drop setting**.

- Adjustable drop size eliminates the need to carry different I.V. sets, allowing for efficient care when you need it most
- Drip chamber is attached directly to the bag spike for easy access, adjustments and drip chamber visibility
- Available with Sure-lok<sup>™</sup> Needle-Free and Pre-pierced Y-sites in an array of configurations
- DEHP-Free, Latex-Free, and Non-Pyrogenic







## AMSafe-3<sup>®</sup> is available with the Sure-Lok<sup>™</sup> Needle-Free Y-site

- Low resistance to luer connection and disconnection for easy administration
- Smooth flat top with tight septum seal to support effective disinfection
- High flow rate for efficient administration
- Clear housing for visual confirmation of flush



Sure-Lok<sup>™</sup>  
Needle-Free Y-site

## AMSafe-3® with Sure-Lok™ Needle-Free Y-site

Item Code	Description	Length	Priming Volume (Approx.)	Qty/Unit
<b>D38301</b>	 I.V. Admin. Set, Roller Clamp, Pre-Pierced and Sure-Lok™ Needle-Free Y-sites, and Rotating Male Luer Lock	83 in. (211 cm)	28 mL	50/Case
<b>D38302</b>	 I.V. Admin. Set, Roller Clamp, Pre-Pierced and Sure-Lok™ Needle-Free Y-sites, Pre-Attached Extension Set, Pre-Pierced Y-site, Slide Clamp, Sure-Lok™ Needle-Free Y-site, and Rotating Male Luer Lock	91 in. (231 cm)	30 mL	50/Case
<b>D38305</b>	 I.V. Admin. Set, Roller Clamp, Pre-Pierced and Sure-Lok™ Needle-Free Y-sites, Backcheck Valve, Pre-Attached Extension Set, Pre-Pierced Y-site, Slide Clamp, Sure-Lok™ Needle-Free Y-site, and Rotating Male Luer Lock	91 in. (231 cm)	30 mL	50/Case
<b>D38307</b>	 I.V. Admin. Set, Roller Clamp, Pre-Pierced Y-site with Pre-Attached Extension Set with Pinch Clamp, Sure-Lok™ Needle-Free Y-site, and Rotating Male Luer Lock	82 in. (208 cm)	27 mL	50/Case

To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)

**AMSINO**

Advancing Healthcare Worldwide®

**Amsino International, Inc.**

708 Corporate Center Drive, Pomona, CA 91768

800.632.6746 • fax: 909.626.3888 • [info@amsino.com](mailto:info@amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSafe-3 and Sure-Lok are registered trademarks of Amsino International, Inc. © Copyright 2018, Amsino International, Inc. All rights reserved.