


AMShield® Disinfecting Cap

The Power of 70% Isopropyl Alcohol (IPA)


In an effort to protect the safety of patients, Amsino is launching the AMShield® Disinfecting Cap (patented design). Most needle-free valves are unprotected against airborne contaminants in the hospital setting and other contaminants that can create risks to patients. With the unique design of a compressed foam sphere, utilizing 70% Isopropyl Alcohol acting as a disinfectant, the needle-free valves are protected against bacterial ingress. The connector surface and threads are also protected by the disinfecting agent.

Follow hospital protocol and utilize standard technique to swab the connector. The AMShield® Disinfecting Cap can be rotated on the connector due to no thread which allowing nurses to “scrub the hub”. This is especially necessary when the connector is heavily contaminated, for example, from blood residuals. Utilizing Disinfecting Caps should become a standard of care.



Features

- Fits all needle-free valves
- Stays securely in place
- Contains 70% Isopropyl Alcohol
- Helps reduce CLABSI (central line-associated bloodstream infections)
- Available individually packaged or as IV pole multi-strips (10caps/strip) to fit an institution's requirements
- Bright gold color is easy to identify when installed on the IV line for visual verification of compliance monitoring
- Sterile, DEHP-free, latex-free, single-use-only



Item Code	Description	Quantity/Unit	
DCF2000	AMShield® Disinfecting Cap, individually packaged	200 Caps/Box	10 Boxes/Case
DCF2500S	AMShield® Disinfecting Cap, 10 Caps mounted on a strip	25 Strips/Box	10 Boxes/Case



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768
800.632.6746 • fax: 909.626.3888 • info@amsino.com

To view our complete product line,
visit us at www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMShield® is a registered trademark of Amsino International, Inc. © Copyright 2020, Amsino International, Inc. All rights reserved.