


AMSure®

Enteral Syringes with ENFit® Tip

AMSure® Enteral Syringes are available with ENFit® Transition Connector if requested on your PO. The ENFit® Transition Connector is compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.



ENFit® Tip



ENFit® Transition Connector

AMSure® Thumb Control Ring Syringe with ENFit® Tip

AMSure® Thumb Control Ring Syringe with ENFit® tip is designed for maximum control allowing for a one handed administration. The syringe tip is compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.



AMSure® Flat Top Piston Syringe with ENFit® Tip

AMSure® Flat Top Piston Syringe with ENFit® tip has visible dosing marks designed for ease of use. The syringe tip is compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.









AMSure® Enteral Feeding Irrigation Kits with ENFit® Syringe

AMSure® Enteral Feeding Irrigation Kit is designed for simplicity and efficiency of care. The non-sterile kit includes a 60 mL thumb control ring syringe with ENFit® tip, patient label, 500 mL graduated container and either a poly pouch or resealable bag. The syringe tip is compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.



For a limited time, Amsino is offering a significant price reduction on these items. Supplies are limited, please contact your Amsino representative for pricing.

Order Information

Item Code	Image	Description	Features	Qty/Unit
ENS015		AMSure® 60 mL Thumb Control Ring Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS015NS ENS016 equivalent		AMSure® 60 mL Thumb Control Ring Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENS115		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS115NS ENS116 equivalent		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENK126		AMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Control Ring Piston Syringe with ENFit® Tip, Patient Label	Sterile Poly Pouch	30/Case
ENK127		AAMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Control Ring Piston Syringe with ENFit® Tip, Patient Label	Non-Sterile Resealable Pole Bag*	30/Case

*Resealable pole bag designed for safe and convenient storage.

All syringes include a tip protector and are latex-free, DEHP-free and BPA-free.

To view our complete product line, visit us at www.amsino.com



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768

800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure® is a registered trademark of Amsino International, Inc. © Copyright 2020, Amsino International, Inc. All rights reserved.