



Amsino International, Inc.  
708 Corporate Center Drive  
Pomona, CA 91768  
USA

Tel: (909) 626-5888 Fax: (909) 626-9006  
Toll Free: 1-800-63-amsino  
E-mail: [amsino@amsino.com](mailto:amsino@amsino.com)  
[www.amsino.com](http://www.amsino.com)

May 15, 2017

Dear Valued Customer,

The AS016 Thumb Control Piston Irrigation Syringe, 60cc is a non-sterile, Class I Medical Device that does not required an Expiration Date



The date printed on the labeling net to the symbol is the date of manufacture.

The recommended shelf life is 5 years (60 months) from the manufacture date.

Sincerely,

Holiven Ji  
Director of QA & RA