



# AMSure® High Flow Oxygen Nasal Cannula

AMSure® high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

## Features

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most High Flow Nasal Cannula devices
- Single-patient-use, Latex free
- Sizes available: S/M/L Adult

## Nasal Cannula



- A Adjustable angle**  
Rotatable silicone clip for easy angle adjustment
- B Slide-in pliable nasal prongs**
- C An additional clip** is provided for the overhead strap to carry the tube on the patient's clothes
- D Swivel connector**  
Allows for moving direction easily



## Order Information

Item Code	Description	Qty/Unit
<b>AS75280</b>	AMSure® Nasal Cannula, Size: Small-Adult	10/Case
<b>AS75281</b>	AMSure® Nasal Cannula, Size: Medium-Adult	10/Case
<b>AS75282</b>	AMSure® Nasal Cannula, Size: Large-Adult	10/Case



## Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • info@amsino.com

To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2021, Amsino International, Inc. All rights reserved.