

AMSure®

Closed Suction Catheter

24-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

Infection Control


- Closed system to prevent cross infections in facility
- Clean easily with needle free port

Safety:

- Avoid interruptions of oxygen supply due to suction
- Safety lock to prevent un-expected suction
- Clear depth mark for accurate insertion
- R-polished catheter tip to reduce respiratory tract injury

Convenience:

- Reduce additional step to separate from ventilator before suction
- Flexible connection to ventilator without moving patient




Clear Mark on Vacuum Control Easy for Nurses to know if suction is Closed or Open (Safe for patients).

Graduation Marks on Catheter Easy to identify the depth when catheter is inserted.



Integrated closed structure Effectively protects patients and medical staff




Color Coded buttons as per ISO 8836:2009

Smooth Catheter tip



Amsino provides an abundant portfolio for closed suction system



MDI Port for Drug Delivery

Double Swivel or T-Piece Options to Connect to Ventilator

Flex Connector Easy to Accommodate Ventilator Circuits position

Endotracheal Kits with MDI Port Adapter

Item Codes	Size (Fr)	Double Swivel Elbow	T-Piece	Endotracheal	Flex Connector	MDI Port	QTY/Case
VS363RC124	10	✓		54cm (21.3 in)	15 mm x 22 mm	✓	20/case
VS364RC124	12	✓		54cm (21.3 in)	15 mm x 22 mm	✓	20/case
VS365RC124	14	✓		54cm (21.3 in)	15 mm x 22 mm	✓	20/case
VS363TP124	10		✓	54cm (21.3 in)	22 mm x 22mm	✓	20/case
VS364TP124	12		✓	54cm (21.3 in)	22 mm x 22 mm	✓	20/case
VS365TP124	14		✓	54cm (21.3 in)	22 mm x 22 mm	✓	20/case

Tracheostomy Kits with MDI Port Adapter

Item Codes	Size (Fr)	Double Swivel Elbow	T-Piece	Tracheostomy	Flex Connector	MDI Port	QTY/Case
VS363TC124	10	✓		30.5cm (12 in)	15 mm x 22 mm	✓	20/case
VS364TC124	12	✓		30.5cm (12 in)	15 mm x 22 mm	✓	20/case
VS365TC124	14	✓		30.5cm (12 in)	15 mm x 22 mm	✓	20/case
VS363TTP124	10		✓	30.5cm (12 in)	22 mm x 22 mm	✓	20/case
VS364TTP124	12		✓	30.5cm (12 in)	22 mm x 22 mm	✓	20/case
VS365TTP124	14		✓	30.5cm (12 in)	22 mm x 22 mm	✓	20/case

***All products are 10 PCS/Box, 2 Boxes/Case**



Amsino International, Inc.
 708 Corporate Center Drive, Pomona, CA 91768
 800.632.6746 • fax: 909.626.3888 • info@amsino.com

For the full Amsino product line, visit us at www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©Copyright 2020, Amsino International, Inc. All rights reserved.