

2018 PRODUCT CATALOG



AMSINO

Advancing Healthcare Worldwide®

TABLE OF CONTENTS

GENERAL ORDERING INFORMATION	ii
AMSure® ENTERAL FEEDING SYSTEMS.....	1
Enteral Feeding Pump.....	2
Enteral Feeding Sets	2
Enteral Syringes.....	3
Enteral Feeding Irrigation Kits and Trays.....	4
AMSafe® MEDICATION DELIVERY SYSTEMS.....	5
Gravity Administration Sets	6
10 Drops/mL	6
15 Drops/mL	8
60 Drops/mL	11
Selectable Drip Chamber	
AMSafe-3® 10, 15 and 60 drops/mL	12
Secondary Piggyback Sets	13
Y-type Blood Sets	13
I.V. Extension Sets	13
Standard Extension Sets	13
Microbore Extension Sets.....	14
Smallbore Extension Sets	15
Flow Controller Extension Sets	15
Pressure Infusion Extension Sets	16
Needle-free Connectors	16
Bag Spike	16
Needle-free Vial Adapter	16
AMSure® UROLOGICAL CARE	17
Urinary Collection Devices	18
Urinary Drainage Bags	18
Urinary Leg Bags	19
Extension Tube for Leg Bags.....	19
Pediatric Urine Collector	19
Foley Catheters	20
Foley Insertion Trays	21
Pre-Connected Drainage Bag Foley Catheter Trays	22
Urethral Catheters.....	23
Urethral Insertion Trays and Kits.....	24
Pre-Connected Drainage Bag Urethral Catheter Kit	25
Self-Catheterization Kits	26
Irrigation Trays and Kits	27
Syringes	28
Specimen Containers	29
AMSure® RESPIRATORY THERAPY SYSTEMS	30
Respiratory Masks	31
Suction Catheters, Trays and Kits.....	33
Nebulizer Kits.....	34
Aerosol Drainage Bag	35
Tracheostomy Clean and Care Trays.....	35
Suction Tubing	36
Suction Yankauers and Kits	37
Suction Device Holster and Kits.....	38
Accessories	39
AMSure® SURGICAL AND NURSING SYSTEMS....	40
Suction Tubing	41
Suction Yankauer and Kits	42
Suction Device Holster and Kits.....	43
Ear/Ulcer Syringes.....	43
Vaginal Speculums	44
Disposable Umbilical Cord Clamp	44
Specimen Containers	44
Cleansing Enema Sets	45
5-in-1 Straight Connector	45
RECEPTAL® Waste Fluid Management Systems	46
RECEPTAL® Canisters	47
RECEPTAL® Liners	47
EZE-VAC® Suction Canister System	48
ACCU-MEASURE® Critical Measurement Devices	48
Suction Tubing	49
Suction Yankauers and Kits	50
Suction Device Holster and Kits	51
Accessories	52
Reusable Hardware	54
Floor Stands and Trolleys	54
Canister Supports	55
Regulator Mounts.....	56
INDEX OF ITEM CODES.....	57

GENERAL ORDERING INFORMATION

CONTACT INFORMATION

- Phone: (909) 626-5888 or (800) MD-AMSINO, M-F – 7:00 am to 5:00 pm, PST
- Fax: (909) 626-3888 – 24 hours a day
- Address: Customer Service, 708 Corporate Center Drive, Pomona, CA 91768
- E-mail: info@amsino.com
- Website: www.amsino.com

TERMS OF PAYMENT & INVOICE BILLING

- Please pay your invoice and make checks payable to Amsino International, Inc.
- Invoices not paid according to terms of sale are subject to a 1.5% late payment penalty per month.

RETURNS

- Authorization to return products must be obtained from Amsino's Customer Service Department prior to the return of any goods. Failure to obtain authorization will result with denial of the return.
- To obtain authorization the following must be provided:
 - Product number
 - Lot number
 - Quantity
 - Reason for the return
- Amsino's Customer Service Department will provide a Return Authorization Number (RGA). This number must be on the packing slip of the returned products in order to receive credit.
- Product returns are subject to the following fees:
 - 20% restocking fee
 - Return freight fee
 - \$100 handling fee per restocking shipment

PRODUCTS NOT ELIGIBLE FOR RETURNS

- Unauthorized returns
- Products that have been mishandled or abused
- Opened cartons
- Custom products
- Discontinued products
- Products purchased more than 90 days prior to return request

MINIMUM ORDER

- \$250.00 minimum
- \$25.00 handling fee for orders less than minimum
- \$2,000.00 minimum first time order for new accounts

FREIGHT

- Contact your sales representative for freight minimums.
- Fees will apply to drop shipment orders.
- Freight & MOQ policies will apply.
- All pump orders must have separate PO and FOB Amsino.
- Freight subject to fuel charge.

DROP SHIPMENT

- Fees will apply to drop shipment orders. Freight & MOQ policies will apply.

NOTES

- AMSINO®, AMSure®, AMSafe® and AMSafe-3® are registered trademarks of Amsino International, Inc.
- Product availability and pricing is subject to change without prior notice.
- Concealed Damaged Shipments must be reported to Amsino along with a copy of the inspection report or waiver within 3 business days from receipt of product. According to the National Motor Freight Classification Supplement 49 4949 U.S.C. 14706 and 49 CFR Part 370, 30013, all Concealed Damages must be reported into the appropriate freight carrier no more than (5) five days from delivery of shipment.

AMSure® ENTERAL FEEDING SYSTEMS

OPTIMIZING NUTRITIONAL DELIVERY

AMSure® Enteral Feeding System enhances delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The system is designed to provide best patient outcomes and includes:

- Sentinel® Enteral Feeding Pump 2
- ENFit® and Non-ENFit® Sets 2
- ENFit® and Non-ENFit® Syringes 3
- ENFit® and Non-ENFit® Irrigation Kits and Trays 4



ENTERAL FEEDING PUMP

The Sentinel® enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel® pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications



ITEM CODE	DESCRIPTION	QUANTITY/UNIT
S-1000-S1	Sentinel® Enteral Feeding Pump	1/Each

ENTERAL FEEDING SETS

AMSure® Enteral Feeding Sets are available with and without ENFit® connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure® enteral feeding pump sets are recommended to use with the Sentinel® enteral pump and are compatible with the Compat™, EntraFlo™ and Kangaroo™ 224, 324 and PET pumps.

- ENFit® connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit® sets include pre-attached stepped transition connector
- Non-sterile, latex-free and BPA-free



ENF1200K

ITEM CODE	ENFit® CONNECTOR	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY			QUANTITY/ UNIT
		Spike Only- No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Kangaroo™ 224, 324 and PET	Compat™/ EntraFlo™	
E-0500M			●			●	●	●		30/Case
E-1004	●					●	●		●	30/Case
E-1200				●		●	●		●	30/Case
E-1200M				●		●	●	●		30/Case
E-1203				●	●					30/Case
ENF1004K	●	●				●	●		●	30/Case
ENF1200K	●			●		●	●		●	30/Case
ENF1203K	●			●	●					30/Case
ENF500K	●		●			●	●		●	30/Case

ENTERAL SYRINGES

AMSure® Enteral Syringes are available with or without an ENFit® tip. The ENFit® syringe tip is compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are latex-free, DEHP-free and BPA-free.

BULB SYRINGE

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination



AS010

THUMB CONTROL RING SYRINGE

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage



AS016 ENS016

FLAT TOP SYRINGE

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use

AS116 ENS115



ITEM CODE	ENFIT® TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QUANTITY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010			●						●		50/Case
AS011			●			●				●	50/Case
AS011P			●				●			●	50/Case
AS015	●			●		●				●	50/Case
AS016	●			●				●			30/Case
AS115	●				●	●				●	50/Case
AS116	●				●			●			30/Case
ENS015	●			●		●					50/Case
ENS016	●			●				●		●	30/Case
ENS115	●				●			●			50/Case
ENS116	●				●					●	30/Case

ENTERAL FEEDING IRRIGATION KITS AND TRAYS

AMSSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60-mL syringe with or without an ENFit® tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60-mL syringe with or without an ENFit® tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- Latex-free, DEHP-free and BPA-free
- Includes syringe tip protector cap

AS121



AS136



KITS

ITEM CODE	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL RING PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QUANTITY/ UNIT
AS121	●	●					●	30/Case
AS126	●			●			●	30/Case
AS127	●			●		●	●	30/Case
ENK126	●	●		●	●		●	30/Case
ENK127	●	●		●		●	●	30/Case

TRAYS

ITEM CODE	ENFIT® SYRINSGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL RING PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QUANTITY/ UNIT
AS130	●	●	●	●		●	●	●	20/Case
AS136	●	●	●		●	●	●	●	20/Case
ENK136	●	●	●		●	●	●	●	20/Case

AMSafe® MEDICATION DELIVERY SYSTEMS

SAFETY FIRST

The AMSafe® I.V. product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our I.V. product selection includes:

• Gravity Administration Sets	6
10 Drops/mL	6
General Primary Sets	6
Flow Controller I.V. Sets	7
Piggyback Sets with Back Check Valve	8
15 Drops/mL	8
General Primary Sets	8
Piggyback Sets with Back Check Valve	10
60 Drops/mL.....	11
General Primary Sets	11
Piggyback Sets with Back Check Valve	11
Burette Sets	12
Selectable Drip Chamber AMSafe-3® 10, 15 and 60 drops/mL	12
Secondary Piggyback Sets	13
Y-type Blood Sets	13
• I.V. Extension Sets	13
Standard Extension Sets	13
Microbore Extension Sets	14
Smallbore Extension Sets	15
Flow Controller Extension Sets.....	15
Pressure Infusion Extension Sets	16
• Needle-free Connectors	16
• Bag Spike	16
• Needle-free Vial Adapter	16

GRAVITY ADMINISTRATION SETS

Amsino's comprehensive line of gravity administration I.V. sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in an array of sets to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Non-pyrogenic
- DEHP-free, latex-free and sterile

10 DROPS/mL

GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
108303	10	83 in. (211 cm)	Sure-Lok™ Needle-free Y-site	16 mL	50/Case
109602	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
108306	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
MRSA10-84	10	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	10.5 mL	50/Case
AC3110	10	105 in. (266 cm)	(2) Sure-Lok™ Needle-free Y-sites	20 mL	50/Case
MRA10E-95	10	95 in. (241 cm)	Vented/Non-vented, Pre-pierced Y-site, Slide Clamp, Pre-attached Ext Set, Slide Clamp, MicroClave® Clear Needle-free Y-site	18.5 mL	50/Case

10 DROPS/mL

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
108305	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
MRA10E-104	10	104 in. (264 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-sites, Pre-pierced Y-site, with Pre-attached Ext Set, Slide Clamp, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	13.4 mL	50/Case
MRSA10E-107	10	107 in. (272 cm)	Vented/Non-vented, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set	13.4 mL	50/Case
108308	10	83 in. (211 cm)	(2) Sure-Lok™ Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case




FLOW CONTROLLER I.V. SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AFS102	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok™ Needle-free Y-site	18.5 mL	50/Case
AFS208	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites	19.2 mL	50/Case
AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe® Needle-free Y-sites, Slide Clamp	19 mL	50/Case

10 DROPS/mL



PIGGYBACK SETS WITH BACK CHECK VALVE

- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
1010001B	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok™ Needle-free Y-site	19.3 mL	50/Case
					
108302	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
					
MRA102-84	10	84 in. (213 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA103-94	10	94 in. (239 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (3) MicroClave® Clear Needle-free Y-sites	17.5 mL	50/Case
					
MRSA10E-110	10	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, (2) MicroClave® Clear with Pre-attached Ext Set	21.1 mL	50/Case
					

15 DROPS/mL

GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AA3101	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
158310	15	83 in. (211 cm)	Vented/Non-vented, Pre-pierced Y-site	17 mL	50/Case
					







15 DROPS/mL

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
158309	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok™ Needle-free Y-site	17 mL	50/Case
MRSA15-100	15	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site	12.4 mL	50/Case
158306	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
MRSA15-91	15	91 in. (231 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	11.4 mL	50/Case
159602	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
MRA15-100	15	100 in. (254 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	13 mL	50/Case
1510501	15	105 in. (267 cm)	Vented/Non-vented, (2) Sure-Lok™ Needle-free Y-sites	20 mL	50/Case
158305	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
MRSA15E-88	15	88 in. (224 cm)	Vented/Non-vented, Pinch Clamp, Pre-pierced Y-site with Pre-attached Ext Set, Pinch Clamp, (2) MicroClave® Clear Needle-free Y-sites	11.1 mL	50/Case
158307	15	89 in. (226 cm)	Sure-Lok™ Needle-free Y-site, Pre-pierced Y-site with Pre-attached Ext Set, Sure-Lok™ Needle-free Y-site, Slide Clamp, Pre-pierced Y-site	18.7 mL	50/Case
MRSA15E-120	15	120 in. (305 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site, Rotating Male Luer Lock, Drop-In Ext Set MicroClave® Clear Y-site	15.8 mL	50/Case

15 DROPS/mL

PIGGYBACK SETS WITH BACK CHECK VALVE

- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a backcheck valve to prevent fluid exchange or backflow into the primary infusion line

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRSA15E-96	15	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set, Slide Clamp	12.1 mL	50/Case
					
MRA15B-105	15	105 in. (267 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites	13.23 mL	50/Case
					
AS3100	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
AS3101	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
AS3102	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
AS3110A	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok™ Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					
AS3110B	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Pre-pierced Y-sites, Slide Clamp	22 mL	50/Case
					

60 DROPS/mL

GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
607201	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608306	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
					
608305	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
MRSA60-84	60	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA60E-88	60	88 in. (224 cm)	Vented/Non-vented, Pre-pierced Y-site with Pre-attached Ext Set, (2) MicroClave® Clear Needle-free Y-sites	11.1 mL	50/Case
					
MRSA60-100	60	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site, Pre-pierced Y-site	12.5 mL	50/Case
					

PIGGYBACK SETS WITH BACK CHECK VALVE

- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

MRA60B-92	60	92 in. (234 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Slide Clamp	11.4 mL	50/Case
					
MRSA602-100	60	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites	12.5 mL	50/Case
					

60 DROPS/mL

BURETTE SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRS60B-106	60	106 in. (269 cm)	150 mL Burette Set (w/MicroClave®, Shut-Off), (2) MicroClave® Clear Needle-free Y-sites	12.7 mL	25/Case
MRS60B-117	60	117 in. (297 cm)	24-hour 150 mL Burette Set (w/15 Micron Filter, Shut-off) with (3) MicroClave® Clear Needle-free Y-sites	13.5 mL	25/Case

SELECTABLE DRIP CHAMBER AMSAFE-3® 10, 15 AND 60 DROPS/mL

AMSafe-3® is a 3-in-1 multi-drip chamber I.V. set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential. Now you can stock one I.V. set for multiple drip rates.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- DEHP-free, latex-free and sterile

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
D38301	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	27.7 mL	50/Case
D38302	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok™ Needle-free Y-site	30.2 mL	50/Case
D38305	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Back Check Valve with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok™ Needle-free Y-site	30.5 mL	50/Case


SECONDARY PIGGYBACK SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AS3515	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
AS3500	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case




Y-TYPE BLOOD SETS


ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRSB-94	10	94 in. (239 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Slide Clamp, Roller Clamp, MicroClave® Y-site	57.2 mL	25/Case



MRSB-102	20	102 in. (259 cm)	200 Micron Filter, Roller Clamps, Hand Pump, MicroClave® Y-site	20.2 mL	25/Case
----------	----	------------------	---	---------	---------



MRSB-105	10	105 in. (267 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Pre-pierced Y-site, Slide Clamp, Roller Clamp, MicroClave® Y-site, Pre-attached Ext Set	59.3 mL	25/Case
----------	----	------------------	--	---------	---------




I.V. EXTENSION SETS

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, latex-free, sterile and include a rotating male luer lock unless otherwise noted.



- Standard Extension Sets
- Smallbore Extension Sets
- Pressure Infusion Extension Sets
- Microbore Extension Sets
- Flow Controller Extension Sets

STANDARD EXTENSION SETS



ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
000606	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Slide Clamp	2.5 mL	50/Case




STANDARD EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MR4T07	6 in. (15 cm)	MicroClave® Clear Needle-free Connector, T-connector	0.86 mL	50/Case
				
AE3107	7 in. (18 cm)	Removeable Sure-Lok™ Needle-free Connector, Slide Clamp	1.76 mL	50/Case
				
MR4107	7 in. (18 cm)	Removable MicroClave® Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
MR4207	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
AE3108	8 in. (20 cm)	Removable Sure-Lok™ Needle-free Connector, Pinch Clamp	1.76 mL	100/Case
				
MR4R10	10 in. (25 cm)	Removable MicroClave® Clear Needle-free Connector, Pinch Clamp	0.98 mL	50/Case
				
003201	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
				
003202	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok™ Needle-free Y-sites	5.3 mL	50/Case
				
MR4T32	32 in. (81 cm)	(2) MicroClave® Clear Needle-free Connectors, (2) Slide Clamps	5.02 mL	50/Case
				



MICROBORE EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AE0106	6 in. (15 cm)	Removeable Sure-Lok™ Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
AE1106	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
				

MICROBORE EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AE0107	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok™ Needle-free Y-site	0.9 mL	50/Case
				
AE1112	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
AE1138	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case
				

SMALLBORE EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MR4307	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.24 mL	50/Case
				
MR4B07	7 in. (18 cm)	Bifuse Ext Set with (2) MicroClave® Clear Needle-free Connector, (2) Slide Clamps	0.35 mL	50/Case
				
MR4TF04	4 in. (10 cm)	Trifuse Ext Set with (3) MicroClave® Clear Needle-free Connectors, (3) Slide Clamps	0.40 mL	50/Case
				

FLOW CONTROLLER EXTENSION SETS

Amsino's I.V. extension sets with flow control devices are designed to control I.V. fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AFS107	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
				
MRSFC-19	19 in. (48 cm)	Flow Controller, MicroClave® Clear Needle-free Connector, Clamp	3.2 mL	50/Case
				

PRESSURE INFUSION EXTENSION SETS

For rapid administration of medication or fluids in an emergency, Amsino's single-patient-use pressure infusion I.V. extension set is made to withstand up to 400 psig. Includes a rotating male luer lock.

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MR4P07	7 in. (18 cm)	Pressure Infusion (400psig) Ext Set with MicroClave® Clear Needle-free Connector, Purple Slide Clamp	0.43 mL	50/Case



NEEDLE-FREE CONNECTORS

Prevent needle stick injuries and reduce bacterial contamination with needle-free connectors.


- Swabbable
- MRI-friendly
- Not limited to 24-hour use
- DEHP-free, latex-free and sterile

ITEM CODE	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MR4001	MicroClave® Clear Needle-free Pressure Infusion (400psig) Connector	0.04 mL	100/Case
AY0200	AMSafe® Needle-free Connector	0.04 mL	100/Case




BAG SPIKE

ITEM CODE	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MBSS001	MicroClave® Needle-free Connector I.V. Bag Spike	0.34 mL	50/Case



NEEDLE-FREE VIAL ADAPTER

ITEM CODE	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MB4V01	MicroClave® Clear Needle-free Connector Single Vial Adapter	0.16 mL	50/Case



AMSure® UROLOGICAL CARE

SOLUTIONS FOR INFECTION PREVENTION

AMSure® Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

• Urinary Collection Devices	18
• Foley Catheters	20
• Foley Insertion Trays	21
• Urethral Catheters.....	23
• Urethral Trays and Kits	23
• Self-catheterization Kits	26
• Irrigation Trays and Kits	27
• Syringes	28
• Specimen Containers	29



URINARY COLLECTION DEVICES

URINARY DRAINAGE BAGS

AMSure® Urinary Drainage Bags are designed with needle-free sampling ports and anti-reflux devices to improve safety and infection prevention. AMSure® urinary drainage bags are intended to collect urine from sheaths or catheters.

- Anti-reflux mechanism helps improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, latex-free and BPA-free

ITEM CODE	BAG SHAPE		ANTI-REFLUX TYPE			NEEDLE-FREE SAMPLING PORT		DRAIN TYPE		HANGER TYPE			QUANTITY/UNIT
	Low-profile	Teardrop	Chamber	Flutter Valve	Device	Pre-pierced Blunt Cannula	Male Luer Lock or Luer Slip	T-tap	Metal Clamp	Single Hook	Double Hook	Rope	
AS312		●	●			●		●			●	●	20/Case
AS322		●	●			●		●			●	●	20/Case
AS32200	●		●			●		●		●		●	20/Case
AS32600*	●		●				●	●	●	●		●	20/Case
AS332		●		●		●		●		●			20/Case
AS326	●				●	●			●		●	●	20/Case
AS302		●			●	●		●			●	●	20/Case

*100% sterile



AS322 Teardrop




AS32600 Low-profile

URINARY LEG BAGS

AMSSure® Urinary Leg Bags provide discrete and comfortable drainage. The bag is available in two sizes with pre-attached leg straps, drainage port and an anti-reflux device. Push-pull and twist-turn drainage port design options available for patient dexterity.

- Soft and breathable pre-attached leg straps provide patient comfort
- Anti-reflux device helps improve patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway, latex-free and BPA-free



ITEM CODE	600 mL LEG BAG	900 mL LEG BAG	PUSH-PULL DRAIN PORT	TWIST-TURN DRAIN PORT	18 IN. (45.72 CM) EXTENSION TUBE WITH ADAPTER AND PROTECTIVE CAP	QUANTITY/UNIT
AS306N	●		●			48/Case
AS307N		●	●			48/Case
AS308N	●			●		48/Case
AS309N		●		●		48/Case
LBE306N	●		●		●	48/Case
LBE307N		●	●		●	48/Case
LBE308N	●			●	●	48/Case
LBE309N		●		●	●	48/Case

EXTENSION TUBE FOR LEG BAGS

The 18 in. leg bag extension tubing is essential for comfort and proper placement of the leg bag. The leg bag can be positioned either on the thigh or calf. The extension tubing can be cut to the length the patient requires. It includes an adapter and protective cap.

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, latex-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case



PEDIATRIC URINE COLLECTOR

AMSSure® Pediatric Urine Collectors are hypoallergenic and provide graduation marks designed for convenient and simple urine collection.

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS409	200 mL, Pediatric Urine Bag	50/Box



FOLEY CATHETERS

2-WAY SILICONIZED LATEX FOLEY CATHETERS

AMSure® 2-Way Siliconized Latex Foley catheters are durable and flexible which help improve patient comfort. The Foley catheters are 16 in. (40.64 cm) long, available in a variety of French and balloon sizes.

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free and BPA-free



2-WAY 100% SILICONE FOLEY CATHETERS

AMSure® 2-way 100% Silicone Foley Catheters are ideal for patients with latex allergies and longer-term catherization.

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, latex-free, DEHP-free and BPA-free



ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS41012	12 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41014	14 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41016	16 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41018	18 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41020	20 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41022	22 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41024	24 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41026	26 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41028	28 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41030	30 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS42012	12 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42014	14 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42016	16 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42018	18 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42020	20 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42022	22 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42024	24 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42026	26 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42028	28 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42030	30 Fr, 16 in. (40.64 cm)	30 cc	10/Box

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS41012S	12 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41014S	14 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41016S	16 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41018S	18 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41020S	20 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41022S	22 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41024S	24 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41026S	26 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41028S	28 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS42016S	16 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42018S	18 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42020S	20 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42022S	22 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42024S	24 Fr, 16 in. (40.64 cm)	30 cc	10/Box

FOLEY INSERTION TRAYS

AMSure® Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter.

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, latex-free, DEHP-free and BPA-free



Trays include:

- 1000 mL collection tray
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Lubrication jelly
- Fenestrated drape

ITEM CODE	10 CC SYRINGE PRE-FILLED WITH STERILE WATER	30 CC SYRINGE PRE-FILLED WITH STERILE WATER	QUANTITY/UNIT
AS880	●		20/Case
AS890		●	20/Case

ADD-A-FOLEY CATHETER TRAY

AMSure® Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter.

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine solution help maintain sterility and reduce the risk of infection
- Sterile and latex-free

Trays include:

- 2000 mL low-profile urinary drainage bag
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone iodine solution
- Forceps
- Lubricating jelly
- 2 oz. specimen container with patient label
- Central supply room wrap and Tyvek header bag



ITEM CODE	10 CC SYRINGE PRE-FILLED WITH STERILE WATER	30 CC SYRINGE PRE-FILLED WITH STERILE WATER	QUANTITY/UNIT
AS89300	●		10/Case
AS89330		●	10/Case

PRE-CONNECTED DRAINAGE BAG FOLEY CATHETER TRAYS

AMSSure® Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to the low-profile drainage bag which helps to minimize contamination.

- Tray includes all necessary components for aseptic catheter insertion
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help maintain sterility and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Tray is available with both siliconized latex and 100% silicone catheters to suit patient needs
- Sterile

Pre-connected drainage bag Foley tray contents include:

- 2000 mL low-profile urinary drainage bag
- Pre-connected 16 in. Foley catheter with 5 cc balloon
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone iodine solution
- Forceps
- Lubricating jelly
- Syringe prefilled with 10 cc sterile water
- 2 oz. specimen container with a patient label
- Central supply room wrap and Tyvek header bag



AS89316

ITEM CODE	FOLEY CATHETER FRENCH SIZE			CATHETER MATERIAL		QUANTITY/ UNIT
	14 French	16 French	18 French	Siliconized Latex	100% Silicone (Latex-free)	
AS89314	●			●		10/Case
AS89316		●		●		10/Case
AS89318			●	●		10/Case
AS89314S	●				●	10/Case
AS89316S		●			●	10/Case
AS89318S			●		●	10/Case

URETHRAL CATHETERS

RED RUBBER LATEX URETHRAL CATHETER

AMSure® Red Rubber Urethral Catheters are pliable and designed for short-term use. The urethral catheters are available in various French sizes as well as female and male lengths.

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free and BPA-free




ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS44008	8 Fr, 11 in. (27.94 cm)	Female	100/Case
AS44010	10 Fr, 11 in. (27.94 cm)	Female	100/Case
AS44012	12 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44014	14 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44016	16 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44018	18 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44020	20 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44022	22 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44024	24 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44028	28 Fr, 16 in. (40.64 cm)	Male	100/Case

PVC INTERMITTENT URETHRAL CATHETERS

AMSure® PVC Intermittent Urethral Catheters provide a smooth surface firm enough for insertion yet soft enough to be slightly pliable. The urethral catheters are designed for short-term use and available in various French sizes as well as female, pediatric and male lengths.

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Clear shaft provides easy viewing as you measure
- PVC material eliminates latex allergy concerns
- Individually packaged, sterile, latex-free and BPA-free



ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS860608	8 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860610	10 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860612	12 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860614	14 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860616	16 Fr, 6 in. (15.24 cm)	Female	50/Case
AS861006	6 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861008	8 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861010	10 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861608	8 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861610	10 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861612	12 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861614	14 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861616	16 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861618	18 Fr, 16 in. (40.64 cm)	Male	50/Case

PVC INTERMITTENT URETHRAL CATHETER WITH R[®] POLISHED[®] EYES

AMSSure[®] PVC Intermittent Catheters with R-Polished[®] technology are designed for short-term use. R-Polished[®] technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort. The intermittent catheters are available in various French sizes as well as female, pediatric and male lengths.

- Advanced R-Polished[®] eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- PVC material eliminates latex allergy concerns
- Sterile, latex-free and BPA-free

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS960610	10 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960612	12 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960614	14 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960616	16 Fr, 6 in. (15.24 cm)	Female	50/Case
AS961010	10 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS961610	10 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961612	12 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961614	14 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961616	16 Fr, 16 in. (40.64 cm)	Male	50/Case



URETHRAL INSERTION TRAYS AND KITS

AMSSure[®] Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization.

- Trays provide clinician flexibility to choose their preferred drainage or leg bag
- Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection

PVC CATHETER

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, latex-free and BPA-free

ITEM CODE	URETHRAL CATHETER TYPE	QUANTITY/UNIT
AS870	14 Fr PVC (Latex-free)	20/Case
AS87114	14 Fr Red Rubber Latex	20/Case
AS87118	18 Fr Red Rubber Latex	20/Case

RED RUBBER CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free and BPA-free

Trays include:

- 1000 mL outer basin tray
- 16 in. urethral catheter
- Waterproof and fenestrated drapes
- Vinyl powder-free gloves
- (3) 10% povidone iodine swab sticks
- Lubricating jelly
- 4 oz. specimen container



PRE-CONNECTED DRAINAGE BAG URETHRAL CATHETER KIT

AMSSure® Pre-connected Drainage Bag Closed Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the low-profile drainage bag which helps to minimize contamination.

- Kits provide clinician flexibility to choose their preferred drainage or leg bag
- Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection
- Sterile

PVC CATHETER

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, latex-free and BPA-free

RED RUBBER CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free and BPA-free

Kits include:

- 1200 mL low-profile urinary drainage bag
- 16 in. (40.64 cm) urethral catheter
- CSR wrap, waterproof and fenestrated drapes
- (2) Nitrile powder-free gloves
- (3) Povidone-iodine swab sticks
- Lubricating jelly
- 2 oz. specimen container with ID label



AS89214R

ITEM CODE	14 FR RED RUBBER CATHETER CONTAINS LATEX	14 FR PVC CATHETER LATEX-FREE	QUANTITY/UNIT
AS89214R	●		20/Case
AS89214		●	20/Case

SELF-CATHETERIZATION KITS

URETHRAL SELF-CATHETERIZATION KIT WITH R[®] POLISHED[®] EYES

AMSure[®] Self Catheterization kits with R-Polished[®] technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polished[®] technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization.

- Advanced R-Polished[®] eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, latex-free and BPA-free

Kits include:

- Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter
- 1100 mL collection bag with insertion tip and protector cap
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Waste disposal bag
- Gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)



ITEM CODE	FRENCH SIZE	QUANTITY/UNIT
AS85008	8 Fr	50/Case
AS85010	10 Fr	50/Case
AS85012	12 Fr	50/Case
AS85014	14 Fr	50/Case
AS85016	16 Fr	50/Case

IRRIGATION TRAYS AND KITS

AMSSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. Latex-free, DEHP-free and BPA-free, the trays are available sterile or non-sterile to suit your clinical need.

BULB SYRINGE

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- Latex-free, DEHP-free and BPA-free
- Includes syringe tip protector



THUMB CONTROL RING SYRINGE

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- Latex-free, DEHP-free and BPA-free
- Includes syringe tip adapter and protector cap



Non-sterile irrigation kits include:

- 500 mL graduated container
- Patient ID label and tip protector

ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QUANTITY/UNIT
AS121		●		30/Case
AS126	●			30/Case
AS127	●		●	30/Case

Sterile irrigation trays include:

- 1000 mL outer tray
- 500 mL graduated container
- Waterproof underpad
- Alcohol prep pad
- Patient ID label and tip protector



ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QUANTITY/UNIT
AS130		●	20/Case
AS136	●		20/Case

SYRINGES

- Latex-free, DEHP-free and BPA-free

BULB SYRINGE


- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Tip protector included

THUMB CONTROL RING SYRINGE

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Includes syringe luer tip adapter and protector cap

FLAT TOP SYRINGE

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- Latex-free, DEHP-free and BPA-free
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



ITEM CODE	60 mL SYRINGE TYPE			PACKAGE TYPE				STERILE	QUANTITY/UNIT
	Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Pack	Resealable Pole Bag	Bulk		
AS116			●			●			30/Case
AS016		●				●			30/Case
AS010	●						●		50/Case
AS011	●			●				●	50/Case
AS011P	●				●			●	50/Case
AS015		●		●				●	50/Case
AS115			●	●				●	50/Case

SPECIMEN CONTAINERS

AMSSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- Latex-free, DEHP-free and BPA-free

ITEM CODE	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	QUANTITY/ UNIT
AS340		●		●		●	100/Case
AS341				●	●		100/Case
AS343				●			500/Case
AS345	●			●*			100/Case

*100% sterile



AS340

AMSure® RESPIRATORY THERAPY SYSTEMS

ENHANCING O₂ DELIVERY

AMSure® Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

• Respiratory Masks	31
• Suction Catheters, Trays and Kits	33
• Nebulizer Kits	34
• Aerosol Drainage Bag	35
• Tracheostomy Clean and Care Trays	35
• Suction Tubing	36
• Suction Yankauers and Kits	37
• Suction Device Holster and Kits	38
• Accessories	39



RESPIRATORY MASKS

One Mask™

One Mask™ is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide FiO₂ ranges from 28-99% by simply changing the flow rate of oxygen. One Mask™ is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO₂
- Access for hydration, oral suction or oral medicine
- Efficient FiO₂ escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Non-sterile and latex-free




ITEM CODE	DESCRIPTION	QUANTITY/UNIT
OM001	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

NASAL OXYGEN CANNULA

AMSSure® Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, latex-free and BPA-free



ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS75080	Nasal Cannula, 7 ft. Star-shaped Lumen Oxygen Tubing	Non-Flared Prongs, Adult	50/Case
AS75085	Nasal Cannula, 7 ft. Star-shaped Lumen Oxygen Tubing	Flared Prongs, Adult	50/Case

OXYGEN MASK

AMSSure® Oxygen Masks provide medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, latex-free and BPA-free



ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS74010	Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	Adult	50/Case

NON-REBREATHER OXYGEN MASK

AMSure® Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, latex-free and BPA-free



AS75020

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS75010	Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	Adult	50/Case
AS75020	Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	Pediatric	50/Case

AEROSOL MASK

AMSure® Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, latex-free and BPA-free



AS75050

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS75050	Aerosol Mask	Pediatric	50/Case

TRACHEOSTOMY MASK

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, latex-free and BPA-free



AS75060

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS75060	Adult Tracheostomy Mask	50/Case

SUCTION CATHETERS, TRAYS AND KITS

SUCTION CATHETERS

AMSSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- Opposing staggered eyelets reduce the potential for blocked suction
- Whistle cap port enhances controlled suction
- Sterile, latex-free and BPA-free



ITEM	DESCRIPTION	FEATURES	QTY/UNIT
AS361	6 Fr, 17 ½ in.	Straight Pack	50/Case
AS361C	6 Fr, 17 ½ in.	Coiled Pack	50/Case
AS362	8 Fr, 17 ½ in.	Straight Pack	50/Case
AS362C	8 Fr, 17 ½ in.	Coiled Pack	50/Case
AS363	10 Fr, 23 in.	Straight Pack	50/Case
AS363C	10 Fr, 23 in.	Coiled Pack	50/Case
AS364	12 Fr, 23 in.	Straight Pack	50/Case
AS364C	12 Fr, 23 in.	Coiled Pack	50/Case
AS365	14 Fr, 23 in.	Straight Pack	50/Case
AS365C	14 Fr, 23 in.	Coiled Pack	50/Case
AS366	16 Fr, 23 in.	Straight Pack	50/Case
AS366C	16 Fr, 23 in.	Coiled Pack	50/Case
AS367	18 Fr, 23 in.	Straight Pack	50/Case
AS367C	18 Fr, 23 in.	Coiled Pack	50/Case

MINI SUCTION CATHETER TRAYS

The AMSSure® Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- Opposing staggered eyelets reduce the potential for blocked suction
- DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients
- Whistle tip enhances controlled suction
- Sterile, latex-free and BPA-free



ITEM CODE	DESCRIPTION	CATHETER SIZE	FEATURES	QUANTITY/UNIT
SCT08	Tray, Vinyl Gloves and 17 ½ in. Suction Catheter	8 Fr	DeLee Tip	100/Case
SCT10	Tray, Vinyl Gloves and 23 in. Suction Catheter	10 Fr	DeLee Tip	100/Case
SCT12	Tray, Vinyl Gloves and 23 in. Suction Catheter	12 Fr	Whistle Tip	100/Case
SCT14	Tray, Vinyl Gloves and 23 in. Suction Catheter	14 Fr	Whistle Tip	100/Case
SCT16	Tray, Vinyl Gloves and 23 in. Suction Catheter	16 Fr	Whistle Tip	100/Case

SUCTION CATHETER KITS

AMSure® Suction Catheter Kits provide hassle-free and convenient extraction of mucus or saliva from the upper respiratory airway.

- Uniquely designed two-eyed whistle tip offers safe and effective suctioning
- Sterile, latex-free and BPA-free



AS371

ITEM	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
AS371	17 1/2 in. Whistle Tip Catheter	6 Fr	1 Vinyl Glove	50/Case
AS381	17 1/2 in. Whistle Tip Catheter	6 Fr	2 Vinyl Gloves	50/Case
AS372	17 1/2 in. Whistle Tip Catheter	8 Fr	1 Vinyl Glove	50/Case
AS382	17 1/2 in. Whistle Tip Catheter	8 Fr	2 Vinyl Gloves	50/Case
AS373	23 in. Whistle Tip Catheter	10 Fr	1 Vinyl Glove	50/Case
AS383	23 in. Whistle Tip Catheter	10 Fr	2 Vinyl Gloves	50/Case
AS374	23 in. Whistle Tip Catheter	12 Fr	1 Vinyl Glove	50/Case
AS384	23 in. Whistle Tip Catheter	12 Fr	2 Vinyl Gloves	50/Case
AS375	23 in. Whistle Tip Catheter	14 Fr	1 Vinyl Glove	50/Case
AS385	23 in. Whistle Tip Catheter	14 Fr	2 Vinyl Gloves	50/Case
AS376	23 in. Whistle Tip Catheter	16 Fr	1 Vinyl Glove	50/Case
AS386	23 in. Whistle Tip Catheter	16 Fr	2 Vinyl Gloves	50/Case
AS387	23 in. Whistle Tip Catheter	18 Fr	2 Vinyl Gloves	50/Case

NEBULIZER KITS

AMSure® Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or an aerosol mask.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, latex-free and BPA-free



AS78010



AS78020

T-MOUTHPIECE

- Side-entry tubing design allows nebulizer to stand on its own while adding medications

AEROSOL MASK

- Soft and comfortable elastic strap reduces facial pressure

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS78010	Nebulizer, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78020	Nebulizer, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol Mask	50/Case

AEROSOL DRAINAGE BAG

AMSSure® Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- Pre-connected self-locking plastic chain secures the system to a stable object
- Convenient drain port allows for easy one-handed drainage
- Non-sterile and latex-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS75035	750 mL Aerosol Drainage Bag	50/Case



AS75035

TRACHEOSTOMY CLEAN AND CARE TRAYS

AMSSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- Sterile, latex-free and BPA-free

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS860	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
AS861	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	20/Case



AS861

SUCTION TUBING

SUCTION TUBING WITH CONNECTORS

AMSSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Collapse-resistant striated tubing improves clinical outcomes
- Molded connectors deliver secure connections to reduce procedure interruptions
- Latex-free, BPA-free, sterile and non-sterile options



ITEM CODE	DESCRIPTION	CONNECTOR COLOR	FEATURES	STERILE	QUANTITY/UNIT
AS820	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green		●	50/Case
AS820B	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green			100/Case
AS821	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green		●	50/Case
AS821B	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green			50/Case
AS822	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green		●	50/Case
AS822B	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green			50/Case
AS823	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green		●	20/Case
AS823B	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green			50/Case
AS824	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green		●	20/Case
AS824B	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green			50/Case
AS825	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue		●	50/Case
AS825B	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue			50/Case
AS826	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue		●	50/Case
AS826B	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue			50/Case
AS826C	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue	Pinch Clamp	●	50/Case
AS827	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue		●	20/Case
AS827B	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue			50/Case
AS828	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue		●	20/Case
AS828B	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue			50/Case
AS829B-01	1/4 in. ID x 100 ft. LN (6 mm x 30.0 m)	Light Blue	Connectors every 6 ft. (1.8 m)		1/Case
AS825BM*	5/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue			50/Case
AS825M	5/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue		●	50/Case
AS826BM	5/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue			50/Case
AS826M	5/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue		●	50/Case
AS827BM	5/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue			50/Case
AS827M*	5/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue		●	50/Case
AS828BM	5/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue			50/Case
AS828M	5/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue		●	50/Case

*Made to order

SUCTION TUBING WITHOUT CONNECTORS

- PVC, non-sterile, disposable, latex-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
43457-25	1/32 in. ID x 24 in. (7 mm x 0.61 m)	25/Case
43455-01	1/32 in. ID x 100 ft. (7 mm x 30 m)	1/Case

43457-25



SUCTION YANKAUERS AND KITS

AMSSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, latex-free, DEHP-free and BPA-free



ITEM CODE	HANDLE		TIP		FEATURES		KIT WITH TUBING	QUANTITY/UNIT
	Rigid	Flexible	Regular	Bulbous	Vented	Non-Vented		
AS830	●			●				50/Case
AS831	●			●				50/Case
AS832	●				●			50/Case
AS833	●				●			50/Case
AS834		●	●					50/Case
AS835		●	●		●			50/Case
AS836		●		●				50/Case
AS837		●		●	●			50/Case
AS825P	●			●		●	By 6 ft. (1.8 m)	50/Case
AS827P*	●			●		●	By 12 ft. (3.66 m)	20/Case

*Made to order

SUCTION DEVICE HOLSTER AND KITS

STERICADI™

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- All kits includes SteriCadi™ are pre-packaged with a sterile yankauer
- Sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
SIH0001	SteriCadi™ Suction Device Holder with Clip	30/Case

STERICADI™ KITS

ITEM CODE	STERICADI™ HOLDER	YANKAUER	TIP			FEATURES	QUANTITY/UNIT
			Rigid	Regular Beveled	Bulbous		
SIH0002	●	●		●			30/Case
SIH0003	●	●		●		●	30/Case
SIH0004	●	●			●		30/Case
SIH0005	●	●			●	●	30/Case



ACCESSORIES

OXYGEN TUBING

7 ft. Star-shaped Lumen Tubing with Standard 22 mm Connector

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, latex-free and BPA-free

ITEM	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS76007	7 ft. Oxygen Tubing	Star-shaped Lumen	50/Case



RESPIRATORY ADAPTERS

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- Non-sterile, latex-free and BPA-free

ITEM	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS730100	Universal Cuff Adapter	22 mm ID Connector	50/Case
AS730200	Multi Adapter	15 mm ID / 22 mm OD Connector	50/Case



AS730100 AS730200

AMSure® SURGICAL AND NURSING SOLUTIONS

EVERYDAY PATIENT CARE SOLUTIONS

AMSure® Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

• Suction Tubing	41
• Suction Yankauers and Kits	42
• Suction Device Holster and Kits	43
• Ear/Ulcer Syringes	43
• Vaginal Speculums	44
• Disposable Umbilical Cord Clamp	44
• Specimen Containers	44
• Cleansing Enema Sets	45
• 5-in-1 Straight Connector	45



SUCTION TUBING

AMSURE® SUCTION TUBING WITH CONNECTORS

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Collapse-resistant striated tubing improves clinical outcomes
- Molded connectors deliver secure connections to reduce procedure interruptions
- Latex-free, BPA-free, sterile and non-sterile options



ITEM CODE	DESCRIPTION	CONNECTOR COLOR	FEATURES	STERILE	QUANTITY/UNIT
AS820	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green		●	50/Case
AS820B	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green			100/Case
AS821	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green		●	50/Case
AS821B	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green			50/Case
AS822	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green		●	50/Case
AS822B	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green			50/Case
AS823	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green		●	20/Case
AS823B	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green			50/Case
AS824	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green		●	20/Case
AS824B	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green			50/Case
AS825	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue		●	50/Case
AS825B	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue			50/Case
AS826	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue		●	50/Case
AS826B	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue			50/Case
AS826C	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue	Pinch Clamp	●	50/Case
AS827	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue		●	20/Case
AS827B	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue			50/Case
AS828	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue		●	20/Case
AS828B	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue			50/Case
AS829B-01	1/4 in. ID x 100 ft. LN (6 mm x 30.0 m)	Light Blue	Connectors every 6 ft. (1.8 m)		50/Case
AS825BM*	9/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue			50/Case
AS825M	9/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue		●	50/Case
AS826BM	9/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue			50/Case
AS826M	9/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue		●	50/Case
AS827BM	9/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue			50/Case
AS827M*	9/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue		●	50/Case
AS828BM	9/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue			50/Case
AS828M	9/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue		●	50/Case

*Made to order

SUCTION YANKAUER AND KITS

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, latex-free, DEHP-free and BPA-free



ITEM CODE	HANDLE		TIP		FEATURES		KIT WITH TUBING	QUANTITY/UNIT
	Rigid	Flexible	Regular	Bulbous	Vented	Non-Vented	1/4 in. ID (6 mm)	
AS830	●		●			●		50/case
AS831	●		●		●			50/case
AS832	●			●		●		50/case
AS833	●			●	●			50/case
AS834		●	●			●		50/case
AS835		●	●		●			50/case
AS836		●		●		●		50/case
AS837		●		●	●			50/case
AS825P	●			●		●	By 6 ft. (1.8 m)	50/case
AS827P*	●			●		●	By 12 ft. (3.66 m)	20/case

*Made to order

SUCTION DEVICE HOLSTER AND KITS

STERICADI™

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- All kits include SteriCadi™ and are pre-packaged with a sterile yankauer
- Sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
SIH0001	SteriCadi™ Suction Device Holder with Clip	30/Case

STERICADI™ KITS

ITEM CODE	STERICADI™ HOLDER	YANKAUER	TIP		FEATURES	QUANTITY/UNIT
			Rigid	Regular Beveled		
SIH0002	●	●		●		30/Case
SIH0003	●	●		●	●	30/Case
SIH0004	●	●			●	30/Case
SIH0005	●	●			●	30/Case



EAR/ULCER SYRINGES

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options and latex-free

ITEM	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS00501	1 oz. syringe	Non-sterile, Bulk	300/Case
AS00502	2 oz. syringe	Non-sterile, Bulk	200/Case
AS00503	3 oz. syringe	Non-sterile, Bulk	150/Case
AS00501S	1 oz. syringe	Sterile, Single Pack	50/Case
AS00502S	2 oz. syringe	Sterile, Single Pack	50/Case
AS00503S	3 oz. syringe	Sterile, Single Pack	50/Case



VAGINAL SPECULUMS

AMSure® Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- Round-edge design helps to increase patient comfort
- Single-use disposable
- Non-sterile and latex-free



AS032M

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
AS032S	Small Vaginal Speculum	100/Case
AS032M	Medium Vaginal Speculum	100/Case
AS032L	Large Vaginal Speculum	100/Case

DISPOSABLE UMBILICAL CORD CLAMP

AMSure® Umbilical Cord Clamp with safety lock helps provide a secure connection.

- Safety lock is designed to decrease risk of accidental disconnection
- Grooved teeth help provide an effective grip
- Sterile, latex-free and BPA-free



UCC100

SPECIMEN CONTAINERS

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- Latex-free, DEHP-free and BPA-free

ITEM CODE	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABEL	QUANTITY/UNIT
AS340		●		●		●	100/Case
AS341			●	●	●		100/Case
AS343			●				500/Case
AS345	●			●*			100/Case

*100% sterile



AS340

CLEANSING ENEMA SETS

AMSure® Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

BAG SET

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, latex-free and BPA-free



BUCKET SET

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, latex-free and BPA-free



ITEM	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS330	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case

5-IN-1 STRAIGHT CONNECTOR

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 7/16 in.) sizes of tubing.

- Latex-free, DEHP-free and BPA-free

ITEM	DESCRIPTION	FEATURES	QUANTITY/UNIT
AS807	5-in-1 Connector	Sterile	100/Case
AS807B	5-in-1 Connector	Non-sterile, Bulk	500/Case



RECEPTAL® Waste Fluid Management Systems

SIMPLE, SAFE AND COMPLETE SOLUTIONS

Amsino is a leading global provider of waste fluid management systems. The breadth of Amsino's RECEPTAL® waste fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

• RECEPTAL® Canisters	47
• RECEPTAL® Liners	47
• EZE-VAC® Suction Canister System	48
• ACCU-MEASURE® Critical Measurement Devices	48
• Suction Tubing	49
• Suction Yankauers and Kits	50
• Suction Device Holster and Kits	51
• Accessories	52
• Reusable Hardware.....	54



RECEPTAL® CANISTERS

RECEPTAL® Canisters are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL® Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, latex-free and DEHP-free



ITEM CODE	DESCRIPTION	QUANTITY/UNIT
43449-05	1000 mL Canister	10/Case
43423-01	1500 mL Canister	10/Case
43445-01	2000 mL Canister	10/Case
43503-01	3000 mL Canister	8/Case

RECEPTAL® LINERS

RECEPTAL® liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL® disposable liners are to be used with the same size RECEPTAL® Reusable Canisters.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, latex-free and BPA-free



ITEM CODE	VOLUME				POUR SPOUT	ORTHOPEDIC PORT	FEATURES	PRE-GELLED SOLIDIFIER	QUANTITY/UNIT
	1000 mL	1500 mL	2000 mL	3000 mL					
43010-05	●				●		●		50/Case
43025-01		●			●				50/Case
43044-01		●			●		●		50/Case
43027-01			●		●				50/Case
43041-01			●			●			50/Case
43042-01			●				●		50/Case
43044-05			●		●		●		50/Case
OL655-01				●	●	●			50/Case
OL656-01				●	●	●	●		50/Case
43010S	●						●	●	50/Case
43015S		●					●	●	50/Case
43020S			●				●	●	50/Case
43030S				●			●	●	50/Case

EZE-VAC® SUCTION CANISTER SYSTEM

EZE-VAC® suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD®, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
43207-01	700 mL Canister with ACCU-MEASURE® Critical Measurements, Lid and 2 Elbows	32/Case
43203-01	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01	3300 mL Canister, Lid and 2 Elbows	12/Case



ACCU-MEASURE® CRITICAL MEASUREMENT

ACCU-MEASURE® critical measurement disposable canisters enable the user to make contact-free critical measurement of smaller quantities of aspirate.

- Canisters are non-sterile, latex-free, DEHP-free and BPA-free

RECEPTAL® ACCU-MEASURE®

RECEPTAL® ACCU-MEASURE® is a 250 mL disposable critical measurement canister with graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL. The canister can be conveniently attached to a RECEPTAL® canister or used independently.



ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
43306-01	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL® Canisters	10/Case

EZE-VAC® ACC-MEASURE®

EZE-VAC® ACC-MEASURE® is a 700 mL disposable critical measurement canister with graduations in 5 mL measure increments from 0 to 100 mL, 10 mL increments measure from 100 mL to 300 mL and 25 mL increments measure from 300 mL to 700 mL.



ITEM CODE	DESCRIPTION	QUANTITY/UNIT
43207-01	700 mL Critical Measurement Canister, Lid and 2 Elbows	32/Case

SUCTION TUBING

SUCTION TUBING WITH CONNECTORS

AMSSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Collapse-resistant striated tubing improves clinical outcomes
- Molded connectors deliver secure connections to reduce procedure interruptions
- Latex-free, BPA-free and sterile and non-sterile options



ITEM CODE	DESCRIPTION	CONNECTOR COLOR	FEATURES	STERILE	QUANTITY/UNIT
AS820	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green		●	50/Case
AS820B	3/16 in. ID x 18 in. LN (5 mm x 46 cm)	Green			100/Case
AS821	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green		●	50/Case
AS821B	3/16 in. ID x 6 ft. LN (5 mm x 1.8 m)	Green			50/Case
AS822	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green		●	50/Case
AS822B	3/16 in. ID x 10 ft. LN (5 mm x 3.0 m)	Green			50/Case
AS823	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green		●	20/Case
AS823B	3/16 in. ID x 12 ft. LN (5 mm x 3.7 m)	Green			50/Case
AS824	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green		●	20/Case
AS824B	3/16 in. ID x 20 ft. LN (5 mm x 6.1 m)	Green			50/Case
AS825	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue		●	50/Case
AS825B	1/4 in. ID x 6 ft. LN (6 mm x 1.8 m)	Light Blue			50/Case
AS826	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue		●	50/Case
AS826B	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue			50/Case
AS826C	1/4 in. ID x 10 ft. LN (6 mm x 3.0 m)	Light Blue	Pinch Clamp	●	50/Case
AS827	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue		●	20/Case
AS827B	1/4 in. ID x 12 ft. LN (6 mm x 3.7 m)	Light Blue			50/Case
AS828	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue		●	20/Case
AS828B	1/4 in. ID x 20 ft. LN (6 mm x 6.1 m)	Light Blue			50/Case
AS829B-01	1/4 in. ID x 100 ft. LN (6 mm x 30.0 m)	Light Blue	Connectors every 6 ft. (1.8 m)		1/Case
AS825BM*	5/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue			50/Case
AS825M	5/32 in. ID x 6 ft. LN (7 mm x 1.8 m)	Medium Blue		●	50/Case
AS826BM	5/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue			50/Case
AS826M	5/32 in. ID x 10 ft. LN (7 mm x 3.0 m)	Medium Blue		●	50/Case
AS827BM	5/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue			50/Case
AS827M*	5/32 in. ID x 12 ft. LN (7 mm x 3.7 m)	Medium Blue		●	50/Case
AS828BM	5/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue			50/Case
AS828M	5/32 in. ID x 20 ft. LN (7 mm x 6.1 m)	Medium Blue		●	50/Case

*Made to order

SUCTION YANKAUERS AND KITS

AMSSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, latex-free, DEHP-free and BPA-free



ITEM CODE	HANDLE		TIP		FEATURES		KIT WITH TUBING	QUANTITY/UNIT
	Rigid	Flexible	Regular	Bulbous	Vented	Non-Vented	1/4 in. ID (6 mm)	
AS830	●		●			●		50/Case
AS831	●		●		●			50/Case
AS832	●			●		●		50/Case
AS833	●			●	●			50/Case
AS834		●	●			●		50/Case
AS835		●	●		●			50/Case
AS836		●		●		●		50/Case
AS837		●		●	●			50/Case
AS825P	●			●		●	By 6 ft. (1.8 m)	50/Case
AS827P*	●			●		●	By 12 ft. (3.66 m)	20/Case

*Made to order

SUCTION DEVICE HOLSTER AND KITS

STERICADI™

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- All kits includes SteriCadi™ are pre-packaged with a sterile yankauer
- Sterile, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
SIH0001	SteriCadi™ Suction Device Holder with Clip	30/Case

STERICADI™ KITS

ITEM CODE	STERICADI™ HOLDER	YANKAUER	TIP		FEATURES	QUANTITY/UNIT
			Rigid	Regular Beveled		
SIH0002	●	●		●		30/Case
SIH0003	●	●		●	●	30/Case
SIH0004	●	●			●	30/Case
SIH0005	●	●			●	30/Case



SIH0002

ACCESSORIES

RECEPTAL® SPECIMEN TRAPS

RECEPTAL® Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL® and EZE-VAC® lids with a pour spout.



43303-01

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	FEATURES		QUANTITY/UNIT
		Suction Port Size	Stiffener	
43303-01	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	⅜ in. (0.95 cm)		24/Case
43303-02	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	⅜ in. (0.95 cm)	●	24/Case
43303-05	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	½ in. (1.27 cm)	●	24/Case

EXTERNAL VAC-GARD® SHUT-OFF VALVE

VAC-GARD® external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, latex-free, DEHP-free and BPA-free

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
43309-01	External Shut-Off Valve	¼ in. (0.64 cm) Male x Female Elbow	10/Case



43309-01

EZE-VAC® TANDEM TUBING

EZE-VAC® Tandem disposable tubing allows two or more EZE-VAC® canisters to work together to meet procedural requirements.

- Non-sterile, latex-free and BPA-free



43317-25

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
43317-25	18 in. (0.45 m) Tandem Tubing	Use with EZE-VAC® Canisters	25/Case

RECEPTAL® ELBOW AND TAPER ADAPTERS

The elbows and taper adapter attach to the suction port to provide more flexibility.

ITEM CODE	DESCRIPTION	FEATURES		QUANTITY/UNIT
		1/4 in. (0.64 cm)	3/8 in. (0.95 cm)	
43494-25	1/4 in. (0.64 cm) Male x Female Elbow	●		25/Case
43316-01	1/4 in. (0.64 cm) Male x Female Elbow	●		100/Case
43473-25	1/4 in. (0.64 cm) Male Taper		●	25/Case



RECEPTAL® CANISTER VACUUM "T"

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
44300-05	1/4 in. (0.64 cm) Male Tapers, O-ring and Screw	5/Case



RECEPTAL® VACUUM TAPER

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
43448-25	1/4 in. (0.64 cm) Male x Female Taper	25/Case



RECEPTAL® ACCU-MEASURE® DEVICE BRACKET

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL® Reusable Canister.

ITEM CODE	DESCRIPTION	FEATURES	QUANTITY/UNIT
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL® Canisters	10/Case



REUSABLE HARDWARE

FLOOR STANDS AND TROLLEYS

Reusable Floor Stand and Trolley configurations provide procedural flexibility.

ITEM CODE	DESCRIPTION	FEATURES		CANISTER SIZE COMPATIBILITY				QUANTITY/UNIT
		Regulator Mounting Plate	Stainless Steel with Conductive Casters	1500 mL	2000 mL	3000 mL	All sizes	
43476-P2	Adjustable Height, Built-in Regulator Manifold Trolley w/(4) Polycarbonate Canister Supports			●	●			1/Case
43476-P3	Adjustable Height, Built-in Regulator Manifold Trolley w/(4) Polycarbonate Canister Supports					●		1/Case
43475-01	Adjustable Height, Small Head Trolley, Holds Up to 8 RECEPTAL® Canisters or 4 EZE-VAC® Canisters		●	●	●	●	●	1/Case
43474-01	Adjustable Height, Large Head Trolley, holds Up to 12 RECEPTAL® or 8 EZE-VAC® Canisters	●	●	●	●	●	●	1/Case
43447-01	In-Line Quad Floor with Rollers	●	●	●	●			1/Case
43432-01	Tandem Floor Stand with Rollers		●	●	●			1/Case
43431-01	Single Canister Floor Stand			●	●			1/Case



43476-P2



43474-01



43475-01



43447-01



43432-01



43431-01

CANISTER SUPPORTS

Reusable Canister Supports are available in single or tandem configurations to accommodate RECEPTAL® and EZE-VAC® suction canisters.

ITEM CODE	SUPPORT TYPE		CANISTER SIZE COMPATIBILITY					MATERIAL	REGULATOR MOUNTING PLATE	QUANTITY/UNIT
			RECEPTAL®			EZE-VAC®				
	Single	Tandem	1000 mL	1500 mL	2000 mL	3000 mL	All Sizes	Polycarbonate	Metal	
43425-05	●			●					●	5/Case
43426-05	●			●	●				●	5/Case
43426P	●			●	●			●		20/Case
43436-05	●					●			●	5/Case
43436P	●					●		●		10/Case
43428-05	●			●	●		●	●	●	5/Case
43429-01		●		●	●			●		1/Case
43439-01		●				●		●		1/Case
43469-12	●						●	●	●	12/Case

Canister Support Wall Brackets are sold separately.



WALL BRACKET

Wall channel brackets are used for mounting canisters.

ITEM CODE	FEATURE	QUANTITY/UNIT
43444-05	Stainless Steel	5/Case
43444P	Polycarbonate	50/Case



BRACKET CLAMP ADAPTER

A Bracket Clamp Adaptor is used with wall bracket (43444-05) for attachment to a vertical pole.



ITEM CODE	DESCRIPTION	QUANTITY/UNIT
44303-10	Stainless Steel	10/Case

CANISTER SUPPORT BRACKET

This is used to attach Canister directly to a vacuum regulator.
Can be used with 1000 mL, 1200 mL and 1500 mL Canisters.

ITEM CODE	DESCRIPTION	QUANTITY/UNIT
44304-05	Stainless Steel	5/Case



44304-05

REGULATOR MOUNTS

Amsino Regulator Mounts facilitate suction tubing connection to accommodate various clinical settings and requirements.

ITEM CODE	DESCRIPTION	MOUNT TYPE		ATTACHMENT		QUANTITY/UNIT
		Standard	DISS	Regulator	Hardware Mount	
43408-05	½ in. (0.32 cm) Male Pipe Thread x ¼ in. (0.64 cm) Male Taper with 1⅛ in. (2.86 cm) Threaded Shank and 2 Locknuts	●		●	●	5/Case
43422-05	½ in. (0.32 cm) Male Pipe Thread x ¼ in. (0.64 cm) Male Taper	●		●		5/Case
43401-05	Male DISS x ¼ in. (0.64 cm) Taper with 1¼ in. (3.18 cm) Shank		●		●	5/Case
43400-05	Male DISS x ¼ in. (0.64 cm) Male Taper		●	●		5/Case



43408-05



43401-05



43422-05



43444-05

INDEX OF ITEM CODES

ITEM CODE	PAGE						
000606.....	13	43408-05.....	56	AS00501	43	AS362	33
003201.....	14	43422-05.....	56	AS00501S	43	AS362C	33
003202.....	14	43423-01	47	AS00502	43	AS363	33
0L655-01	47	43425-05	55	AS00502S	43	AS363C	33
0L656-01	47	43426-05	55	AS00503	43	AS364	33
1010001B.....	8	43426P	55	AS00503S	43	AS364C	33
107201.....	6	43428-05	54	AS010.....	3, 28	AS365	33
108301.....	6	43429-01	55	AS011	3, 28	AS365C	33
108302.....	8	43431-01	54	AS011P	3, 28	AS366	33
108303.....	6	43432-01	54	AS015.....	3, 28	AS366C	33
108305.....	7	43436-05	55	AS016.....	3, 28	AS367	33
108306.....	6	43436P	55	AS032L	44	AS367C	33
108308.....	7	43439-01	55	AS032M	44	AS371	34
109602.....	6	43444-05.....	55	AS032S	44	AS372	34
1510501	9	43444P	55	AS115.....	3, 28	AS373	34
158305.....	9	43445-01	47	AS116.....	3, 28	AS374	34
158306.....	9	43447-01	54	AS121	4, 27	AS375	34
158307.....	9	43448-25.....	53	AS126.....	4, 27	AS376	34
158309.....	9	43449-05	47	AS127	4, 27	AS381	34
158310.....	8	43455-01	37	AS130	4, 27	AS382	34
159602.....	9	43457-25.....	37	AS136	4, 27	AS383	34
43010-05	47	43469-12	55	AS302	18	AS384	34
43010S	47	43473-25.....	53	AS306N	19	AS385	34
43015S	47	43474-01	54	AS307N	19	AS386	34
43020S	47	43475-01	54	AS308N	19	AS387	34
43025-01	47	43476-P2	54	AS309N	19	AS409	19
43027-01	47	43476-P3	54	AS3100	10	AS41012	20
43030S	47	43494-25.....	53	AS3101	10	AS41012S	20
43041-01	47	43503-01	47	AS3102	10	AS41014	20
43042-01	47	44300-05.....	53	AS311	19	AS41014S	20
43044-01	47	44303-10.....	55	AS3110A	10	AS41016	20
43044-05	47	44304-05.....	56	AS3110B	10	AS41016S	20
43203-01	48	607201.....	11	AS312	18	AS41018	20
43206-01	48	608305.....	11	AS322	18	AS41018S	20
43207-01	48	608306.....	11	AS32200	18	AS41020	20
43207-01	48	AA3101.....	8	AS326	18	AS41020S	20
43213-01	48	AC3110.....	6	AS32600	18	AS41022	20
43303-01	52	AE0106	14	AS330	45	AS41022S	20
43303-02	52	AE0107	15	AS332	18	AS41024	20
43303-05	52	AE1106	14	AS333	45	AS41024S	20
43305-01	48	AE1112	15	AS340	29, 44	AS41026	20
43305-10	53	AE1138	15	AS341	29, 44	AS41026S	20
43306-01	48	AE3107	14	AS343	29, 44	AS41028	20
43309-01	52	AE3108	14	AS345	29, 44	AS41028S	20
43316-01	53	AFS102	7	AS3500	13	AS41030	20
43317-25	52	AFS107	15	AS3515	13	AS42012	20
43400-05	56	AFS208	7	AS361	33	AS42014	20
43401-05	56	AFS208B.....	7	AS361C	33	AS42016	20

INDEX OF ITEM CODES

ITEM CODE	PAGE	ITEM CODE	PAGE	ITEM CODE	PAGE	ITEM CODE	PAGE
AS42016S	20	AS825B	36, 41, 49	AS870.....	24	MBSS001	16
AS42018	20	AS825BM	36, 41, 49	AS87114.....	24	MR4001	16
AS42018S	20	AS825M	36, 41, 49	AS87118.....	24	MR4107	13
AS42020	20	AS825P	37, 42, 50	AS880.....	21	MR4207	14
AS42020S	20	AS826	36, 41, 49	AS890.....	21	MR4307	15
AS42022	20	AS826B	36, 41, 49	AS89214.....	25	MR4B07	15
AS42022S	20	AS826BM	36, 41, 49	AS89214R	25	MR4P07	16
AS42024	20	AS826C	36, 41, 49	AS89300.....	21	MR4R10	14
AS42024S	20	AS826M	36, 41, 49	AS89314	22	MR4T07	14
AS42026	20	AS827	36, 41, 49	AS89314S	22	MR4T32	14
AS42028	20	AS827B	36, 41, 49	AS89316	22	MR4TF-04	15
AS42030	20	AS827BM	36, 41, 49	AS89316S	22	MRA102-84	8
AS44008	23	AS827M*	36, 41, 49	AS89318	22	MRA10E-104	7
AS44010	23	AS827P	37, 42, 50	AS89318S	22	MRA10E-95	6
AS44012	23	AS828	36, 41, 49	AS89330.....	21	MRA15-100	9
AS44014	23	AS828B	36, 41, 49	AS960610	24	MRA15B-105	10
AS44016	23	AS828BM	36, 41, 49	AS960612	24	MRA60B-92	11
AS44018	23	AS828M	36, 41, 49	AS960614	24	MRS60B-106	12
AS44020	23	AS829B-01 ..	36, 41, 49	AS960616	24	MRS60B-117	12
AS44022	23	AS830	37, 42, 50	AS961010	24	MRSA10-84	6
AS44024	23	AS831	37, 42, 50	AS961012	24	MRSA103-94	8
AS44028	23	AS832	37, 42, 50	AS961014	24	MRSA10E-107	7
AS730100	39	AS833	37, 42, 50	AS961016	24	MRSA10E-110	8
AS730200	39	AS834	37, 42, 50	AY0200	16	MRSA15-100	9
AS74010.....	31	AS835	37, 42, 50	D38301	12	MRSA15-91	9
AS75010.....	32	AS836	37, 42, 50	D38302	12	MRSA15E-120	9
AS75020.....	32	AS837	37, 42, 50	D38305	12	MRSA15E-88	9
AS75035.....	35	AS85008.....	26	E-0500M.....	2	MRSA15E-96	10
AS75050.....	32	AS85010.....	26	E-1004.....	2	MRSA60-100	11
AS75060.....	32	AS85012.....	26	E-1200.....	2	MRSA60-84	11
AS75080.....	31	AS85014.....	26	E-1200M.....	2	MRSA602-100	11
AS75085.....	31	AS85016.....	26	E-1203.....	2	MRSA60E-88	11
AS76007.....	39	AS860.....	35	ENF1004K.....	2	MRSB-94	13
AS78010.....	34	AS860608	23	ENF1200K.....	2	MRSB-102	13
AS78020.....	34	AS860610	23	ENF1203K.....	2	MRSB-105	13
AS807	45	AS860612	23	ENF500K.....	2	MRSFC-19	15
AS807B	45	AS860614	23	ENK126.....	4	OM001	31
AS820	36, 41, 49	AS860616	23	ENK127.....	4	S-1000-S1	2
AS820B	36, 41, 49	AS861.....	35	ENK136.....	4	SCT08	33
AS821	36, 41, 49	AS861006	23	ENS015	3	SCT10	33
AS821B	36, 41, 49	AS861008	23	ENS016	3	SCT12	33
AS822	36, 41, 49	AS861010	23	ENS115	3	SCT14	33
AS822B	36, 41, 49	AS861608	23	ENS116	3	SCT16	33
AS823	36, 41, 49	AS861610	23	LBE306N.....	19	SIH0001	38, 43, 51
AS823B	36, 41, 49	AS861612	23	LBE307N.....	19	SIH0002	38, 43, 51
AS824	36, 41, 49	AS861614	23	LBE308N.....	19	SIH0003	38, 43, 51
AS824B	36, 41, 49	AS861616	23	LBE309N.....	19	SIH0004	38, 43, 51
AS825	36, 41, 49	AS861618	23	MB4V01	16	SIH0005	38, 43, 51
						UCC100	44



Amsino International, Inc.

708 Corporate Center Drive
Pomona, CA 91768

P: 800.632.6746
F: 909.626.3888

info@amsino.com

www.amsino.com