

PRODUCT CATALOG





AMSure® Enteral Feeding Systems..... 1

PUGGLE® Enteral Feeding Pump.....	2
Sentinel™ Enteral Feeding Pump.....	3
ENFit® and Non-ENFit Sets.....	3
ENFit and Non-ENFit Syringes.....	3
ENFit and Non-ENFit Irrigation Kits and Trays.....	5

AMSafe® Medication Delivery Systems 6

Gravity Administration Sets.....	7
10 Drops/mL.....	7
15 Drops/mL.....	9
60 Drops/mL.....	11
AMSafe-3® Selectable Drip Chamber	
10, 15 & 60 drops/mL.....	11
Secondary Piggyback Sets.....	13
IV Extension Sets.....	13
Standard Extension Sets.....	13
Microbore Extension Sets.....	14
Smallbore Extension Sets.....	14
Needle-Free Connectors.....	15
IV Pump Administration Sets.....	16
Sure-Lok® Needle-Free IV Pump Administration Sets.....	16

AMSure Urological Care 17

Foley Catheter Trays & Collection.....	18
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	18
Add-a-Foley Catheter Trays.....	19
Pre-Connected Foley Catheter Urine Meter Kits.....	20
Add-a-Foley Urine Meter Trays.....	21
Foley Insertion Trays.....	21
Urinary Drainage Bags.....	22
Urine Meter.....	23
ACCU-MEASURE® Pediatric Urine Meter.....	23
Foley Catheters.....	24
Hydrophilic Latex Foley Catheters.....	24
Silicone-Coated Latex Foley Catheters.....	25
100% Silicone Foley Catheters.....	26
Intermittent Urinary Catheters.....	27
Urethral Insertion Trays and Kits.....	27
Closed System Urethral Drainage Kits.....	28
Self-Catheterization Kits with R-Polished™ Eye Catheters.....	29
PVC Intermittent Urethral Catheter with R-Polished Eyes.....	30
Red Rubber Latex Urethral Catheter.....	31
Irrigation Sets.....	32
Cystoscopy & TUR Irrigation Sets.....	32
Ancillary & Patient Care Products.....	33
Urinary Leg Bags.....	33
Extension Tube for Leg Bags.....	33
Pediatric Urine Collector.....	33
Irrigation Trays and Kits.....	34
Syringes.....	35
Specimen Containers.....	36

Inhalation, Irrigation, and Injection Solutions..... 37

Sterile Water for Inhalation, USP.....	38
Sterile Water Prefilled Humidifiers, USP.....	38
Nebulizer Kits.....	63
Saline Products.....	39
Prefilled Saline Flush Syringes.....	39
Sodium Chloride Inhalation Solution, (0.9%) USP.....	39
Device and Wound Irrigation.....	40

Waste Fluid Management Systems 41

EZE-VAC® Suction Canister System.....	42
EZE-VAC 2.0 Suction Canister System.....	43
RECEPTAL Liners.....	44
RECEPTAL Twist Liners.....	49
RECEPTAL Outer Shells.....	46
ACCU-MEASURE® Critical Measurement Devices.....	47
RECEPTAL Solidifiers.....	47
AMSure® Suction Tubing with Connectors.....	48
Suction Yankauers and Kits.....	49
SteriCadi™ Suction Device Holster System.....	50
Accessories.....	50
Reusable Hardware.....	54
iRECEPTAL® Surgical Waste Fluid Management System.....	57

AMSure Respiratory Therapy Systems..... 59

Respiratory Masks and Nasal Cannulas.....	60
Nebulizer Kits.....	63
Accessories.....	64
Aerosol Drainage Bag.....	64
Heat and Moisture Exchangers.....	65
Tracheostomy Clean and Care Trays.....	65
Sterile Water and Saline.....	66
Suction Catheters, Trays and Kits.....	67
Suction Tubing with Connectors.....	70
Suction Yankauers and Kits.....	71

AMSure Surgical & Nursing Systems..... 72

Specimen Containers.....	73
Vaginal Speculums.....	73
Disposable Umbilical Cord Clamp.....	73
Ear/Ulcer Syringes.....	74
Cleansing Enema Sets.....	74
Connectors.....	76

Index of REFs 77

Notes 78



Optimizing Nutritional Delivery

AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The systems are designed to provide best patient outcomes and includes:

PUGGLE® Enteral Feeding Pump	2
Sentinel™ Enteral Feeding Pump.....	3
ENFit & Non-ENFit Enteral Feeding Sets.....	3
ENFit & Non-ENFit Enteral Syringes	4
ENFit & Non-ENFit Irrigation Kits & Trays	5

PUGGLE® Enteral Feeding Pump

The versatile pump is designed for pediatric and adult ambulatory applications. Simplified program settings make feeding on-the-go a snap. **PUGGLE is only available the US.**

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours

Dimensions: 5.4 x 5.0 x 2.5 inches (13.7 x 12.7 x 6.4 cm)
Weight: 1.43 lbs (649 g)
Operating Range: 1 to 400 mL/h in 1 mL increments
Delivered Volume: 1 – 9999 mL in 1 mL increments
Accuracy: +/- 5%



Scan the code below to view the PUGGLE instructional video



REF	DESCRIPTION	QTY/UNIT
PG001	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
PE90ST	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit® and Transition Connector	30/Case
PE90B5T	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
PE90B12T	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case



Sentinel™ Enteral Feeding Pump

The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each



ENF1200K

ENFit® and Non-ENFit Enteral Feeding Sets

AMSure Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel™ enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free options, not manufactured with natural rubber latex

REF	ENFit CONNECTOR	STANDARD BORE (DEHP FREE)	LARGE BORE (DEHP FREE)	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
				Spike Only No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/EntraFlo™	
E-1004				•				•	•	•	30/Case
E-1200						•		•	•	•	
E-1203						•	•				
ENF1004K ©	•			•				•	•	•	
ENF1200K ©	•					•		•	•	•	
ENF1203K ©	•					•	•				
ENF1203 KD	•	•				•	•				
ENF1203 KLD	•		•			•	•				
ENF500K ©	•				•			•	•	•	

EntraFlo™ is a registered trademark of Medline Industries, Inc.

Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.

ENFit® is a registered trademark of GEDSA.

ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

Bulb Syringe

- Flexible soft bulb promotes one-handed administration

Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage

Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



AS010



ENS015NS



AS116

ENS115

ENS115T

REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010 ©			•						•		50/Case
AS011 ©			•			•				•	50/Case
AS115 ©		•			•	•				•	50/Case
AS116		•			•			•			30/Case
ENS015 ©	•			•		•				•	50/Case
ENS015NS	•			•				•			30/Case
ENS115 ©	•				•	•				•	50/Case
ENS115NS	•				•			•			30/Case
ENS115T*	•				•			•			30/Case

*Includes Transition Adapter

ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	
AS127		•		•		•	•	
ENK127	•	•		•		•	•	

Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	
ENK136	•	•	•		•	•	•	•	

AMSafe®

MEDICATION DELIVERY SYSTEMS



Safety First


The AMSafe IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:

Gravity Administration IV Sets.....	7
10 Drops/mL.....	7
15 Drops/mL.....	9
60 Drops/mL.....	11
AMSafe-3® Selectable Drip Chamber, 10, 15 & 60 drops/mL.....	11
Secondary Piggyback Sets	13
Standard Extension Sets	13
Microbore Extension Sets.....	14
Smallbore Extension Sets	14
Needle-Free Connectors	15
IV Pump Administration Sets	16
Sure-Lok® Needle-Free IV Pump Administration Sets.....	16

Gravity Administration IV Sets








Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex



108310

10 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108303	10	83 in. (211 cm)	Sure-Lok® Needle-free Y-site	16 mL	50/Case
					
108304	10	89 in. (226 cm)	Non-vented, Roller Clamp, (2) Sure-Lok Needle-free Y-sites with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
109602	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
					
108306	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
108310	10	83 in. (211 cm)	Vented/Non-vented Spike, Roller Clamp, (2) Sure-Lok Needle-free Y-sites	16 mL	50/Case
					
AC3110	10	105 in. (266 cm)	(2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					

10 Drops/mL

General Primary Sets (continued)

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
108305	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case




108308	10	83 in. (211 cm)	(2) Sure-Lok Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
---------------	----	-----------------	--	---------	---------




10 Drops/mL

Flow Controller Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS102	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok Needle-free Y-site	18.5 mL	50/Case



AFS208 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) AMSafe® Needle-free Y-sites	19.2 mL	50/Case
-----------------	----	------------------	--	---------	---------



AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe Needle-free Y-sites, Slide Clamp	19 mL	50/Case
----------------	----	-----------------	---	-------	---------




Gravity Administration IV Sets

10 Drops/mL


General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
1010001B	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case



108302	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
---------------	----	-----------------	--	---------	---------




15 Drops/mL

General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AA3101	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case



158309	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok® Needle-free Y-site	17 mL	50/Case
---------------	----	-----------------	---	-------	---------



158306	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
---------------	----	-----------------	---	-------	---------




159602	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
---------------	----	-----------------	-----------------------	---------	---------



1510501	15	105 in. (266 cm)	Vented/Non-vented, (2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
----------------	----	------------------	---	-------	---------







158305	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
---------------	----	-----------------	---	---------	---------



15 Drops/mL

Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line






REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AS3100	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
AS3101	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
AS3102	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
AS3110A	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					



Gravity Administration IV Sets

60 Drops/mL

General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608304	60	89 in. (226 cm)	Non-vented, Roller Clamp, (2) Sure-Lok® Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
608310 ©	60	83 in. (211 cm)	Vented/Non-vented Spike, Roller Clamp, (2) Sure-Lok® Needle-free Y-sites	16 mL	50/Case
					



AMSafe-3® Selectable Drip Chamber

10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
D38301 ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site	27.7 mL	50/Case
					
D38302 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.2 mL	50/Case
					



Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AS3515 ©	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
AS3500 ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case









IV Extension Sets

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.







AE3108

Standard Extension Sets



REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
000606	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok® Needle-free Y-site, Slide Clamp	2.5 mL	50/Case
				
AE3107 ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
				
AE3108 ©	8 in. (20 cm)	Removable Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
				
AE3108B	8 in. (20 cm)	Bonded Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
				
003201	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
				
003202 ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
				

Microbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AE0106 ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
AE1106 ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
				
AE0107 ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok Needle-free Y-site	0.9 mL	50/Case
				
AE1112	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
AE1138 ©	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case
				

Flow Controller Extension Sets

Amsino's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.


REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS107 ©	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
				
AFS109	19 in. (51 cm)	Flow Controller, NeuFlo® Needle-free Injection Site, Slide Clamp	2.83 mL	50/Case
				

Needle-Free Connectors

Our AMSafe® needle free valves are manufactured to provide both clinicians and patients save procedural outcomes.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AY0200 	Sure-Lok® Needle-free Connector	0.33 mL	100/Case
AY0206	NeuFlo® Neutral Pressure Valve	0.08 mL	100/Case



AMSafe® Prefilled Normal Saline Flush Syringe

- 0.9% Sodium Chloride Injection, USP
- 3 mL, 5 mL, 10 mL and 20 mL
- Sterile Field Ready (SFR) and Sterile Fluid Path packaging configurations







For more details and ordering information, please refer to page 39.



IV Pump Administration Sets

Sure-Lok® Needle-Free IV Pump Administration Sets

- FDA cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP092	10	92 in. (234 cm)	Roller Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP092Y1	10	92 in. (234 cm)	1 Needle-free Valve, Roller Clamp, 1 Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100FY1	10	100 in. (254 cm)	1 Needle-free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100F12Y1	10	100 in. (254 cm)	Sure-lok Needle-Free Valve, Slide and Roller clamp, 1.2 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
					



UROLOGICAL CARE

Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

Foley Catheter Trays & Collection	18
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	18
Add-a-Foley Catheter Trays	19
Pre-Connected Foley Catheter Urine Meter Kits	20
Add-a-Foley Urine Meter Trays	21
Foley Insertion Trays	21
Urinary Drainage Bags	22
Urine Meter	23
ACCU-MEASURE® Pediatric Urine Meter	23
Foley Catheters	24
Hydrophilic Latex Foley Catheters	24
Silicone-Coated Latex Foley Catheters	25
100% Silicone Foley Catheters	26
Intermittent Urinary Catheters	27
Urethral Insertion Trays and Kits	27
Closed System Urethral Drainage Kits	28
Self-Catheterization Kits with R-Polished™ Eye Catheters	29
PVC Intermittent Urethral Catheter with R-Polished™ Eyes.....	30
Red Rubber Latex Urethral Catheter	31
Irrigation Sets	32
Cystoscopy & TUR Irrigation Sets	32
Ancillary & Patient Care Products.....	33
Urinary Leg Bags	33
Extension Tube for Leg Bags	33
Pediatric Urine Collector	33
Irrigation Trays and Kits	34
Syringes	35
Specimen Containers	36

Foley Catheter Trays & Collection



AS89316

AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure® Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

Trays include: 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S ©	14 Fr	5cc	10/Case
	AS89316S ©	16 Fr		
	AS89318S ©	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89314 ©	14 Fr	5cc	10/Case
	AS89316 ©	16 Fr		
	AS89318 ©	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr	5cc	10/Case
	AS89316H	16 Fr		
	AS89318H	18 Fr		

Foley Catheter Trays & Collection



AS89300

AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10cc	10/Case
AS89330 ©	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30cc	10/Case

AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs). (HCPCS codes: See product chart below)

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag



CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89814S ©	14 Fr	5cc	10/Case
	AS89816S ©	16 Fr		
	AS89818S ©	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89814 ©	14 Fr	5cc	10/Case
	AS89816 ©	16 Fr		
	AS89818 ©	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89814H ©	14 Fr	5cc	10/Case
	AS89816H ©	16 Fr		
	AS89818H ©	18 Fr		

Foley Catheter Trays & Collection

AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. (HCPCS Code A4354)



AS89800

- Sterile, not manufactured with natural rubber latex

Trays include: 400 ml Urine Meter Preconnected with 5cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30cc	10/Case

AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)



AS880

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	AS880 ©	10cc	20/Case
	AS890 ©	30cc	
BZK Swab Sticks	AS880K ©	10cc	20/Case
	AS890K ©	30cc	

AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube. (HCPCS code A4358)

- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex



AS322 Pear



AS326 Oval



AS32600 Low-Profile

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312 ©	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322 ©	Pear	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	
AS332 ©	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	
AS302 ©	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	
AS32200 ©	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600* ©	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	
AS326* ©	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	

*100% sterile packed in peel pouch

Foley Catheter Trays & Collection

AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. (HCPCS code A4358)



- **Anti-Reflux Chamber helps improve patient safety by reducing urine backflow**
- **Twist valve drains at the bottom of the urine meter for easy drainage**
- **Leak-proof seams ensure fluid containment**
- **Needle-free sampling port reduces the risk of needle stick injuries**
- **Air vents help to maintain consistent urine flow delivering safer patient outcomes**
- **Sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS328 ©	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case

ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)



- **Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag**
- **Needle-free sampling port with luer lock access**
- **Leak-proof seams ensure fluid containment**
- **Meter bottle capacity 150 mL**
- **Adjustable bed sheet clamp for securing drainage bag tubing**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS327100	ACCU-MEASURE® Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case

Foley Catheters

AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
5cc	AS41014H©	14 Fr	16 in. (40.6 cm)	10/Box
	AS41016H©	16 Fr		
	AS41018H©	18 Fr		
30cc	AS42014H	14 Fr	16 in. (40.6 cm)	10/Box
	AS42016H©	16 Fr		
	AS42018H©	18 Fr		

Foley Catheters

AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
3cc	AS40008 ©	8Fr	10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr		
5cc	AS41012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS41014 ©	14 Fr		
	AS41016 ©	16 Fr		
	AS41018 ©	18 Fr		
	AS41020 ©	20 Fr		
	AS41022 ©	22 Fr		
	AS41024 ©	24 Fr		
	AS41026 ©	26 Fr		
	AS41028 ©	28 Fr		
	AS41030 ©	30 Fr		
30cc	AS42012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS42014 ©	14 Fr		
	AS42016 ©	16 Fr		
	AS42018 ©	18 Fr		
	AS42020 ©	20 Fr		
	AS42022 ©	22 Fr		
	AS42024 ©	24 Fr		
	AS42026 ©	26 Fr		
	AS42028 ©	28 Fr		
	AS42030 ©	30 Fr		

3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	AS43016 ©	16 Fr	16 in. (40.6 cm)	10/Box
	AS43018 ©	18 Fr		
	AS43020 ©	20 Fr		
	AS43022 ©	22 Fr		
	AS43024 ©	24 Fr		
	AS43026 ©	26 Fr		

AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
5cc	AS41012S ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS41014S ©	14 Fr		
	AS41016S ©	18 Fr		
	AS41018S ©	18 Fr		
	AS41020S ©	20 Fr		
	AS41022S ©	22 Fr		
	AS41024S ©	24 Fr		
	AS41026S ©	26 Fr		
30cc	AS42016S ©	16 Fr	16 in. (40.6 cm)	10/Box
	AS42018S ©	18 Fr		
	AS42020S ©	20 Fr		

3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	AS43018S ©	18 Fr	16 in. (40.6 cm)	10/Box
	AS43020S ©	20 Fr		
	AS43022S ©	22 Fr		
	AS43024S ©	24 Fr		
	AS43026S ©	26 Fr		

Intermittent Urinary Catheters



AS870

AMSure® Urethral Insertion Trays and Kits

AMSure® Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection

PVC Catheter Tray

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex

Red Rubber Catheter Trays

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
AS870 ©	14 Fr PVC*	20/Case
AS87114 ©	14 Fr Red Rubber Latex	
AS87118 ©	18 Fr Red Rubber Latex	

*Not manufactured with natural rubber latex.

AMSure® Closed System Urethral Drainage Kits

AMSure® Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection**
- **Sterile**

Red Rubber Catheter Kits

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

PVC Catheter Kits

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS89214R

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	LENGTH	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	AS89214R ©	14 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	20/Case
	AS89214	14 Fr	16 in. (40.6 cm)	PVC Latex- Free	
BZK Swab Sticks	AS89214RK ©	14 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	20/Case
	AS89216RK ©	16 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	

Intermittent Urinary Catheters



AS85010

AMSure® Self-Catheterization Kits with R-POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Introducer tip reduces the potential for urinary tract infection (UTI)**
- **Collection bag with built-in thumb holes allows for ease of use**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

SWAB STICK	REF	SIZE	LENGTH	QTY/UNIT
PVP Swab Sticks	AS85008 ©	8 Fr	16 in. (40.6 cm)	50/Case
	AS85010 ©	10 Fr		
	AS85012 ©	12 Fr		
	AS85014 ©	14 Fr		
	AS85016	16 Fr		
BZK Swab Sticks	AS85008K ©	8 Fr	16 in. (40.6 cm)	50/Case
	AS85010K ©	10 Fr		
	AS85012K ©	12 Fr		
	AS85014K ©	14 Fr		
	AS85016K ©	16 Fr		

AMSure® PVC Intermittent Urethral Catheter with **R**POLISHED™ Eyelets

AMSure PVC Intermittent Catheters with R-Polish technology use radio frequency to create ultra-smooth eyelets for optimal patient comfort. They are designed to provide a smooth surface in a catheter that is firm enough for insertion, yet soft enough to be slightly pliable. These urethral catheters are designed for intermittent or short-term use and are available in various French sizes, as well as female, pediatric and male lengths. A Coudé tip is also available to assist in more difficult catheterizations where strictures or other anomalies may be present in the urethra. (HCPCS code A4351)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- PVC material eliminates latex allergy concerns
- Clear shaft provides easy viewing as you measure
- Short-term catheter use helps decrease infections
- Individually packaged, sterile
- BPA-free, not manufactured with natural rubber latex



FEATURES	REF	SIZE	LENGTH	QTY/UNIT
Female	AS860608 ©	8 Fr	6 in. (15.24 cm)	50/Case
	AS860610 ©	10 Fr		
	AS860612 ©	12 Fr		
	AS860614 ©	14 Fr		
	AS860616 ©	16 Fr		
Pediatric	AS861006 ©	6 Fr	10 in. (25.40 cm)	50/Case
	AS861008 ©	8 Fr		
	AS861010 ©	10 Fr		
Male	AS861608 ©	8 Fr	16 in. (40.64 cm)	50/Case
	AS861610 ©	10 Fr		
	AS861612 ©	12 Fr		
	AS861614 ©	14 Fr		
	AS861616 ©	16 Fr		
	AS861618 ©	18 Fr		

Intermittent Urinary Catheters

AMSure Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free



REF	SIZE	LENGTH	QTY/UNIT
AS44008 ©	8 Fr	11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr		
AS44012 ©	12 Fr	16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr		
AS44016 ©	16 Fr		
AS44018 ©	18 Fr		
AS44020 ©	20 Fr		
AS44022 ©	22 Fr		
AS44024 ©	24 Fr		
AS44028 ©	28 Fr		

Irrigation Sets

AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
ABL400 ©	Single Lead	20/Case



REF	TYPE	QTY/UNIT
ABL402 ©	Two Leads	20/Case



REF	TYPE	QTY/UNIT
ABL404 ©	Four Leads	20/Case

Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex

SIZE	REF	DRAIN PORT STYLE	QTY/UNIT
600 mL	AS306N ©	Push-Pull	48/Case
	AS308N ©	Twist-Turn	
900 mL	AS307N ©	Push-Pull	48/Case
	AS309N ©	Twist-Turn	



AS307N

Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS311 ©	18 in. (45.72 cm) Extension Tube	50/Case



Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS409 ©	200 mL Pediatric Urine Bag	50/Box



Irrigation Trays and Kits

AMSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need. (HCPCS code A4320)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector



AS121

Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap



AS126

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector

Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			
AS127	•		•	

Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		



AS136

Ancillary & Patient Care Products

Syringes

(HCPCS code A4322)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap

AS010

AS016NS

AS115

TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS011 ©	Poly Pouch	Sterile (EO)	•		50/Case
	AS010 ©	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	AS015 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016NS	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	AS115 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116	Resealable Pole Bag	Non-Sterile	•	•	30/Case

AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.
(HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
AS340 ©		•		•		•	•	100/Case
AS341 ©			•	•	•	•	•	100/Case
AS343 ©			•				•	500/Case
AS345 ©	•			•			•	100/Case



AS340





Inhalation, Irrigation, and Injection Solutions

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Sterile Water for Inhalation, USP	38
Sterile Water Prefilled Humidifiers, USP	38
Nebulizer Kits	63
Saline Products	39
Device and Wound Irrigation	40

AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.

- **Manufactured in USA**

Sterile Water for Inhalation (Bottles)

- 500 mL and 1000 mL bottles Sterile Water, USP
- For use with nebulizers and humidifiers
- Hanging Spikeable bottle helps solve the problem of look-alike solution “bags” and prevents mishaps from occurring – Use with vented spikes only



AS1065

AS1000

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	—	500 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	

Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case (10 Boxes – 100 per Box)
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	

Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

Saline Products

AMSafe® Prefilled Saline Flush Syringe

- 3 mL 5 mL, 10 mL and 20 mL syringes, 0.9% Sodium Chloride Injection, USP
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared
- Sterile, non-pyrogenic, DEHP-free, PVC-free, not manufactured with natural rubber latex



12 mL and 20 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1203 ☺	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205 ☺	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210 ☺	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020 ☺	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
IVF1203SS ☺	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205SS ☺	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210SS ☺	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020SS ☺	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

10 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)



REF	DESCRIPTION	QTY/UNIT
IVF1003 ☺	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1005 ☺	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1010 ☺	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1010TM	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

Manufactured in
North America

MEDXL
An **AMSINO** Company

Pre-Filled Saline Flush Syringes



Praxiject SP Flush Syringe

REF	DESCRIPTION	QTY/UNIT
3705-3US1	Praxiject Sterile Path, 10 mL of 0.9% Sodium Chloride Injection USP in a 10 mL syringe – sterile fluid pathway	180/Case (6 Boxes – 30 per Box)

Praxiject SF Flush Syringe

REF	DESCRIPTION	QTY/UNIT
3705CUS	Praxiject Sterile Field, 10 mL of 0.9% Sodium Chloride Injection USP in a 10 mL syringe – sterile field compatible	600/Case (6 Boxes – 100 per Box)

AMSure® Sterile Sodium Chloride Inhalation Solution (vials)

- 3 mL, 5 mL and 15 mL Sterile Sodium Chloride Inhalation Solution (0.9%), USP
- Single-use vials with graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case (10 Boxes – 100 per Box)
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case (10 Boxes – 100 per Box)
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case (3 Boxes – 48 per Box)



ASR0039



ASR0059



ASR0159

AMSure Device and Wound Irrigation

- 100 mL sterile water or saline bottle with screw cap and inner seal
- Easy open twist screw top container
- Used for device and wound irrigation

REF	DESCRIPTION	QTY/UNIT
AS399 ©	100 mL Sterile Water for Device and Wound Irrigation	48/Case
AS398 ©	100 mL Sterile Saline 0.9% for Device and Wound Irrigation	48/Case



AS399



AS398



AMSure Sal-Sense™ Single Use Saline Topical Solution 0.9% Sodium Chloride

- Single patient use avoids cross contamination
- Save time from standard practice of drawing solution
- Twist design to open effortlessly
- Easy to irrigate area in any position



REF	DESCRIPTION	QTY/UNIT
AS0159	15 mL Single Use Saline Topical Solution	144/Case (3 Boxes – 48 per Box)
AS0359	35 mL Single Use Saline Topical Solution	120/Case (6 Boxes – 20 per Box)
AS0559	50 mL Single Use Saline Topical Solution	120/Case (6 Boxes – 20 per Box)

RECEPTAL®

WASTE FLUID MANAGEMENT SYSTEMS



Waste Fluid Management Solutions


Amsino is a leading global provider of suction fluid management systems. The breadth of Amsino's RECEPTAL suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

EZE-VAC® Suction Canister System	42
EZE-VAC 2.0 Suction Canister System.....	43
RECEPTAL Liners	44
RECEPTAL Twist Liners	45
RECEPTAL Outer Shells	46
ACCU-MEASURE® Critical Measurement Devices	47
RECEPTAL Solidifiers	47
AMSure® Suction Tubing with Connectors.....	48
Suction Yankauers and Kits	49
SteriCadi™ Suction Device Holster System	50
Accessories.....	51
Reusable Hardware	54
iRECEPTAL® Surgical Waste Fluid Management System	57

EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01 ©	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01 ©	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01 ©	3300 mL Canister, Lid and 2 Elbows	12/Case

EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
43210-01	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID x 6 ft. (5 mm x 1.8 m) + 3/16 in ID x 18 in. (5 mm x 46 cm)	30/Case

RECEPTAL Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
1000	BULK PACKAGED		
	43010-05 ©	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	
1500	BULK PACKAGED		
	43025-01 ©	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	
	43015S ©	Liner with VAC-GARD and 45g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	
	0L654-01 ©	Liner with VAC-GARD and Pour Spout	
2000	BULK PACKAGED		
	43027-01 ©	Liner with Pour Spout,	50/Case
	43041-01 ©	Liner with Ortho Port	
	43042-01 ©	Liner with VAC-GARD	
	43044-05 ©	Liner with VAC-GARD and Pour Spout	
	43020S ©	Liner with VAC-GARD and 60g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	
	0G905-01 ©	Liner with Ortho Port	
	0P507-01 ©	Liner with Pour Spout	
3000	BULK PACKAGED		
	0L656-01*	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	0L655-01*	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	
	43030S* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	

* Dual port allows for dual suctioning during procedures.
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID

RECEPTAL® Twist Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.


- **Premeasured, preloaded and contact-free Solidifier Dispenser Cap**
- **Integrated solidifier eliminates the need for open ports**
- **All-in-one solidifier and liner solution reduces costs associated with extra inventory**
- **100% leak tested seamless liners minimize risk of contamination**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**



43121S



43120S shown above
in canister 43445-01
(canister not included)



REF	DESCRIPTION	QTY/UNIT
43110S	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43115S ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
43116S ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	
43120S ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
43121S ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	

RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex



Color-coded by size for easy selection.

REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	
43445-01 ©	2000 mL Outer Shell	
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

ACCU-MEASURE® Critical Measurement Devices

ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL® ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

RECEPTAL
ACCU-MEASURE
Device



REF	DESCRIPTION	FEATURES	QTY/UNIT
43306-01	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case

RECEPTAL® Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- No shaking required
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4320SP

RECEPTAL Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
4305SP	500 Solidifier Pouch, solidifies 500 mL	96/Case
4312SP	1200 Solidifier Pouch, solidifies 1200 mL	
4315SP	1500 Solidifier Pouch, solidifies 1500 mL	
4320SP	2000 Solidifier Pouch, solidifies 2000 mL	
4330SP	3000 Solidifier Pouch, solidifies 3000 mL	

RECEPTAL Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
4305SB	500 Solidifier Bottle, solidifies 500 mL	96/Case
4312SB	1200 Solidifier Bottle, solidifies 1200 mL	
4315SB	1500 Solidifier Bottle, solidifies 1500 mL	
4320SB	2000 Solidifier Bottle, solidifies 2000 mL	
4330SB	3000 Solidifier Bottle, solidifies 3000 mL	
43160SB	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case



4315SB

4320SB

AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex




AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 [®] (50/Case)	AS820B [®] (100/Case)
	6 ft. (1.8 m)	AS821 [®] (50/Case)	AS821B [®] (50/Case)
	10 ft. (3.0 m)	AS822 [®] (50/Case)	AS822B [®] (50/Case)
	12 ft. (3.7 m)	AS823 [®] (20/Case)	AS823B [®] (50/Case)
	20 ft. (6.1 m)	AS824 [®] (20/Case)	AS824B [®] (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 [®] (50/Case)	AS825B [®] (50/Case)
	10 ft. (3.0 m)	AS826 [®] (50/Case)	AS826B [®] (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C [®] (50/Case)	—
	12 ft. (3.7 m)	AS827 [®] (20/Case)	AS827B [®] (50/Case)
	20 ft. (6.1 m)	AS828 [®] (20/Case)	AS828B [®] (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 [®] (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M [®] (50/Case)	AS825BM [®] (50/Case)
	10 ft. (3.0 m)	AS826M [®] (50/Case)	AS826BM [®] (50/Case)
	12 ft. (3.7 m)	AS827M [®] (20/Case)	AS827BM [®] (50/Case)
	20 ft. (6.1 m)	AS828M [®] (20/Case)	AS828BM [®] (50/Case)

Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



Yankauers

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			AS830 © 50/Case	AS830B 500/Case
	•		•	AS831 © 50/Case	AS831B 500/Case
		•		AS832 © 50/Case	AS832B 500/Case
		•	•	AS833 © 50/Case	AS833B 500/Case
Flexible	•			AS834 © 50/Case	AS834B 500/Case
	•		•	AS835 © 50/Case	AS835B 500/Case
		•		AS836 © 50/Case	AS836B 500/Case
		•	•	AS837 © 50/Case	AS837B 500/Case

Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
AS825P ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
AS827P* ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case


*Made to order

SteriCadi™ Suction Device Holster System

SteriCadi is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to support improved infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder, Sterile	30/Case



Accessories



RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
43303-01	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
43303-02	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	
43303-05	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	



External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
43309-01	External Shut-Off Valve, ¼ in. (0.64 cm) Male × Female Elbow	10/Case



RECEPTAL® Accu-Measure® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
43305-10	Canister Attachment Bracket for 250 mL Canister	10/Case

Suction Tubing

For permanent Outer Shell-to-vacuum source,
or Outer Shell-to-Outer Shell connections.

REF	DESCRIPTION	QTY/UNIT
43455-01	100 ft. (30.48 m) PVC Non-Sterile SuctionTubing 9/32 in. (7 mm) I.D.	1/Case
43457-25	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case



EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.

- Non-sterile, BPA-free, not manufactured with natural rubber latex


REF	DESCRIPTION	QTY/UNIT
43317-25	18 in. (0.45 m) Tandem Tubing	25/Case



RECEPTAL Outer Shell Connection Tubing

RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum “T” of the RECEPTAL canisters between them for easy tandem connection.

REF	DESCRIPTION	QTY/UNIT
CTM50B ©	1/4 in. x 20 in. (6.35 mm x 500 mm), Female Connector	100/Case
CTM20B ©	1/4 in. x 8 in. (6.35 mm x 200 mm), Female Connector	300/Case



Tubing Manifold

For connection of 2, 3 or 4 canisters.

REF	DESCRIPTION	QTY/UNIT
TMF001	Tubing Manifold	10/Case



Accessories

RECEPTAL® Elbow Adapters and Taper

Use with Patient Port.



REF	DESCRIPTION	QTY/UNIT
43318-01	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
43316-01	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
43494-25	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
43473-25	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case

RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.



REF	DESCRIPTION	QTY/UNIT
44300-05	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	100/Case
43448-25	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case

5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS807 ©	5-in-1 Connector, Sterile	100/Case
AS807B	5-in-1 Connector, Non-sterile	500/Case

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case

Reusable Hardware

Roll Stands

REF	DESCRIPTION	QTY/UNIT
43474-01	Large head, stainless steel Adjustable Roll Stand, conductive casters, holds up to 12 RECEPTAL® or 8 EZE-VAC® Canisters, use with 1500 mL, 2000 mL and 3000 mL canisters	1/Case
43475-01	Small head, stainless steel Adjustable Roll Stand, conductive casters, holds multiple RECEPTAL and EZE-VAC Canisters	1/Case
43476-P2	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 1500 mL and 2000 mL canisters	1/Case
43476-P3	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 3000 mL canisters	1/Case



Single and In-Line Floor Stands

REF	DESCRIPTION	QTY/UNIT
43431-01	Single Canister Floor Stand with tip-resistant base, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43432-01	Dual Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43447-01	In-Line Quad Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case



Reusable Hardware

RECEPTAL® Liner System Tandem Canister Supports



REF	DESCRIPTION	QTY/UNIT
43429-01*	For 1500 mL and 2000 mL canisters	1/Case
43439-01	For 3000 mL canisters	

RECEPTAL Liner System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43425-05*	For 1000 mL canister, metal	5/Case
43424-10	For 1000 mL canister support, with mounting plate, metal	5/Case
43426-05	For 1500 mL and 2000 mL canisters, metal	5/Case
43438-10	For 3000 mL canisters, with mounting plate for vacuum regulator, metal	10/Case
43428-05*	For 1500 mL and 2000 mL canisters, with mounting plate for vacuum regulator, metal	5/Case
43436-05	For 3000 mL canisters, metal	5/Case
43426P	For 1000 mL, 1500 mL, and 2000 mL canisters, plastic	20/Case
43436P	For 3000 mL canisters, plastic	10/Case
43424-V	For 1000/1500 mL canister support ring with valve	10/Case
43428-V	For 2000 mL canister support ring with valve	10/Case
43438-V	For 3000 mL canister support ring with valve	10/Case

EZE-VAC® Rigid System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43467-01	For EZE-VAC 2.0 800 mL and 1200 mL canisters, metal, small: ID=4.4 in. (11.7 cm), wall to outer edge = 8.5 in. (21.6 cm)	12/Case
43469-12*	For EZE-VAC 1200 mL, 2000 mL and 3300 mL canisters, metal, large: ID=5.6 in. (14.2 cm), wall to outer edge = 9.4 in. (23.9 cm)	12/Case
43469-V	For EZE-VAC 1200 mL, 2000 mL, and 3300 mL with valve	10/Case

*Shown on Wall Bracket (REF 43444-05)

Canister Brackets

REF	DESCRIPTION	QTY/UNIT
43444-05	Wall channel bracket for mounting canister, stainless steel	5/Case
43444P	Wall channel bracket for mounting canister, polycarbonate	50/Case
44303-10	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
44304-05	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
RCB001	Metal Rail Bracket with Universal Holder Clamp	30/Case



iRECEPTAL® Surgical Waste Fluid Management System

The iReceptal Digital Surgical Suction System is an all-in-one, closed suction, container, and disposal system. This advanced, easy-to-use surgical suction solution is completely self-contained to provide for independent collection and disposal of waste fluids, along with a docking station for automatic and closed emptying of the rover.

iReceptal is currently available in EMEA and China markets. iReceptal is not currently cleared by the US FDA.

- Constantly fully closed system, from collection to disposal
- Staff is never exposed to potentially infectious fluid waste
- Helps optimize surgical turnaround times
- Manifold provides multiple suction lines at the same time with backflow protection
- Separates fluid from tissue and bone fragments
- Waste fluid discharges directly into facility's drain system
- High volume capacity of 20L
- Volume gradation indications
- Easy to maneuver with large wheels
- Lowers disposal costs by eliminating full suction canisters in the waste stream
- Automatically cleans when docked in under 5 minutes

The system is comprised of:

- Mobile Rover
- Docking Station
- Disposable Suction Manifold
- Single-Use Smoke Filter
- Cleaning Detergent




Mobile Rover



Docking Station

REF	DESCRIPTION	QTY/UNIT
iRR100-01 ©	iRECEPTAL 1.0 Rover	1/Each
iRR150-01 ©	iRECEPTAL 1.5 Rover with 7cm (2.7") LCD Screen	
iRD102-01 ©	iRECEPTAL Docking Station (220VAC)	
iRD101-01	iRECEPTAL Docking Station (110VAC)	
iRM100-01 ©	4-Port Manifold (disposable)	
iRC001-01	Detergent for Docking Station	
iRC002-01	Detergent for Rover	



Manifold Storage

- 4 ports for suction tubing with backflow control



Intelligent Control Panel

- Fluid read-out (Version 1.5 only)



Electromagnetic Couplers

- Allow for easy docking



Direct to Drain Systems

- Hands-free disposal and self-cleaning

Completely Closed System

- Avoid direct touch of waste fluid

20L Volume

- Suitable for large amounts of fluid





Enhancing Oxygen Delivery

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO₂ levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Respiratory Masks and Nasal Cannulas.....	60
Nebulizer Kits	63
Accessories.....	64
Aerosol Drainage Bag	64
Heat and Moisture Exchangers	65
Tracheostomy Clean and Care Trays	65
Sterile Water and Saline	66
Suction Catheters, Trays and Kits	67
AMSure Suction Tubing with Connectors.....	70
Suction Yankauers and Kits	71

Respiratory Masks and Nasal Cannulas



AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS74010	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

Respiratory Masks and Nasal Cannulas

High Flow Nasal Oxygen Cannula

Amsino's AMSure® high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most high flow nasal cannula devices
- Single patient use, not manufactured with natural rubber latex

AS75280




REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75280 ©	High Flow Nasal Oxygen Cannula	Small, Adult	10/Case
AS75285	High Flow Nasal Oxygen Cannula	Medium, Adult	
AS75282 ©	High Flow Nasal Oxygen Cannula	Large, Adult	

Nasal Oxygen Cannula

AMSure Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

AS75080




REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75080 ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
AS75085 ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	

Cushion Nasal Oxygen Cannula

AMSure Cushion Nasal Oxygen Cannulas with soft foam ear covers provide low-to-medium flow oxygen delivery. Soft, pliable prongs are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75181 ©	Soft Nasal Oxygen Cannula	Curved Prongs with Cushion, Adult	25/Case

Soft Nasal Oxygen Cannula

AMSure® Soft Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75180 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Adult	50/Case
AS75190 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Pediatric	
AS75199 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Infant	



AS75180

Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75010	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
AS75020	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	



AS75010

Pediatric Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75050	Pediatric Aerosol Mask	50/Case



AS75050

Respiratory Masks and Nasal Cannulas




Tracheostomy Mask

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060	Adult Tracheostomy Mask	50/Case



AS78010



AS78015

Nebulizer Kits

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78015 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece	
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	
AS78025 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece, Reusable	

Accessories

Respiratory Adapters

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS730100	Universal Cuff Adapter	22 mm ID Connector	50/Case
AS730200	Multi Adapter	15 mm ID / 22 mm OD Connector	



Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- Pre-connected self-locking plastic chain secures the system to a stable object
- Convenient drain port allows for easy one-handed drainage
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75035	750 mL Aerosol Drainage Bag	50/Case



Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
HME001 ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
HME002 ©	Heat and Moisture Exchangers, Sterile	

Tracheostomy Clean and Care Trays



AS861

AMSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS860 ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
AS861 ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	

AMSure® Sterile Water and Saline

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.

- **Manufactured in USA**

Sterile Water for Inhalation (Bottles)

- 500 mL and 1000 mL bottles Sterile Water, USP
- For use with nebulizers and humidifiers
- Hanging Spikeable bottle helps solve the problem of look-alike solution “bags” and prevents mishaps from occurring – Use with vented spikes only



AS1065

AS1000

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	—	500 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	

Unit Dose Vials for Inhalation

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation
- Available in both sterile water and sterile sodium chloride inhalation solution (0.9%), USP



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case (10 Boxes – 100 per Box)
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case (3 Boxes – 48 per Box)

Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing




AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

Suction Catheters, Trays and Kits



AS361

Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- Opposing staggered eyelets reduce the potential for blocked suction
- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	STRAIGHT PACK	COILED PACK	QTY/UNIT
AS361 ©	6 Fr	17½ in.	•		50/Case
AS361C ©	6 Fr			•	
AS362 ©	8 Fr		•		
AS362C ©	8 Fr			•	
AS363 ©	10 Fr	23 in.	•		
AS363C ©	10 Fr			•	
AS364 ©	12 Fr		•		
AS364C ©	12 Fr			•	
AS365 ©	14 Fr		•		
AS365C ©	14 Fr			•	
AS366 ©	16 Fr		•		
AS366C ©	16 Fr			•	
AS367 ©	18 Fr	23 in.	•		
AS367C ©	18 Fr			•	

Suction Catheter Kits

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. Kits include coiled packed suction catheters in various French sizes, a pop-up solution cup with sterile inner wrap, and one or two vinyl gloves.

- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex



AS371

REF	SIZE	LENGTH	VINYL GLOVES		QTY/UNIT
			ONE GLOVE	TWO GLOVES	
AS371 ©	6 Fr	17 ½ in.	•		50/Case
AS381 ©	6 Fr			•	
AS372 ©	8 Fr		•		
AS382 ©	8 Fr			•	
AS373 ©	10 Fr	23 in.	•		
AS383 ©	10 Fr			•	
AS374 ©	12 Fr		•		
AS384 ©	12 Fr			•	
AS375 ©	14 Fr		•		
AS385 ©	14 Fr			•	
AS376 ©	16 Fr		•		
AS386 ©	16 Fr			•	
AS387 ©	18 Fr			•	



SCT10

Mini Suction Catheter Trays

The AMSure® Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- Opposing staggered eyelets reduce the potential for blocked suction
- DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients
- Whistle tip enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	DESCRIPTION	DELEE TIP	WHISTLE TIP	QTY/UNIT
SCT08	8 Fr	17 ½ in.	Suction Catheter, Tray and Vinyl Gloves	•		100/Case
SCT10	10 Fr	23 in.		•		
SCT12	12 Fr	23 in.			•	
SCT14	14 Fr	23 in.			•	
SCT16	16 Fr	23 in.			•	

Closed Suction Catheter Trays

72-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

- Closed system to prevent cross infections in facility
- Clean easily with needle free port
- Avoid interruptions of oxygen supply due to suction



VS365TP124

TYPE	REF	SIZE	LENGTH	ADAPTER TYPE		FLEX CONNECTOR		QTY/UNIT
				Double Swivel Elbow	T-Piece	15 mm x 22 mm	22 mm x 22 mm	
Endotracheal Kits	VS363RC124	10 Fr	21.3 in. (54 cm)	•		•		20/Case
	VS364RC124	12 Fr		•		•		
	VS365RC124	14 Fr		•		•		
	VS363TP124	10 Fr			•		•	
	VS364TP124	12 Fr			•		•	
	VS365TP124	14 Fr			•		•	
Tracheostomy Kits	VS363TC124	10 Fr	12 in. (30.5 cm)	•		•		
	VS364TC124	12 Fr		•		•		
	VS365TC124	14 Fr		•		•		
	VS363TTP124	10 Fr			•		•	
	VS364TTP124	12 Fr			•		•	
	VS365TTP124	14 Fr			•		•	

Pediatric Closed Suction Catheter Trays

72-hour Catheter from 5 FR to 8 FR Catheter Kits with MDI Port Adapter.

- Closed system to prevent cross infections in facility
- Clean easily with needle free port
- Avoid interruptions of oxygen supply due to suction

REF	SIZE	LENGTH	DESCRIPTION	QTY/UNIT
VS360RC72	5 Fr	12 in. (30.5 cm)	Endotracheal Y-Adapters (2.0 and 2.5 mm)	20/Case
VS361RC72	6 Fr		Endotracheal Y-Adapters (2.5, 3.0 and 3.5 mm)	
VS3615RC72	7 Fr		Endotracheal Y-Adapters (3.0, 3.5 and 4.0 mm)	
VS362RC72	8 Fr		Endotracheal Y-Adapters (3.0, 3.5 and 4.0 mm)	

Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 © (50/Case)	AS820B © (100/Case)
	6 ft. (1.8 m)	AS821 © (50/Case)	AS821B © (50/Case)
	10 ft. (3.0 m)	AS822 © (50/Case)	AS822B © (50/Case)
	12 ft. (3.7 m)	AS823 © (20/Case)	AS823B © (50/Case)
	20 ft. (6.1 m)	AS824 © (20/Case)	AS824B © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 © (50/Case)	AS825B © (50/Case)
	10 ft. (3.0 m)	AS826 © (50/Case)	AS826B © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C © (50/Case)	—
	12 ft. (3.7 m)	AS827 © (20/Case)	AS827B © (50/Case)
	20 ft. (6.1 m)	AS828 © (20/Case)	AS828B © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 © (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M © (50/Case)	AS825BM © (50/Case)
	10 ft. (3.0 m)	AS826M © (50/Case)	AS826BM © (50/Case)
	12 ft. (3.7 m)	AS827M © (20/Case)	AS827BM © (50/Case)
	20 ft. (6.1 m)	AS828M © (20/Case)	AS828BM © (50/Case)

Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS835

AS832

Yankauers

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			AS830 © 50/Case	AS830B 500/Case
	•		•	AS831 © 50/Case	AS831B 500/Case
		•		AS832 © 50/Case	AS832B 500/Case
		•	•	AS833 © 50/Case	AS833B 500/Case
Flexible	•			AS834 © 50/Case	AS834B 500/Case
	•		•	AS835 © 50/Case	AS835B 500/Case
		•		AS836 © 50/Case	AS836B 500/Case
		•	•	AS837 © 50/Case	AS837B 500/Case

Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
AS825P ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
AS827P* ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

*Made to order

Suction Device Holster System

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder with Clip	30/Case





Everyday Patient Care Solutions

AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

Specimen Containers..... 73

Vaginal Speculums..... 73

Disposable Umbilical Cord Clamp..... 73

Ear/Ulcer Syringes 74

Cleansing Enema Sets..... 74

Connectors..... 76

Specimen Containers

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**



AS340

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
AS340 ©		•		•		•	•	100/Case
AS341 ©			•	•	•	•	•	100/Case
AS343 ©			•				•	500/Case
AS345 ©	•			•			•	100/Case

Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- **Round-edge design helps to increase patient comfort**
- **Single-use disposable**
- **Non-sterile, not manufactured with natural rubber latex**



AS032M

REF	DESCRIPTION	QTY/UNIT
AS032S	Small Vaginal Speculum	100/Case
AS032M	Medium Vaginal Speculum	100/Case
AS032L	Large Vaginal Speculum	100/Case

Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- **Safety lock is designed to decrease risk of accidental disconnection**
- **Grooved teeth help provide an effective grip**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



REF	DESCRIPTION	QTY/UNIT
UCC100	Umbilical Cord Clamp, sterile	50/Case

Ear/Ulcer Syringes

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS00501 ©	1 oz. syringe	Non-sterile, Bulk	300/Case
AS00502 ©	2 oz. syringe	Non-sterile, Bulk	200/Case
AS00503 ©	3 oz. syringe	Non-sterile, Bulk	150/Case
AS00501S ©	1 oz. syringe	Sterile, Single Pack	50/Case
AS00502S ©	2 oz. syringe	Sterile, Single Pack	50/Case
AS00503S ©	3 oz. syringe	Sterile, Single Pack	50/Case



Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS330 ©	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case



AS330



AS333




Connectors

5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3 /16 in., ¼ in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS807 	5-in-1 Connector	Sterile	100/Case
AS807B	5-in-1 Connector	Non-sterile, Bulk	500/Case



AS807

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606.....	13	43317-25.....	52	AS032L.....	73	AS374.....	68
003201.....	13	43423-01.....	46	AS032M.....	73	AS375.....	68
003202.....	13	43424-10.....	55	AS032S.....	73	AS376.....	68
0G901-01.....	44	43424-V.....	55	AS0352.....	38, 66	AS381.....	68
0G902-01.....	44	43425-05.....	55	AS0352E.....	38, 66	AS382.....	68
0G903-01.....	44	43426-05.....	55	AS0359.....	40	AS383.....	68
0G904-01.....	44	43426P.....	55	AS0552.....	38, 66	AS384.....	68
0G905-01.....	44	43428-05.....	55	AS0552E.....	38, 66	AS385.....	68
0L654-01.....	44	43428-V.....	55	AS0559.....	40	AS386.....	68
0L655-01.....	44	43429-01.....	55	AS1000.....	38, 66	AS387.....	68
0L656-01.....	44	43431-01.....	54	AS1065.....	38, 66	AS398.....	40
0P507-01.....	44	43432-01.....	54	AS1065E.....	38, 66	AS399.....	40
0P508-01.....	44	43436-05.....	55	AS115.....	4, 39	AS409.....	33
1010001B.....	9	43436P.....	55	AS116.....	4, 39	AS40008.....	25
107201.....	7	43438-10.....	55	AS121.....	5, 34	AS40010.....	25
108301.....	7	43438-V.....	55	AS126.....	5, 34	AS41012.....	25
108302.....	9	43439-01.....	55	AS127.....	5, 34	AS41012S.....	26
108303.....	7	43444-05.....	56	AS130.....	5, 34	AS41014.....	25
108304.....	7	43444P.....	56	AS136.....	5, 34	AS41014H.....	24
108305.....	8	43445-01.....	46	AS1565.....	38, 66	AS41014S.....	26
108306.....	7	43447-01.....	54	AS302.....	22	AS41016.....	25
108308.....	8	43448-25.....	53	AS306N.....	33	AS41016H.....	24
108310.....	7	43449-05.....	46	AS307N.....	33	AS41016S.....	26
109602.....	7	43455-01.....	52	AS308N.....	33	AS41018.....	25
1510501.....	9	43457-25.....	52	AS309N.....	33	AS41018H.....	24
158305.....	9	43467-01.....	55	AS43016.....	25	AS41018S.....	26
158306.....	9	43469-12.....	55	AS43018.....	25	AS41020.....	25
158309.....	9	43469-V.....	55	AS43018S.....	26	AS41020S.....	26
159602.....	9	43473-25.....	53	AS43020.....	25	AS41022.....	25
3705-3US1.....	39	43474-01.....	54	AS43020S.....	26	AS41022S.....	26
3705CUS.....	39	43475-01.....	54	AS43022.....	25	AS41024.....	25
43010-05.....	44	43476-P2.....	54	AS43022S.....	26	AS41024S.....	26
43010S.....	44	43476-P3.....	54	AS43024.....	25	AS41026.....	25
43015S.....	44	43494-25.....	53	AS43024S.....	26	AS41026S.....	26
43020S.....	44	43503-01.....	46	AS43026.....	25	AS41028.....	25
43025-01.....	44	44300-05.....	53	AS43026S.....	26	AS41030.....	25
43027-01.....	44	44303-10.....	56	AS3100.....	10	AS42012.....	25
43030S.....	44	44304-05.....	56	AS3101.....	10	AS42014.....	25
43041-01.....	44	607201.....	11	AS3102.....	10	AS42014H.....	24
43042-01.....	44	608304.....	11	AS3110A.....	10	AS42016.....	25
43044-01.....	44	608305.....	11	AS312.....	22	AS42016H.....	24
43044-05.....	44	608306.....	11	AS322.....	22	AS42016S.....	26
43055B.....	47	608310.....	11	AS32200.....	22	AS42018.....	25
43055P.....	47	AA3101.....	9	AS326.....	22	AS42018H.....	24
43110S.....	45	ABL400.....	32	AS328.....	23	AS42018S.....	26
43115S.....	45	ABL402.....	32	AS32600.....	22	AS42020.....	25
43116S.....	45	ABL404.....	32	AS327100.....	23	AS42020S.....	26
43120S.....	45	AC3110.....	7	AS330.....	81	AS42022.....	25
43121S.....	45	AE0106.....	17	AS333.....	81	AS42024.....	25
43125B.....	47	AE0107.....	17	AS340.....	36, 73	AS42026.....	25
43125P.....	47	AE1106.....	17	AS341.....	36, 73	AS42028.....	25
43155B.....	47	AE1112.....	17	AS343.....	36, 73	AS42030.....	25
43155P.....	47	AE1138.....	17	AS345.....	36, 73	AS44008.....	31
43160SB.....	47	AE3107.....	13	AS3500.....	13	AS44010.....	31
43202-01.....	43	AE3108.....	13	AS3515.....	13	AS44012.....	31
43202-820821.....	43	AE3108B.....	13	AS361.....	67	AS44014.....	31
43203-01.....	42	AFS102.....	8	AS361C.....	67	AS44016.....	31
43206-01.....	42	AFS107.....	17	AS362.....	67	AS44018.....	31
43205B.....	47	AFS109.....	17	AS362C.....	67	AS44020.....	31
43205P.....	47	AFS208.....	8	AS363.....	67	AS44022.....	31
43210-01.....	43	AFS208B.....	8	AS363C.....	67	AS44024.....	31
43213-01.....	42	AS00501.....	74	AS364.....	67	AS44028.....	31
43303-01.....	51	AS00501S.....	74	AS364C.....	67	AS730100.....	64
43303-02.....	51	AS00502.....	74	AS365.....	67	AS730200.....	64
43303-05.....	51	AS00502S.....	74	AS365C.....	67	AS74010.....	60
43305-10.....	47	AS00503.....	74	AS366.....	67	AS75010.....	62
43305-10.....	51	AS00503S.....	74	AS366C.....	67	AS75020.....	62
43306-01.....	47	AS010.....	4, 35	AS367.....	67	AS75035.....	64
43309-01.....	51	AS011.....	4, 35	AS367C.....	67	AS75050.....	62
43305B.....	47	AS015.....	35	AS371.....	68	AS75060.....	63
43305P.....	47	AS0159.....	40	AS372.....	68	AS75080.....	61
43316-01.....	53	AS016NS.....	35	AS373.....	68	AS75085.....	61

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
AS75180	62	AS860612	30	ENS015NS	4		
AS75181	61	AS860614	30	ENS115	4		
AS75190	62	AS860616	30	ENS115NS	4		
AS75199	62	AS861	65	FP092	16		
AS75280	61	AS861006	30	FP092Y1	16		
AS75282	61	AS861008	30	FP100BY2	16		
AS75285	61	AS861010	30	FP100FY1	16		
AS78010	63	AS861608	30	FP100F12Y1	16		
AS78015	63	AS861610	30	FP106FY2	16		
AS78020	63	AS861612	30	FP109BY3	16		
AS78025	63	AS861614	30	HME001	65		
AS807	53, 76	AS861616	30	HME002	65		
AS807B	53, 76	AS861618	30	iRC001-01	57		
AS820	48, 70	AS870	27	iRC002-01	57		
AS820B	48, 70	AS87114	27	iRD101-01	57		
AS821	48, 70	AS87118	27	iRD102-01	57		
AS821B	48, 70	AS880	21	iRM100-01	57		
AS822	48, 70	AS880K	21	iRR100-01	57		
AS822B	48, 70	AS890	21	iRR150-01	57		
AS823	48, 70	AS890K	21	IVF1003	39		
AS823B	48, 70	AS89214	28	IVF1005	39		
AS824	48, 70	AS89214R	28	IVF1010	39		
AS824B	48, 70	AS89214RK	28	IVF1010TM	39		
AS825	48, 70	AS89216RK	28	IVF1203	39		
AS825B	48, 70	AS89300	19	IVF1203SS	39		
AS825BM	48, 70	AS89314	18	IVF1205	39		
AS825M	48, 70	AS89314H	18	IVF1205SS	39		
AS825P	49, 71	AS89314S	18	IVF1210	39		
AS826	48, 70	AS89316	18	IVF1210SS	39		
AS826B	48, 70	AS89316H	18	IVF2020	39		
AS826BM	48, 70	AS89316S	18	IVF2020SS	39		
AS826C	48, 70	AS89318	18	OM001	60		
AS826M	48, 70	AS89318H	18	PE90B12T	2		
AS827	48, 70	AS89318S	18	PE90B5T	2		
AS827B	48, 70	AS89330	19	PE90ST	2		
AS827BM	48, 70	AS89800	21	PG001	2		
AS827M	48, 70	AS89814	20	RCB001	62		
AS827P	49, 71	AS89814H	20	S-1000-SI	3		
AS828	48, 70	AS89814S	20	SCT08	68		
AS828B	48, 70	AS89816	20	SCT10	68		
AS828BM	48, 70	AS89816H	20	SCT12	68		
AS828M	48, 70	AS89816S	20	SCT14	68		
AS829B-01	48, 70	AS89818	20	SCT16	68		
AS830	49, 71	AS89818H	20	SIH0001	50, 71		
AS830B	49, 71	AS89818S	20	TMF001	52		
AS831	49, 71	AS89830	21	UCC100	73		
AS831B	49, 71	ASR0035	38, 66	UYC001	53, 76		
AS832	49, 71	ASR0035E	38, 66	UYC001S	53, 76		
AS832B	49, 71	ASR0039	40	VS360RC72	69		
AS833	49, 71	ASR0039E	40	VS361RC72	69		
AS833B	49, 71	ASR0055	38, 66	VS3615RC72	69		
AS834	49, 71	ASR0055E	38, 66	VS362RC72	69		
AS834B	49, 71	ASR0059	40	VS363RC124	69		
AS835	49, 71	ASR0059E	40	VS363TC124	69		
AS835B	49, 71	ASR0159	40	VS363TP124	69		
AS836	49, 71	ASR0159E	40	VS363TTP124	69		
AS836B	49, 71	AY0200	15	VS364RC124	69		
AS837	49, 71	AY0206	15	VS364TC124	69		
AS837B	49, 71	CTM20B	52	VS364TP124	69		
AS85008	33	CTM50B	52	VS364TTP124	69		
AS85008K	33	D38301	11	VS365RC124	69		
AS85010	33	D38302	11	VS365TC124	69		
AS85010K	33	E-1004	3	VS365TP124	69		
AS85012	33	E-1200	3	VS365TTP124	69		
AS85012K	33	E-1203	3				
AS85014	33	ENF1004K	3				
AS85014K	33	ENF1200K	3				
AS85016	33	ENF1203K	3				
AS85016K	33	ENF500K	3				
AS860	65	ENK127	5				
AS860608	30	ENK136	5				
AS860610	30	ENS015	4				

Notes



Amsino Medical Group Co Ltd.

708 Corporate Center Drive
Pomona, CA 91768
USA

P: +1 800.632.6746

info@amsino.com

www.amsino.com

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required. CE mark in the catalog is only based on latest registration status. The CE status may be changed after the latest revision of the catalog.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO, ACCU-MEASURE, AMSafe, AMSafe-3, AMShield, AMSure, AMSure HG, EZE-VAC, iRECEPTAL, MedXL, One Mask, Praxiject, Puggle, R-Polished, Receptal, Safgard, Sal-Sense, SENTINEL, SteriCadi, Sure-Lok, THE-POLE-SYRINGE and VAC-GARD are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2025, Amsino Medical Group Co. Ltd. All rights reserved.