



*AMSafe*<sup>®</sup>

Gravity IV Sets for  
Surgery Centers

A cost effective alternative  
without sacrificing quality  
or reliability.

For over 25 years, Amsino  
has provided high value  
intravenous products to  
the US and global markets.

***AMSINO***

Advancing Healthcare Worldwide<sup>®</sup>

## AMSafe® Gravity IV Sets for Surgery Centers

Amsino's **AMSafe®** brand of quality intravenous products offers healthcare providers a top value choice.

- Primary IV sets (15 and 60 drops/mL)
- Secondary IV set (15 drops/mL)
- Microbore and standard bore extension sets
- Needle-free PRN connectors

Item	Drops/mL	Length	Priming Volume	Qty/Case
<b>158309</b> Vented/Non-Vented Basic Set, Roller Clamp, 1 <b>AMSafe®</b> Needle-Free Y Site 6" from Distal End, Rotating Male Luer Lock, PE Poly Pouch	15	83"	17mL	50
<b>1510501</b> Vented Basic Set, Roller Clamp, 2 <b>AMSafe®</b> Needle-Free Y Site, Rotating Male Luer Lock, PE Poly Pouch	15	105"	21mL	50
<b>AS3101</b> Vented/Non-Vented, Back Check Valve, 1 Slide Clamp, 2 <b>AMSafe®</b> Needle-Free Y Sites, Roller Clamp, Rotating Male Luer Lock, Tyvek Form Fill Seal Pkg.	15	110"	22mL	50
<b>AS3110A</b> Vented/Non-vented, Back Check Valve, 3 <b>AMSafe®</b> Needle-Free Y Sites, 1 Slide Clamp, Roller Clamp, Rotating Male Luer Lock, Paper Poly Pouch	15	110"	22mL	50
<b>IVA3160</b> Vented/Non-Vented, Back Check Valve, 3 <b>AMSafe®</b> Needle-Free Y Sites, 1 Slide Clamp, Roller Clamp, Rotating Male Luer Lock, Paper Poly Pouch	60	110"	22mL	50
<b>AS3515</b> Vented/Non-Vented Secondary Set, Roller Clamp, Rotating Male Luer Lock, Hanger, PE Poly Pouch	15	40"	10mL	50

Item	Length	Priming Volume	Qty/Case
<b>AE1106</b> Microbore Tubing, Female Luer Lock, 1 Slide Clamp, Rotating Male Luer Lock, Tyvek Poly Pouch	6"	0.1mL	50
<b>AE0106</b> Microbore tubing, 1 Removable <b>AMSafe®</b> Needle-Free Connector, 1 Slide Clamp, Rotating Male Luer Lock, Tyvek Poly Pouch	6"	0.1mL	50
<b>AE3107</b> Standard Bore Tubing, 1 <b>AMSafe®</b> Needle-Free Connector, 1 Slide Clamp, Rotating Male Luer Lock, Tyvek Poly Pouch	7"	1.2mL	50
<b>003202</b> Standard Bore Tubing, Female Luer Lock, 2 Slide Clamps, 2 <b>AMSafe®</b> Needle-Free Y Sites, Rotating Male Luer Lock, Tyvek Poly Pouch	27"	4.9mL	50
<b>000606</b> Standard Bore Tubing, Female Luer Lock, 1 Pre-Pierced Y Injection Site, 1 <b>AMSafe®</b> Needle-Free Y Site, Rotating Male Luer Lock, PE Poly Pouch	6"	1.5mL	50
<b>AY0200</b> <b>AMSafe®</b> Needle-Free PRN Connector	100/Case		

## Features and Benefits

- Cost-effective alternative without sacrificing safety or quality
- 2 Needle-free Y-site options to meet all needs
- Easy-open packaging designed for compact storage and easy product access
- IV expertise: Amsino is the #1 provider of IV administration sets in Emergency Medical Services



Amsino is a leading medical device company dedicated to advancing healthcare worldwide by developing, manufacturing and marketing medical products that improve the safety and effectiveness of patient care.

Call Customer Service at **800-632-6746**

Visit our online catalog at **[www.Amsino.com](http://www.Amsino.com)**

**AMSINO**

Advancing Healthcare Worldwide®

**Amsino International, Inc.**

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • [info@amsino.com](mailto:info@amsino.com)