

AMSafe-3.® Now you can stock one IV set for all drip rates.




Provides the flexibility to select 10, 15 or 60 drops/mL.

- Features AMSafe® Needle-Free Valve.
- Chamber spikes directly into bag.
- No need for clamps above and below chamber.
- Includes male rotating luer lock.
- Available with or without pre-attached extension set

Amsino was the first to introduce IV sets with both Luer-activated and pre-pierced Y-sites specifically designed for EMS. Our compact packaging is uniquely designed to open fast. And our Needle-Free valve avoids concerns with positive displacement valves, is swabbable, MRI-friendly, and not limited to 24-hour use.



D38301
AMSafe-3 I.V. Set, 10,15, 60 Drops, Roller Clamp, 1 Pre-Pierced Y Site, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock, 83"



D38302
D38301 with pre-attached extension set including additional Pre-Pierced Y Site, AMSafe® Needle-Free Y Site, Slide Clamp, Rotating Luer Lock, 91"

AMSINO
Advancing Healthcare Worldwide®

Amsino International, Inc.
855 Towne Center Drive, Pomona, CA 91767
1-800-632-6746
Fax: 909-626-3888
info@amsino.com

www.amsino.com

AMSINO

Advancing Healthcare Worldwide®





The Specialists in EMS IV Sets










AMSafe-3® Multi-Drop IV Sets – 10, 15 or 60 drops/mL

Latex-Free, DEHP-Free, Packaged in PE Poly Pouch, 50 each per case

Item Number	Description	Length	Priming Volume
D38301	AMSafe-3 IV Set, Roller Clamp, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock 	83"	24mL
D38302	AMSafe-3 IV Set, Roller Clamp, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, Pre-Attached Extension Set with 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock 	91"	24mL




AMSure® EMS IV Sets

Latex-Free, DEHP-Free, Packaged in PE Poly Pouch, 50 each per case

Item Number			Description	Length	Priming Volume
10 Drops/mL	15 Drops/mL	60 Drops/mL			
108306	158306	608306	Non-Vented, Roller Clamp, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock 	83"	17mL
108302	NA	NA	Vented/Non-vented, Back Check Valve, Roller Clamp, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock 	83"	17mL
108305	158305	608305	Non-Vented, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site with Pre-Attached Extension Set, Roberts Clamp, Rotating Luer Lock 	89"	16mL
NA	158307	NA	Non-Vented, Roller Clamp, 2 AMSafe® Needle-Free Y Sites, 2 Pre-Pierced Y Injection Sites with Pre-Attached Extension Set, 1 Slide Clamp, Rotating Luer Lock 	89"	16mL
107201	NA	607201	Non-Vented, Roller Clamp, 1 Y Injection Site, Rotating Luer Lock 	72"	15mL
NA	AA3101	NA	Vented/Non-Vented, 1 Slide Clamp, Roller Clamp, 1 Y Injection Site, Rotating Luer Lock 	78"	16mL
AFS102	NA	NA	Vented/Non-vented, Flow Controller, 1 Slide Clamp, 1 AMSafe® Needle-Free Y Site, Rotating Luer Lock 	100"	20mL

AMSure® EMS Extension Sets and Connector

Latex-Free, DEHP-Free, Packaged in PE Poly Pouch

Item Number	Description	Length	Priming Volume	Units/Case
AE3108	1 AMSafe® Needle-Free Connector, Roberts Clamp, Rotating Luer Lock, PE Poly Pouch. Priming volume: 1.2mL 	8"	1.2mL	100
000606	Female Luer Lock, 1 Pre-Pierced Y Injection Site, 1 AMSafe® Needle-Free Y Site, 1 Slide Clamp, Rotating Luer Lock, PE Poly Pouch. Priming volume: 1.5mL 	6"	1.5mL	50
AY0200	AMSafe® Needle-Free PRN Connector 	0.7"	0.01mL	100

AMSafe-3™

To order **AMSafe-3** samples,
email AMSafe-3@amsino.com
or call 1-800-632-6746.

www.amsino.com