

Introducing the new ENFit Connector

A safer way to provide enteral feeding

The new ENFit connector has a unique design that:

- Does not allow connectivity with any non-enteral connectors
- Requires a push turn motion to firmly connect with feeding tubes
- Promotes a safe enteral feeding connectivity reducing misconnections and improving patient safety



A new global design change for enteral connectors will be phased into use and applied to all enteral patients.

AMSINO

Advancing Healthcare Worldwide®

AMSure® ENTERAL FEEDING SYSTEMS

from **AMSINO**

AMSure® Enteral Feeding Administration Kits

Item	Description (For use with the Sentinel® Enteral Feeding Pump)
ENF1004K	Pump Spike Set with ENFit Connector & Transition Connector
ENF500K	500 mL Bag Pump Set with ENFit Connector & Transition Connector
ENF1200K	1200 mL Bag Pump Set with ENFit Connector & Transition Connector
ENF1203K	1200 mL Bag Gravity Set with ENFit Connector & Transition Connector

AMSure® Enteral Feeding Syringes

Item	Description
ENS016	60cc Thumb Control Ring Syringe with ENFit Tip
ENS116	60cc Flat Top Piston Syringe with ENFit Tip

AMSure® Enteral Feeding Irrigation Kits

Item	Description
ENK126	500cc Graduated Container, 60cc Thumb Control Ring Piston Syringe with ENFit Tip, patient I.D. label
ENK127	500cc Graduated Container, 60cc Thumb Control Ring Piston Syringe with ENFit Tip, patient I.D. label, resealable pole bag
ENK136	1000cc Outer Tray, 500cc Graduated Container, 60cc Thumb Control Ring Piston Syringe with ENFit Tip, tip protector, alcohol prep pad, large moisture-proof underpad, latex-free, sterile

To facilitate a smooth conversion to the new ENFit Connector Amsino has added a Transition Connector to the Administration Kits which works with traditional feeding tubes. Once the ENFit Feeding Tubes are fully established, the Transition Connector will be discontinued.

During this transition phase Amsino is collaborating with the Global Enteral Device Supplier Association (GEDSA), a nonprofit trade association formed to introduce international standards for medical device tubing connectors. As an associate member of GEDSA, Amsino's goal is to enhance patient safety by reducing the risk of misconnections.



Rx only.

About Stay Connected

The Stay Connected initiative for using safer connectors is a communications program created to facilitate a successful transition from medical device applications that allow connection between unrelated delivery systems to new, safer connectors that ensure compatibility and consistency while reducing the likelihood of tubing misconnections. Visit www.stayconnected.org.

Your Enteral Feeding Partners

AMSINO

Advancing Healthcare Worldwide®

Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 • 800.632.6746 • fax: 909.626.3888 • info@amsino.com • Online catalog at www.AmsinoCatalog.com

All trademarks are the property of their respective owners.
© Copyright 2016, Amsino International, Inc. All rights reserved.

EN-US-2017-06-09C