

# EZE-VAC<sup>®</sup> 2Go<sub>≡</sub>

## Goes Anywhere, Safe Everywhere

NEW!

### Safe and convenient fluid collection in every healthcare setting.


The EZE-VAC 2Go 800 mL disposable suction canister is universally designed to fit portable suction machine holders and restricted space environments. This single-use canister is engineered for convenient, safe, and reliable waste fluid collection.

- **Integrated mechanical float valve** automatically stops suction when full, helping to minimize the risk of fluid overflow and contamination.
- **Clearly marked patient and vacuum ports** on lid helps prevent misconnection.
- **Conveniently designed** with pre-attached caps on lid patient port and pour spout to support caregiver safety.
- **Canisters are nestable** for compact storage.



### Ordering Information

Item Code	Description	Qty/Unit
43208-01	800 mL Disposable Suction Canister	48/Case



### Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • info@amsino.com

For the full Amsino product line, visit us at [www.amsino.com](http://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. EZE-VAC is a registered trademarks of Amsino International, Inc. © Copyright 2018, Amsino International, Inc. All rights reserved.