

EZE-VAC[®] 2.0 Canisters

NEW

The **EZE-VAC 2.0** Suction Canister System is the second generation of the well-known and reliable EZE-VAC Suction Canister System. The new EZE-VAC 2.0 is a single-use disposable canister that does not need a secondary disposable liner.

This new canister system is very economical and is compatible with most hardware in the market. As with many other of Amsino's Waste Fluid Management Systems, EZE-VAC 2.0 also includes an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

EZE-VAC 2.0 Canisters support best practices by integrating VAC-GARD[®], a non-mechanical shut-off valve that prevents back flow of aerosols and fluids into the vacuum system and includes a 0.3 micron bacterial eliminating filter effective in removing 99.9% of pseudomonas.



- **Lid attaches firmly** to the canister, maximizing vacuum power and reducing risk of leakage.
- **Easy-to-read volume marks** help with patient documentation.
- **Clearly marked** patient and vacuum ports on lid helps prevent misconnection.
- **Canisters are nestable** for compact storage.
- **Multiple volume sizes** to match your needs.
- **Compatible with most hardware** currently being used in the market; this allows for ease of conversion.

Ordering information

Item code	Description	Qty/Unit
43210-01	800 mL EZE-VAC 2.0 Canister and Lid	72/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	60/Case
43216-01	2000 mL EZE-VAC 2.0 Canister and Lid	40/Case
43215-01	3300 mL EZE-VAC 2.0 Canister and Lid	40/Case


Easy-to-read graduated volume marks help with patient documentation

Generous patient label for documentation



Clearly marked patient and vacuum ports

VAC-GARD[®] non-mechanical shut-off valve



Pour Spout

EZE-VAC® 2.0 800ml Canister

NEW

Goes Anywhere. Safe Everywhere.

Safe and convenient fluid collection in every healthcare setting.

The EZE-VAC 2.0 800 mL disposable suction canister is universally designed to fit portable suction machine holders and restricted space environments. This single-use canister is engineered for convenient, safe, and reliable waste fluid collection.

- **Integrated VAC-GARD®** automatically stops suction when full, helping to minimize the risk of fluid overflow and contamination.
- **Clearly marked patient and vacuum ports** on lid help to prevent misconnection.
- **Conveniently designed** with easy-to-read volume marks.
- **Canisters are nestable** for compact storage.



Ordering Information

Item code	Description	Qty/Unit
43210-01	800 mL EZE-VAC 2.0 Canister and lid	72/Case



Ideal for portable suction pumps

AMSINO
Advancing Healthcare Worldwide®

Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768
800.632.6746 • fax: 909.626.3888 • info@amsino.com

For the full Amsino product line, visit us at www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. EZE-VAC is a registered trademarks of Amsino International, Inc. ©Copyright 2018, Amsino International, Inc. All rights reserved.

WFM-GLO-2020-02-08A