

ENgineered with Love

Peace-of-mind with thought-through, intuitive design and reliable performance.



Item Code	Description	Dimensions	Weight	Flow Rate Range	Volume Range	Accuracy*
PG001	PUGGLE™ Enteral Feeding Pump – AC Adapter and Pole Clamp Included	5.4"H x 5.0"W x 2.5"D	1.43 lbs.	1 to 400 mL/h in 1 mL increments	1 to 9,999 mL in 1 mL increments	± 5%

T100

Item Code	Description	Qty/Unit
PE90ST	PUGGLE™ Enteral Pump Spike Set, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B5T	PUGGLE™ Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B12T	PUGGLE™ Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case



Amsino International, Inc.
 708 Corporate Center Drive, Pomona, CA 91768 • Amsino.com
 Phone: 800.632.6746 • Fax: 909.626.3888
 Email: info@amsino.com

800.632.6746



Please refer to the Instructions for Use for additional information.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. ENFit® is a registered trademark of GEDSA. The symbol ® indicates the trademark is registered in the

© Copyright 2021, Amsino International, Inc. All rights reserved.

EN-US-2021-4-30-DS