

NEW!



High-Profile Gastrostomy Feeding Tube

- 3-Port feeding tube with clearly marked feeding, medication and balloon ports for easy administration
- Wavy retention ring with a smooth surface to promote patient comfort and provide ventilation of the stoma site
- Inflatable silicone internal retention balloon for secure placement
- Graduated centimeter markings for tube position verification
- Radiopaque stripe to aid in g-tube visualization during medical imaging

Item Code	French Size	Balloon Size	Qty/Unit
PLG12-10	12 Fr	10 mL	20/Case
PLG14-10	14 Fr	10 mL	20/Case
PLG14-20	14 Fr	20 mL	20/Case
PLG16-10	16 Fr	10 mL	20/Case
PLG16-20	16 Fr	20 mL	20/Case
PLG18-20	18 Fr	20 mL	20/Case
PLG20-10	20 Fr	10 mL	20/Case
PLG20-20	20 Fr	20 mL	20/Case
PLG22-20	22 Fr	20 mL	20/Case
PLG24-20	24 Fr	20 mL	20/Case



For more information, contact Amsino Customer Service at
1.800.632.6746 or cssupport@amsino.com



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768
 800.632.6746 • fax: 909.626.3888 • info@amsino.com • www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use / Operator's Manual for a complete listing of the indications, contraindications, warnings, and precautions. The symbol * indicates that the trademark is registered in the U.S. Patent and Trademark Office. AMSure® and PUGGLE™ are registered trademarks of Amsino International, Inc. © 2018, Amsino International, Inc. All rights reserved.