

# STERiCADI™

## Suction Device Holster System

### A preferred infection control safeguard

SteriCadi™ is a disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Secretion reservoir with safety tip promotes suction device drainage
- Clip design with tubing clasp provides a secure attachment

**Kits are available with a variety of suction yankauer devices**



Vented and Non Vented




Bulbous and Beveled



# SteriCadi™ Suction Device Holster System

Designed to create convenient, consistent placement for suction device storage to help support better infection control measures for patients and caregivers.



## Ordering Information

Item Code	Description	Units/Case
SIH0001	SteriCadi™ Suction Device Holder	30

## SteriCadi™ Kits

Item Code	Description	Units/Case
SIH0002	SteriCadi™ with Regular Beveled Tip, Non-vented, Rigid Suction Yankauer Device, Sterile	30
SIH0003	SteriCadi™ with Regular Beveled Tip, Vented, Rigid Suction Yankauer Device, Sterile	30
SIH0004	SteriCadi™ with Bulbous Tip, Non-vented, Rigid Suction Yankauer Device, Sterile	30
SIH0005	SteriCadi™ with Bulbous Tip, Vented, Rigid Suction Yankauer Device, Sterile	30



Advancing Healthcare Worldwide®

### Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768  
800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use / Operator's Manual for a complete listing of the indications, contraindications, warnings, and precautions. The symbol ® indicates that the trademark is registered in the U.S. Patent and Trademark Office. SteriCadi is a trademark of Amsino International, Inc. © 2016, Amsino International, Inc. All rights reserved.