

AMSure® Enteral Feeding Systems. Safety through simplicity.



ALCOR® AMSure® Enteral Feeding Administration Kits*

Item	Description
E-0500M	500 mL Bag Pump Set
E-1200	1200 mL Bag Pump Set
E-1200M	1200 mL Bag Pump Set includes Magnet
E-1203	1200 mL Gravity Set
E-1004	Pump Spike Set

Purple Tips with IV Warning Tag

Helps avoid hazardous connection to IV



*Use with ALCOR® Sentinel® Enteral Feeding Pump

ALCOR® Sentinel® Enteral Feeding Pump

Item	Description
S-1000-SI	ALCOR® Sentinel® Enteral Feeding Pump

Pump Safety Features & Benefits

- Large, Dual Display** Shows infusion rate & volume delivered simultaneously from a distance
- Patented Drop Detection** Ensures reliable delivery while helping to minimize false alarms
- Advanced Memory** Retains pump history until manually erased, even in power failures
- Audible & Visual Alarms** Volume adjustable only on non-vital alarms; prominent & clear visual alarms
- Simple 3-Step Operation**



1. Press ON 2. Set infusion rate 3. Begin feeding

Amsino Customer Care Center: 800.632.6746

AMSINO

Advancing Healthcare Worldwide®

AMSure® Enteral Feeding Systems with the new ENFit Connector.

A safer way to provide enteral feeding.



ENF1200K

ALCOR® AMSure® Enteral Feeding Administration Kits*

Item	Description
ENF1004K	Pump Spike Set with ENFit Connector & Transition Connector
ENF500K	500 mL Bag Pump Set with ENFit Connector & Transition Connector
ENF1200K	1200 mL Bag Pump Set with ENFit Connector & Transition Connector
ENF1203K	1200 mL Bag Gravity Set with ENFit Connector & Transition Connector

A new global design change for enteral connectors will be phased into use and applied to all enteral patients.

To facilitate a smooth conversion to the new ENFit Connector, Amsino has added a Transition Connector to the Administration Kits which works with traditional feeding tubes. Once the ENFit Feeding Tubes are fully established, the Transition Connector will be discontinued.

The new ENFit connector has a unique design that promotes a safe enteral feeding connectivity reducing misconnections and improving patient safety.



*Use with ALCOR® Sentinel® Enteral Feeding Pump

Enteral Syringes



ENS015
60 cc Piston Irrigation Syringe, Thumb Control Ring, ENFit Tip, Poly Pouch, Sterile



ENS115
60 cc Piston Irrigation Syringe, Flat Top, ENFit Tip, Poly Pouch, Sterile



AS016
60 cc Piston Irrigation Syringe, Thumb Control Ring, Luer Tip Adapter, Resealable IV Pole Bag, Non-sterile



AS116
60 cc Piston Irrigation Syringe, Flat Top, Luer Tip Adapter, Resealable IV Pole Bag, Non-sterile



Advancing Healthcare Worldwide®

Your Enteral Feeding Partner