

AMSure[®]
UROLOGICAL CARE



Comprehensive and Cost-Effective Urological
Solutions Engineered to Achieve Superior Results

SELF-CATHETERIZATION KITS

URETHRAL SELF-CATHETERIZATION KIT WITH R-POLISHED® EYES

AMSure® R-Polished technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort. The AMSure Self-Catheterization kit allows the patient, even with limited dexterity, to perform incontinence care with confidence.

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, latex-free and BPA-free

Kits include:

- Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter
- 1100 mL collection bag with insertion tip and protector cap
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Waste disposal bag
- Gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

HPCS code A4353

ITEM CODE	FRENCH SIZE	QTY/UNIT
AS85008	8 Fr	50/Case
AS85010	10 Fr	50/Case
AS85012	12 Fr	50/Case
AS85014	14 Fr	50/Case
AS85016	16 Fr	50/Case



URETHRAL CATHETERS

RED RUBBER LATEX URETHRAL CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Sterile, DEHP-free and BPA-free



HPCS code A4351

ITEM CODE	DESCRIPTION	FEATURES	QTY/UNIT
AS44008	8 Fr, 11 in. (27.94 cm)	Female	100/Case
AS44010	10 Fr, 11 in. (27.94 cm)	Female	100/Case
AS44012	12 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44014	14 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44016	16 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44018	18 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44020	20 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44022	22 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44024	24 Fr, 16 in. (40.64 cm)	Male	100/Case
AS44028	28 Fr, 16 in. (40.64 cm)	Male	100/Case

PVC INTERMITTENT URETHRAL CATHETERS

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Sterile, latex-free and BPA-free



HPCS code A4351

ITEM CODE	DESCRIPTION	FEATURES	QTY/UNIT
AS860608	8 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860610	10 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860612	12 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860614	14 Fr, 6 in. (15.24 cm)	Female	50/Case
AS860616	16 Fr, 6 in. (15.24 cm)	Female	50/Case
AS861006	6 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861008	8 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861010	10 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS861608	8 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861610	10 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861612	12 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861614	14 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861616	16 Fr, 16 in. (40.64 cm)	Male	50/Case
AS861618	18 Fr, 16 in. (40.64 cm)	Male	50/Case



PVC INTERMITTENT URETHRAL CATHETER WITH R-POLISHED EYES

- Advanced R-Polished technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- Sterile, latex-free and BPA-free

HPCS code A4351



ITEM CODE	DESCRIPTION	FEATURES	QTY/UNIT
AS960610	10 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960612	12 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960614	14 Fr, 6 in. (15.24 cm)	Female	50/Case
AS960616	16 Fr, 6 in. (15.24 cm)	Female	50/Case
AS961010	10 Fr, 10 in. (25.40 cm)	Pediatric	50/Case
AS961610	10 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961612	12 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961614	14 Fr, 16 in. (40.64 cm)	Male	50/Case
AS961616	16 Fr, 16 in. (40.64 cm)	Male	50/Case

URETHRAL INSERTION TRAYS AND KITS

URETHRAL INSERTION TRAY

- Trays provide clinician flexibility to choose their preferred drainage or leg bag

PVC CATHETER

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Sterile, latex-free and BPA-free

RED RUBBER CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Sterile, DEHP-free and BPA-free

Trays include:

- 1000 mL outer basin tray
- 16 in. urethral catheter
- Waterproof and fenestrated drapes
- Vinyl powder-free gloves
- (3) 10% povidone-iodine swab sticks
- Lubricating jelly
- 4 oz. specimen container

HPCS code A4353



ITEM CODE	URETHRAL CATHETER TYPE	QTY/UNIT
AS870	14 Fr PVC (Latex-free)	20/Case
AS87114	14 Fr Red Rubber Latex	20/Case
AS87118	18 Fr Red Rubber Latex	20/Case



PRE-CONNECTED DRAINAGE BAG URETHRAL CATHETER KIT

- Sterile kits provide clinician flexibility to choose their preferred drainage or leg bag
- Gloves and povidone-iodine swab sticks help maintain sterility and reduce the risk of infection

Kits include:

- 1200 mL low-profile urinary drainage bag
- 16 in. (40.64 cm) urethral catheter
- CSR wrap, waterproof and fenestrated drapes
- (2) Nitrile powder-free gloves
- (3) Povidone-iodine swab sticks
- Lubricating jelly
- 2 oz. specimen container with ID label

HPCS code A4353



ITEM CODE	14 FR RED RUBBER CATHETER	14 FR PVC CATHETER	QTY/UNIT
AS89214R*	●		20/Case
AS89214		●	20/Case

*Contains latex

FOLEY CATHETERS

2-WAY SILICONIZED LATEX FOLEY CATHETERS

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage and aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free and BPA-free



HPCS code A4338



ITEM CODE	DESCRIPTION	FEATURES	QTY/UNIT
AS41012	12 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41014	14 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41016	16 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41018	18 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41020	20 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41022	22 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41024	24 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41026	26 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41028	28 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41030	30 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS42012	12 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42014	14 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42016	16 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42018	18 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42020	20 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42022	22 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42024	24 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42026	26 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42028	28 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42030	30 Fr, 16 in. (40.64 cm)	30 cc	10/Box

2-WAY 100% SILICONE FOLEY CATHETERS

- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material reduces encrustation and calcifications, making it ideal for longer-term catheterization
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, latex-free, DEHP-free and BPA-free



HPCS code A4344

ITEM CODE	DESCRIPTION	FEATURES	QTY/UNIT
AS41012S	12 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41014S	14 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41016S	16 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41018S	18 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41020S	20 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41022S	22 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41024S	24 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41026S	26 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS41028S	28 Fr, 16 in. (40.64 cm)	5 cc	10/Box
AS42016S	16 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42018S	18 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42020S	20 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42022S	22 Fr, 16 in. (40.64 cm)	30 cc	10/Box
AS42024S	24 Fr, 16 in. (40.64 cm)	30 cc	10/Box

FOLEY INSERTION TRAYS

FOLEY INSERTION TRAY

- Deep tray basin designed to prevent contamination
- Gloves and povidone-iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, latex-free, DEHP-free and BPA-free

Foley Insertion Trays include:

- 1000 mL collection tray
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Lubrication jelly
- Fenestrated drape
- Pre-filled syringe with sterile water

HPCS code A4310



ITEM CODE	10 cc SYRINGE	30 cc SYRINGE	QTY/UNIT
AS880	●		20/Case
AS890		●	20/Case

ADD-A-FOLEY CATHETER TRAY

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone-iodine solution help maintain sterility and reduce the risk of infection
- Sterile and latex-free

Add-A-Foley Trays include:

- 2000 mL low-profile urinary drainage bag
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone-iodine solution
- Forceps
- Lubricating jelly
- 2 oz. specimen container with patient label
- Central supply room wrap and tyvek header bag
- Pre-filled syringe with sterile water

HPCS code A4354



ITEM CODE	10 cc SYRINGE	30 cc SYRINGE	QTY/UNIT
AS89300	●		10/Case
AS89330		●	10/Case

FOLEY INSERTION TRAYS

PRE-CONNECTED DRAINAGE BAG FOLEY CATHETER TRAYS

- Sterile tray includes all necessary components for aseptic catheter insertion
- Tray improves clinical efficacy through cost-effective inventory management
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Tray is available with both siliconized latex and 100% silicone catheters to suit patient needs

Pre-connected drainage bag Foley tray contents include:

- 2000 mL low-profile urinary drainage bag
- Pre-connected 16 in. Foley catheter with 5 cc balloon
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone-iodine solution
- Forceps
- Lubricating jelly
- Syringe prefilled with 10 cc sterile water
- 2 oz. specimen container with a patient label
- Central supply room wrap and tyvek header bag



HPCS code A4314

ITEM CODE	FOLEY CATHETER FR SIZE			CATHETER MATERIAL		QTY/UNIT
	14 Fr	16 Fr	18 Fr	Siliconized Latex	100% Silicone (Latex-free)	
AS89314	•			•		10/Case
AS89316		•		•		10/Case
AS89318			•	•		10/Case
AS89314S	•				•	10/Case
AS89316S		•			•	10/Case
AS89318S			•		•	10/Case

SYRINGES, IRRIGATION TRAYS AND KITS

BULB SYRINGE

- Designed to stand on the bulb end to help minimize contamination
- Tip protector included


THUMB CONTROL RING SYRINGE

- Re-sealable pole bag designed for safe and convenient storage
- Includes syringe luer tip adapter and protector cap

FLAT TOP SYRINGE

- Visible dosing marks designed for ease of use
- Includes syringe luer tip adapter and protector cap

HPCS code A4322



ITEM CODE	60 mL SYRINGE TYPE			PACKAGE TYPE				STERILE	QTY/UNIT
	Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Pack	Resealable Pole Bag	Bulk		
AS116			•			•			30/Case
AS016		•				•			30/Case
AS010	•						•		50/Case
AS011	•			•				•	50/Case
AS011P	•				•			•	50/Case
AS015		•		•				•	50/Case
AS115			•	•				•	50/Case

Note: All syringes are latex-free, DEHP-free and BPA-free

IRRIGATION TRAYS AND KITS

Non-sterile Kits include:

- 500 mL graduated container
- Patient ID label and tip protector

HPCS code A4320



ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

Sterile Trays include:

- 1000 mL outer tray
- 500 mL graduated container
- Waterproof underpad
- Alcohol prep pad
- Patient ID label and tip protector

HPCS code A4320



ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130		•	20/Case
AS136	•		20/Case

URINARY COLLECTION DEVICES

URINARY DRAINAGE BAG

- Anti-reflux mechanism helps improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, latex-free and BPA-free



AS322 Teardrop



AS32600 Low-profile

HPCS code A4358

ITEM CODE	BAG SHAPE		ANTI-REFLUX TYPE			NEEDLE-FREE SAMPLING PORT		DRAIN TYPE		HANGER TYPE			QTY/UNIT
	Low Profile	Teardrop	Chamber	Flutter Valve	Device	Pre-pierced Blunt Cannula	Male Luer Lock or Luer Slip	T-tap	Metal Clamp	Single Hook	Double Hook	Rope	
AS312		•	•			•		•			•	•	20/Case
AS322		•	•			•		•			•	•	20/Case
AS32200	•		•			•		•		•		•	20/Case
AS32600*	•		•				•	•		•		•	20/Case
AS332		•		•		•		•		•			20/Case
AS326	•				•	•			•		•	•	20/Case
AS302		•			•	•		•			•	•	20/Case

*100% sterile

PEDIATRIC URINE COLLECTOR

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, latex-free, DEHP-free and BPA-free

HPCS code A4328

ITEM CODE	DESCRIPTION	QTY/UNIT
AS409	200 mL, Pediatric Urine Bag	50/Box



AS409

URINARY COLLECTION DEVICES

URINARY LEG BAGS

- Soft and breathable pre-attached leg straps provide patient comfort
- Anti-reflux device helps improve patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway, latex-free and BPA-free




AS307N

HPCS code A4358

ITEM CODE	600 mL LEG BAG	900 mL LEG BAG	PUSH-PULL DRAIN PORT	TWIST-TURN DRAIN PORT	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER AND PROTECTIVE CAP	QTY/UNIT
AS306N	•		•			48/Case
AS307N		•	•			48/Case
AS308N	•			•		48/Case
AS309N		•		•		48/Case
LBE306N	•		•		•	48/Case
LBE307N		•	•		•	48/Case
LBE308N	•			•	•	48/Case
LBE309N		•		•	•	48/Case

EXTENSION TUBE FOR LEG BAGS

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, latex-free and BPA-free




AS311

HPCS code A4331

SPECIMEN CONTAINERS

- 4 oz. leak resistant specimen container
- Latex-free, DEHP-free, BPA-free

ITEM CODE	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	QTY/UNIT
AS340		•		•		•	100/Case
AS341			•	•	•		100/Case
AS343			•				500/Case
AS345	•			•*			100/Case



AS340

*100% sterile

Your Local Distributor



Amsino International, Inc.

708 Corporate Center Drive
Pomona, CA 91768

P: 800.632.6746

F: 909.626.3888

info@amsino.com

www.amsino.com



Advancing Healthcare Worldwide®