

AMSafe® IV Hydration




Amsino manufactures both IV administration sets and pre-filled flush syringes to support various IV hydration treatments.

Gravity Administration IV Sets:


Amsino has a comprehensive line of gravity administration IV sets that are designed to improve safety and reduce infection for both caregivers and patients.

- Available with needle-free, pre-pierced and injection Y-sites in an array of configurations
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex


General Primary Sets

REF	DROPS/mL	LENGTH	DESCRIPTION	APPROXIMATE PRIMING VOL	QTY/UNIT
108301 ©	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108306 ©	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
AA3101 ©	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					

Piggyback Set with Back Check Valve

REF	DROPS/mL	LENGTH	DESCRIPTION	APPROXIMATE PRIMING VOL	QTY/UNIT
AS3101 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					

Standard Extension Set

REF	LENGTH	DESCRIPTION	APPROXIMATE PRIMING VOL	QTY/UNIT
AE3108 ©	8 in. (20 cm)	Removable Sure-Lok Needle-free Connector, Pinch Clamp	1.76 mL	100/Case
				

Microbore Extension Set

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AE0106 ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				

AMSafe® Pre-Filled Flush Syringe

0.9% Sodium Chloride Injection, USP

Features

- DEHP-free, PVC-free, not manufactured with natural rubber latex
- Sterile, non-pyrogenic
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared



10 mL Syringes: Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1010 ©	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1010TM	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1210 ©	10 mL USP Sodium chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)

Also available in 3 mL, 5 mL and Sterile configurations.

See our catalog for more, [click here.](#)

To view our complete IV product line, visit us at [Amsino.com](https://www.amsino.com), **[click here.](#)**



Advancing Healthcare Worldwide®

Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA

+1 800.632.6746 • fax: +1 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSafe is a registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2023, Amsino Medical Group Co. Ltd. All rights reserved.