

IV Pump Administration Sets

For use with Baxter Spectrum, Flo-Gard® 6201 and 6301 series pumps.

Features

- 510K approved IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and Pre-Pierced Y-injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex



Sure-Lok® IV Pump Administration Sets

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
FP092	10	92 in. (234 cm)	Roller Clamp, Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
FP092Y1	10	92 in. (234 cm)	1 Needle-Free Valve, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
FP169BY3	60	109 in. (277 cm)	Back check Valve, 3 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case

Sure-Lok® IV Pump Administration Sets with Filter


Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 Needle-Free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
FP100FY1	10	100 in. (254 cm)	1 Needle-Free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
FP100F12Y1	10	100 in. (254 cm)	Sure-lok Needle-Free Valve, Slide and Roller clamp, 1.2 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case

For use with Baxter Spectrum, Flo-Gard® 6201 and 6301 series pumps.




Pre-Pierced Y-Injection Site IV Pump Administration Sets

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
FP100BI2	10	100 in. (254 cm)	Back check Valve, 2 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case




FP109BI3	10	109 in. (277 cm)	Back check Valve, 3 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
----------	----	------------------	--	-------	---------



Pre-Pierced Y-Injection Site IV Pump Administration Set with Filter

Item Code	Drops/ mL	Length	Description	Approx Priming Volume	Qty/Unit
FP106FI2	10	106 in. (269 cm)	Back check Valve, 2 Pre-Pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case



To view our complete product line, visit us at www.amsino.com



Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA
 +1 800.632.6746 • fax: +1 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Flo-Gard, Sure-Lok and TUBETECH are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2024, Amsino Medical Group Co. Ltd. All rights reserved.