

# AMSure® ENTERAL FEEDING SYSTEMS



**AMSINO**  
Advancing Healthcare Worldwide®

# AMSure®

## ENTERAL FEEDING SYSTEMS



## Optimizing Nutritional Delivery

AMSure enteral feeding systems enhance delivery from nutrition source to patient feeding tube and are designed to provide best patient outcomes. The systems include both stationary and ambulatory options.

AMSure enteral systems include the Sentinel® stationary and PUGGLE® ambulatory enteral feeding pumps, enteral gravity bag feeding sets, pump feeding sets, feeding syringes and irrigation kits and trays.


The comprehensive Enteral Therapy Portfolio is available with standard and ENFit® connections. AMSure is compliant with the new ENFit ISO 80369-3 design standard and supports the goal to improve patient safety by reducing the risk of misconnections.

# AMBULATORY ENTERAL PUMP

## PUGGLE® Enteral Feeding Pump

The versatile pump is designed for pediatric and adult ambulatory applications. Simplified program settings make feeding on-the-go a snap.

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours



PE90B12T



DESCRIPTION	DIMENSIONS	WEIGHT	OPERATING RANGE	DELIVERED VOLUME	ACCURACY*
PUGGLE Enteral Feeding Pump, AC Adapter and Pole Clamp	5.4 x 5.0 x 2.5 inches	1.43 lbs.	1 to 400 mL/h in 1 mL increments	1 – 9999 mL in 1 mL increments	+/- 5%

REF	DESCRIPTION	QTY/UNIT
<b>PG001</b>	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
<b>PE90ST</b>	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
<b>PE90B5T</b>	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
<b>PE90B12T</b>	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case

\*Accuracy of feed at the specified rate over one hour.

# STATIONARY ENTERAL PUMP

## Sentinel® Enteral Feeding Pump

The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each



# ENFit® and Non-ENFit Enteral Feeding Gravity & Pump Sets

AMSure Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel® enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free, not manufactured with natural rubber latex



ENF1200K



ENF500K

REF	ENFit CONNECTOR	STANDARD BORE (DEHP FREE)	LARGE BORE (DEHP FREE)	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
				Spike Only No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/EntraFlo™	
E-1004				•				•	•	•	30/Case
E-1200						•		•	•	•	
E-1203						•	•				
ENF1004K ©	•			•				•	•	•	
ENF1200K ©	•					•		•	•	•	
ENF1203K ©	•					•	•				
ENF1203 KD	•	•				•	•				
ENF1203 KLD	•		•			•	•				
ENF500K ©	•				•			•	•	•	

EntraFlo™ is a registered trademark of Medline Industries, Inc.  
 Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.  
 ENFit® is a registered trademark of GEDSA.

# ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

## Bulb Syringe

- Flexible soft bulb promotes one-handed administration



## Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage



## Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010 ©			•						•		50/Case
AS011 ©			•				•			•	50/Case
AS115 ©		•			•		•			•	50/Case
AS116		•			•			•			30/Case
ENS015 ©	•			•			•			•	50/Case
ENS015NS	•			•				•			30/Case
ENS115 ©	•				•		•			•	50/Case
ENS115NS	•				•			•			30/Case
ENS115T*	•				•			•			30/Case

\*Includes Transition Adapter

# ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



## Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	30/Case
AS127		•		•		•	•	30/Case
ENK127 ©	•	•		•		•	•	30/Case

## Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	20/Case
ENK136 ©	•	•	•		•	•	•	•	20/Case

# Amsino offers a full line of enteral products including both Stationary and Ambulatory Enteral Feeding Pumps, pump and gravity sets, enteral syringes, syringe kits and trays.

To see our complete line of Enteral products, [click here](#).



## Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA  
+1 800.632.6746 • fax: +1 909.626.3888 • [info@amsino.com](mailto:info@amsino.com)

To view our complete product line, visit us at [Amsino.com](https://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure, Puggle and Sentinel are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2025, Amsino Medical Group Co. Ltd. All rights reserved.