

AMSure®

Enteral Syringes & Kits with ENFit® Tip



ENFit® Tip

AMSure® Enteral Syringes are compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.

AMSure® Thumb Controlled Piston Syringe with ENFit® Tip

AMSure® Thumb Controlled Piston Syringe with ENFit® tip is designed for maximum control allowing for a one-handed administration.



AMSure® Flat Top Piston Syringe with ENFit® Tip

AMSure® Flat Top Piston Syringe with ENFit® tip has visible dosing marks designed for ease of use. Popular choice to deliver viscous fluid with ease.










AMSure® Enteral Feeding Irrigation Kits with ENFit® Syringe

AMSure® Enteral Feeding Irrigation Kit is designed for simplicity and efficiency of care. The non-sterile kit includes a 60 mL thumb controlled piston syringe with ENFit® tip, patient label, 500 mL graduated container. Packaged in a poly pouch or resealable bag.



Order Information

Item Code	Image	Description	Features	Qty/Unit
ENS015		AMSure® 60 mL Thumb Controlled Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS015NS <i>ENS016 equivalent</i>		AMSure® 60 mL Thumb Controlled Piston Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENS115T		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip, transition adapter included	Non-Sterile Resealable Pole Bag*	30/Case
ENS115		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS115NS <i>ENS116 equivalent</i>		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENK126		AMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Controlled Piston Syringe with ENFit® Tip, Patient Label	Sterile Poly Pouch	30/Case
ENK127		AMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Controlled Piston Syringe with ENFit® Tip, Patient Label	Non-Sterile Resealable Pole Bag*	30/Case

*Resealable pole bag designed for safe and convenient storage.

All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

To view our complete product line,
visit us at www.amsino.com

AMSINO

Advancing Healthcare Worldwide®

Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768

800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure® is a registered trademark of Amsino Medical Group Co. Ltd. © Copyright 2023, Amsino Medical Group Co. Ltd. All rights reserved.