



AMSure® PVC Intermittent Urethral Catheter with R))POLISHED™ Eyes

AMSure PVC Intermittent Catheters are designed for short-term, intermittent use and are firm enough for easy insertion, yet soft enough to be slightly pliable. The AMSure Intermittent Catheters utilize Amsino's proprietary radio frequency R-Polish technology to create ultra-smooth eyelets and a smooth surface for optimal patient comfort. They are available in various French sizes as well as female, pediatric and male lengths. **(HCPCS code A4351)**

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Clear shaft provides easy viewing as you measure
- PVC material eliminates latex allergy concerns
- Designed to minimize trauma to help reduce risk of infections
- Individually packaged, sterile, BPA-free, not manufactured with natural rubber latex



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	AS860608	8 Fr, 6 in. (15.24 cm)	50/Case
	AS860610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS860612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS860614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS860616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS861006	6 Fr, 10 in. (25.40 cm)	50/Case
	AS861008	8 Fr, 10 in. (25.40 cm)	50/Case
	AS861010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS861608	8 Fr, 16 in. (40.64 cm)	50/Case
	AS861610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS861612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS861614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS861616	16 Fr, 16 in. (40.64 cm)	50/Case
	AS861618	18 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	AS861610C	10 Fr, 16" (40.65 cm)	30/Case
	AS861612C	12 Fr, 16" (40.65 cm)	30/Case
	AS861614C	14 Fr, 16" (40.65 cm)	30/Case
	AS861616C	16 Fr, 16" (40.65 cm)	30/Case
	AS861618C	18 Fr, 16" (40.65 cm)	30/Case