

AMSINO GENERAL CATALOG

AMSINO

Advancing Healthcare Worldwide®





AMSure® Enteral Feeding Systems 1	AMSure® Sterile Water and Saline 39
PUGGLE® Enteral Feeding Pump 2	Inhalation Solutions 40
Sentinel™ Enteral Feeding Pump 3	Irrigation Solutions 41
ENFit® and Non-ENFit Feeding Sets 3	Nebulizer Kits 73
ENFit and Non-ENFit Enteral Syringes 4	Device and Wound Irrigation 42
ENFit and Non-ENFit Irrigation Kits and Trays 5	Saline Products 43
AMSafe® Medication Delivery Systems 6	MedXL Saline Flush Syringes 45
Gravity Administration IV Sets 7	Prefilled Syringes 47
10 Drops/mL 7	PraxiMed Syringe Filling and Labeling System 49
15 Drops/mL 9	RECEPTAL® Fluid Management Systems 51
60 Drops/mL 11	EZE-VAC® Suction Canister System 52
AMSafe-3® Selectable Drip Chamber,	EZE-VAC 2.0 Suction Canister System 53
10, 15 & 60 drops/mL 12	RECEPTAL Liners 54
Secondary Piggyback Sets 13	RECEPTAL Twist Liners 55
IV Extension Sets 13	RECEPTAL Outer Shells 56
Standard Extension Sets 13	ACCU-MEASURE® Critical Measurement Devices 57
Microbore Extension Sets 14	RECEPTAL Solidifiers 57
Flow Controller Extension Sets 14	AMSure Suction Tubing with Connectors 58
Pressure Infusion Extension Sets 14	Suction Yankauers and Kits 59
Needle-Free Connectors 15	Pre-Connected Yankauer and Tubing 59
Vial Access Device 15	SteriCadi™ Suction Device Holster System 60
IV Start Kit 15	Accessories 61
Sure-Lok® Needle-Free IV Pump Administration Sets 16	Reusable Hardware 64
AMSure® Urological Care 17	iRECEPTAL® Surgical Waste Fluid Management System 67
Foley Catheter Trays & Collection 18	AMSure® Respiratory Therapy Systems 69
Pre-Connected Foley Catheter Urine Drainage Bag Trays 18	Respiratory Masks and Nasal Cannulas 70
Add-a-Foley Catheter Trays 20	Nebulizer Kits 73
Pre-Connected Foley Catheter Urine Meter Kits 21	Accessories 74
Add-a-Foley Urine Meter Trays 23	Aerosol Drainage Bag 74
Foley Insertion Trays 24	Heat and Moisture Exchangers 75
Urinary Drainage Bags 25	Tracheostomy Clean and Care Trays 75
Urine Meters 26	Sterile Water and Saline 76
Foley Catheters 27	Suction Catheters, Trays and Kits 77
Hydrophilic Latex Foley Catheters 27	Suction Tubing with Connectors 80
Silicone-Coated Latex Foley Catheters 28	Suction Yankauers and Kits 81
100% Silicone Foley Catheters 29	Pre-Connected Yankauer and Tubing 81
FlexGRIP® Foley Securement Device 30	SteriCadi™ Suction Device Holster System 82
Intermittent Urinary Catheters 30	AMSure® Surgical & Nursing Systems 83
Urethral Insertion Trays and Kits 30	Specimen Containers 84
Closed System Urethral Drainage Kits 31	Vaginal Speculums 84
Self-Catheterization Kits with R-Polished™ Eye Catheters 32	Disposable Umbilical Cord Clamp 84
PVC Intermittent Urethral Catheter with R-Polished Eyes 33	Ear/Ulcer Syringes 85
Red Rubber Latex Urethral Catheter 34	Cleansing Enema Sets 85
Irrigation Sets 35	Connectors 86
Irrigation Trays and Kits 36	AMSINO Kit Components 87
Ancillary & Patient Care Products 37	BZK Antiseptic 87
Syringes and Specimen Containers 37	Alcohol Prep Pads 87
Urinary Leg Bags 38	Providone-Iodine 88
Extension Tube for Leg Bags 38	Wound and Device Irrigation 88
Pediatric Urine Collector 38	Prefilled Syringes with Enhanced Packaging 88
	Index of REFs 89



ENTERAL FEEDING SYSTEMS

Optimizing Nutritional Delivery



AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and *ENFit*® connections. The systems are designed to provide best patient outcomes and includes:

<i>PUGGLE</i> ® Enteral Feeding Pump	2
<i>Sentinel</i> ™ Enteral Feeding Pump	3
<i>ENFit</i> & Non- <i>ENFit</i> Feeding Sets	3
<i>ENFit</i> & Non- <i>ENFit</i> Enteral Syringes	4
<i>ENFit</i> & Non- <i>ENFit</i> Irrigation Kits & Trays	5

PUGGLE® Enteral Feeding Pump

The versatile pump is designed for pediatric and adult ambulatory applications. Simplified program settings make feeding on-the-go a snap.

PUGGLE is only available the US and greater China, including the People’s Republic of China and the Hong Kong SAR, the Macau SAR and Taiwan.

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours

Dimensions: 5.4 x 5.0 x 2.5 inches (13.7 x 12.7 x 6.4 cm)
Weight: 1.43 lbs (649 g)
Operating Range: 1 to 400 mL/h in 1 mL increments
Delivered Volume: 1 – 9999 mL in 1 mL increments
Accuracy: +/- 5%



Scan the code below to view the **PUGGLE** instructional video



REF	DESCRIPTION	QTY/UNIT
PG001	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
PE90ST	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit® and Transition Connector	30/Case
PE90B5T	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
PE90B12T	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case

Sentinel™ Enteral Feeding Pump



The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The *AMSure*® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each

ENFit® and Non-ENFit Enteral Feeding Sets



ENF1200K

AMSure Enteral Feeding Sets are available with and without *ENFit* connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. *AMSure* enteral feeding pump sets are recommended to use with the *Sentinel* enteral pump and are compatible with the *Compat*™ and *EntraFlo*™ pumps.

- *ENFit* connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- *ENFit* sets include pre-attached stepped transition connector
- Non-sterile, BPA-free options, not manufactured with natural rubber latex

REF	ENFit CONNECTOR	STANDARD BORE (DEHP FREE)	LARGE BORE (DEHP FREE)	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
				Spike Only No Bag	500 mL	1200 mL	Gravity Set	Pump Set	<i>Sentinel</i> ®	<i>Compat</i> ™/ <i>EntraFlo</i> ™	
E-1004				•				•	•	•	30/Case
E-1200						•		•	•	•	
E-1203						•	•				
ENF1004K ©	•			•				•	•	•	
ENF1200K ©	•					•		•	•	•	
ENF1203K ©	•					•	•				
ENF1203 KD ©	•	•				•	•				
ENF1203 KLD ©	•		•			•	•				
ENF500K ©	•				•			•	•	•	

EntraFlo™ is a registered trademark of Medline Industries, Inc.

Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.

ENFit® is a registered trademark of GEDSA.

AMSure® ENFit® and Non-ENFit Enteral Syringes

AMSure Enteral Syringes are available with or without a transition adapter *ENFit* tip. The *ENFit* adapter syringe tip is compliant with the *ENFit* ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.



Bulb Syringe


- Flexible soft bulb promotes one-handed administration

Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage

Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



REF	ENFit TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010 ©			•						•		50/Case
AS011 ©			•			•				•	50/Case
AS011P ©			•				•			•	50/Case
AS115 ©		•			•	•				•	50/Case
AS116 ©		•			•			•			30/Case
ENS015 ©	•			•		•				•	50/Case
ENS015NS	•			•				•			30/Case
ENS115 ©	•				•	•				•	50/Case
ENS115NS ©	•				•			•			30/Case
ENS115T ©*	•				•			•			30/Case

*Includes Transition Adapter

MANUFACTURED IN CANADA											
AS115CA		•			•	•				•	50/Case
AS116CA		•			•			•			30/Case
ENS115CA	•				•	•				•	50/Case
ENS115NSCA	•				•			•			30/Case

ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



Enteral Feeding Irrigation Kits

REF	ENFit SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121 ©		•	•				•	30/Case
AS126 ©		•		•			•	
AS127 ©		•		•		•	•	
ENK127	•	•		•		•	•	

Enteral Feeding Irrigation Trays

REF	ENFit SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	
ENK136	•	•	•		•	•	•	•	

AMSafe®

MEDICATION DELIVERY SYSTEMS

Safety First




The *AMSafe* IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:

- Gravity Administration IV Sets.....7
 - 10 Drops/mL.....7
 - 15 Drops/mL.....9
 - 60 Drops/mL..... 11
 - AMSafe-3®* Selectable Drip Chamber, 10, 15 & 60 drops/mL 12
 - Secondary Piggyback Sets 13
- IV Extension Sets 13
 - Standard Extension Sets..... 13
 - Microbore Extension Sets..... 14
 - Smallbore Extension Sets 14
- Needle-Free Connectors 15
- Vial Access Device..... 15
- IV Start Kit 15
- Sure-Lok®* Needle-Free IV Pump Administration Sets..... 16

AMSafe® Gravity Administration IV Sets

AMSINO's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a roller clamp and a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex



108310

10 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201 ©	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301 ©	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
108303 ©	10	83 in. (211 cm)	<i>Sure-Lok</i> ® Needle-free Y-site	16 mL	50/Case
1010001 ©	10	100 in. (254 cm)	Vented/Non-vented Spike, Roller Clamp, <i>Sure-Lok</i> Needle-free Y-site	26.7 mL	50/Case
108304	10	89 in. (226 cm)	(2) <i>Sure-Lok</i> Needle-free Y-sites with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
109602 ©	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
108306 ©	10	83 in. (211 cm)	Pre-pierced Y-site, <i>Sure-Lok</i> Needle-free Y-site	16 mL	50/Case
108310	10	83 in. (211 cm)	Vented/Non-vented Spike, (2) <i>Sure-Lok</i> Needle-free Y-sites	16 mL	50/Case
AC3110 ©	10	105 in. (266 cm)	(2) <i>Sure-Lok</i> Needle-free Y-sites	20 mL	50/Case

AMSafe® Gravity Administration IV Sets

10 Drops/mL

General Primary Sets (continued)

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
108305	10	89 in. (226 cm)	Pre-pierced Y-site, <i>Sure-Lok</i> ® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case



108308	10	83 in. (211 cm)	(2) <i>Sure-Lok</i> Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
---------------	----	-----------------	---	---------	---------




108309	10	83 in. (211 cm)	<i>Sure-Lok</i> Needle-free Y-site, Injection Y-site	21.8 mL	50/Case
---------------	----	-----------------	--	---------	---------




10 Drops/mL

Flow Controller Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS100	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp	26.8 mL	50/Case



AFS102 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, <i>Sure-Lok</i> Needle-free Y-site	18.5 mL	50/Case
-----------------	----	------------------	---	---------	---------



AFS208 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) <i>SSure-Lok</i> Needle-free Y-sites	19.2 mL	50/Case
-----------------	----	------------------	---	---------	---------



AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) <i>Sure-Lok</i> Needle-free Y-sites, Slide Clamp	19 mL	50/Case
----------------	----	-----------------	--	-------	---------




AMSafe® Gravity Administration IV Sets

10 Drops/mL


General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
1010001B ©	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case




108302 ©	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16.4 mL	50/Case
-----------------	----	-----------------	--	---------	---------




15 Drops/mL

General Primary Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AA3101 ©	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case



158309 ©	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok Needle-free Y-site	17 mL	50/Case
-----------------	----	-----------------	--	-------	---------



158306 ©	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
-----------------	----	-----------------	---	-------	---------



159602 ©	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
-----------------	----	-----------------	-----------------------	---------	---------



1510501 ©	15	105 in. (266 cm)	Vented/Non-vented, (2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
------------------	----	------------------	---	-------	---------



158305 ©	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
-----------------	----	-----------------	---	---------	---------



AMSafe® Gravity Administration IV Sets

15 Drops/mL

Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AS3100 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case




AS3101 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok® Needle-free Y-sites	20.8 mL	50/Case
-----------------	----	------------------	---	---------	---------



AS3102 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites	20.4 mL	50/Case
-----------------	----	------------------	---	---------	---------





AS3110A ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Slide Clamps, (3) Sure-Lok Needle-free Y-sites	22.3 mL	50/Case
------------------	----	------------------	---	---------	---------



AMSafe® Gravity Administration IV Sets

60 Drops/mL

General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608301	60	83 in. (211 cm)	Injection Y-site	22.8 mL	50/Case
					
608304	60	89 in. (226 cm)	(2) Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site, with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
608310	60	83 in. (211 cm)	Vented/Non-vented Spike, (2) Sure-Lok Needle-free Y-sites	16 mL	50/Case
					




AMSafe-3® Selectable Drip Chamber

10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes a roller clamp and a rotating male luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
D38301 ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, <i>Sure-Lok</i> ® Needle-free Y-site	27.7 mL	50/Case
					
D38302 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, <i>Sure-Lok</i> Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, <i>Sure-Lok</i> Needle-free Y-site	30.2 mL	50/Case
					
D38303 ©	10, 15 or 60	89 in. (226 cm)	Backcheck Valve, Pre-pierced Y-site, <i>Sure-Lok</i> Needle-free Y-site, Pre-Attached Ext Set, Pinch Clamp	24 mL	50/Case
					
D38304 ©	10, 15 or 60	103 in. (262 cm)	Backcheck Valve, Pre-pierced Y-site with Pre-attached EXT Set with Pinch Clamp, <i>Sure-lok</i> Needle-free Y-site, Pre-pierced Y-site	34 mL	50/Case
					
D38307	10, 15 or 60	82 in. (208 cm)	Pre-pierced Y-site with Pre-Attached Ext Set with <i>Sure-Lok</i> Needle-Free Connector, Pinch Clamp	22 mL	50/Case
					
D38310	10, 15 or 60	83 in. (210 cm)	(2) <i>Sure-Lok</i> Needle-free Y-site	42.4 mL	50/Case
					

AMSafe® Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AS3515 ©	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
AS3500 ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case



AMSafe® IV Extension Sets

AMSINO's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.



AE3108

Standard Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
000606 ©	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Slide Clamp, Sure-Lok® Needle-free Y-site™	2.5 mL	50/Case
AE3107 ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
AE3108 ©	8 in. (20 cm)	Removable Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
AE3108B	8 in. (20 cm)	Bonded Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
003201 ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Injection Y-sites	4.9 mL	50/Case
003202 ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
IVE2500	100 in. (254 cm)	Female Luer Lock, Slide Clamp	27 mL	50/Case

AMSafe® IV Extension Sets

AMSafe Microbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AE0106 ©	6 in. (15 cm)	Removeable <i>Sure-Lok</i> ® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
AE1106 ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
				
AE0107 ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, <i>Sure-Lok</i> Needle-free Y-site	0.9 mL	50/Case
				
AE1112 ©	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
				
AE1138 ©	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case
				
IVE2501	100 in. (254 cm)	Female Luer Lock, Slide Clamp	9 mL	50/Case
				

AMSafe Flow Controller Extension Sets

AMSINO's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS107 ©	20 in. (51 cm)	Flow Controller, Slide Clamp, Injection Y-site	3.5 mL	50/Case
				
AFS109	19 in. (51 cm)	Flow Controller, Slide Clamp, <i>NeuFlo</i> ® Needle-free Y-site	2.83 mL	50/Case
				

AMSafe Pressure Infusion Extension Sets

AMSINO's single-patient-use pressure infusion IV extension set is made to withstand up to 400 psig for rapid administration of medication or fluids in an emergency. Includes a rotating male luer lock.


REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AE0108-HP	8 in. (20 cm)	Neutral High-Pressure Extension Set w/ Bonded Needle-Free Valve, Slide Clamp, (400psi)	0.30 mL	50/Case
				

AMSafe® Needle-Free Connectors

Our AMSafe needle free valves are manufactured to provide both clinicians and patients safe procedural outcomes.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AY0200 ©	Sure-Lok® Needle-free Connector	0.33 mL	100/Case
AY0206 ©	NeuFlo® Neutral Pressure Needle-free Connector	0.08 mL	100/Case



AMSafe® Vial Access Device

- Microbial barrier swabable silicone stem with a split septum
- Flat seal to allow enhanced disinfection
- Compatible with ISO 11040-4, ISO 80369-7, Luer lock, Luer slip and prefilled glass syringes
- Transparent housing to allow clear view into the fluid pathway
- Compatible with TPN solutions, DMSO, blood and oncology drugs
- ETO Sterilized
- Latex-free, PVC-free, DEHP-free



REF	DESCRIPTION	QTY/UNIT
AY0100	20 mm, Sterile	50/Box

AMSINO IV Start Kits

Amsino IV Start Kits offer efficient, organized solutions for initiating IV therapy.

- Designed for IV initiation
- Sequentially packaged components
- Latex-free tourniquet included
- Convenient all-in-one kit
- Customizable configurations available








REF	DESCRIPTION	QTY/UNIT
IVSK005	10cc PFIV Flush, 8" Extension Set, Gauze, Alcohol Prep Pad, Tape, Tegaderm, Tourniquet, Patient Label, PE pouch, Non- Sterile	50/Case

AMSafe® IV Pump Administration Sets

Sure-Lok® Needle-Free IV Pump Administration Sets

- FDA cleared IV pump sets for use with *Baxter Spectrum* and *Flo-Gard®* 6201 and 6301 series pumps
- *Sure-Lok* Needle-free connector and injection site options
- Includes a roller clamp and a rotating male luer lock
- 0.22 and 1.2 micron filter options available
- *TUBETECH®* technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP092	10	92 in. (234 cm)	IV Pump Set	14 mL	50/Case
					
FP092Y1	10	92 in. (234 cm)	<i>Sure-lok</i> Needle-free Valve, Slide Clamp	14 mL	50/Case
					
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 <i>Sure-lok</i> Needle-free Valves, 2 Slide Clamps	16 mL	50/Case
					
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 <i>Sure-lok</i> Needle-free Valves, 2 Slide Clamps	18 mL	50/Case
					
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 <i>Sure-lok</i> Needle-free Valves, 2 Slide Clamps, 0.22 Micron Filter	22 mL	50/Case
					
FP100FY1 ©	10	100 in. (254 cm)	<i>Sure-lok</i> Needle-free Valve, 2 Slide Clamps, 0.22 Micron Filter	22 mL	50/Case
					
FP100F12Y1	10	100 in. (254 cm)	<i>Sure-lok</i> Needle-Free Valve, Slide clamp, 1.2 Micron Filter	22 mL	50/Case
					



UROLOGICAL CARE

Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

Foley Catheter Trays & Collection	18
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	18
Add-a-Foley Catheter Trays	20
Pre-Connected Foley Catheter Urine Meter Kits	21
Add-a-Foley Urine Meter Trays	23
Foley Insertion Trays.....	24
Urinary Drainage Bags	25
Urine Meters.....	26
Foley Catheters	27
Hydrophilic Latex Foley Catheters	27
Silicone-Coated Latex Foley Catheters	28
100% Silicone Foley Catheters	29
Intermittent Urinary Catheters	30
Urethral Insertion Trays and Kits	30
Closed System Urethral Drainage Kits	31
Self-Catheterization Kits with <i>R-Polished™</i> Eye Catheters.....	32
PVC Intermittent Urethral Catheter with <i>R-Polished</i> Eyes	33
Red Rubber Latex Urethral Catheter	34
Irrigation Sets.....	35
Irrigation Trays and Kits.....	36
Ancillary & Patient Care Products.....	37
Syringes	37
Specimen Containers	37
Urinary Leg Bags	38
Extension Tube for Leg Bags	38
Pediatric Urine Collector	38

AMSure® Foley Catheter Trays & Collection



AS89316

AMSure Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

Trays include: 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S ©	14 Fr	5cc	10/Case
	AS89316S ©	16 Fr		
	AS89318S ©	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89314 ©	14 Fr	5cc	10/Case
	AS89316 ©	16 Fr		
	AS89318 ©	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr	5cc	10/Case
	AS89316H	16 Fr		
	AS89318H	18 Fr		

AMSure® Foley Catheter Trays & Collection



AS89416

AMSure Pre-Connected Foley Catheter Urine Drainage Bag Trays with FlexGRIP®

Now designed for enhanced comfort, safety, and infection prevention with the addition of FlexGRIP Securement Device, enhanced skin prep, and updated IFU.

- **Infection Control Focus:** Features sterile gloves, PVP swab sticks, a BZK wipe, and our AMSure SteriGuide IFU to support proper aseptic catheter insertion
- **FlexGRIP Securement Device:** Improves patient comfort and reduces catheter-related complications
- **Pre-Attached Catheter with Tamper-Evident Seal:** Helps prevent disconnection and maintains a sterile fluid pathway

Trays include: 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, 3 PVP swab sticks, 1 BZK wipe, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag, 1 FlexGRIP securement device, SteriGuide IFU.

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89414S	14 Fr	5cc	10/Case
	AS89416S	16 Fr		
	AS89418S	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89414	14 Fr	5cc	10/Case
	AS89416	16 Fr		
	AS89418	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89414H	14 Fr	5cc	10/Case
	AS89416H	16 Fr		
	AS89418H	18 Fr		

AMSure® Foley Catheter Trays & Collection

AMSure Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex



AS89300

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PREFILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10cc	10/Case
AS89330 ©	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30cc	10/Case

AMSure Add-a-Foley Catheter Trays with FlexGRIP®

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, 3 PVP swab sticks, 1 BZK wipe, lubricating jelly, syringe prefilled with 10cc or 30cc sterile water, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag, 1 FlexGRIP securement device, SteriGuide IFU.



AS89300

REF	DESCRIPTION	PREFILLED SYRINGE	QTY/UNIT
AS89400	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10cc	10/Case
AS89430	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30cc	10/Case

AMSure® Foley Catheter Trays & Collection

AMSure Pre-Connected Foley Catheter Urine Meter Kits

The *AMSure* pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs). (HCPCS codes: See product chart below)

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag



AS89814S

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89814S ©	14 Fr	5cc	10/Case
	AS89816S ©	16 Fr		
	AS89818S ©	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89814 ©	14 Fr	5cc	10/Case
	AS89816 ©	16 Fr		
	AS89818 ©	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89814H ©	14 Fr	5cc	10/Case
	AS89816H ©	16 Fr		
	AS89818H ©	18 Fr		

AMSure® Foley Catheter Trays & Collection



AMSure Pre-Connected Foley Catheter Urine Meter Kits with FlexGRIP®

AS89914S

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 3 PVP swab sticks, 1 BZK wipe, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, 1 FlexGRIP securement device, SteriGuide IFU.

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89914S	14 Fr	5cc	10/Case
	AS89916S	16 Fr		
	AS89918S	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89914	14 Fr	5cc	10/Case
	AS89916	16 Fr		
	AS89918	18 Fr		
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89914H	14 Fr	5cc	10/Case
	AS89916H	16 Fr		
	AS89918H	18 Fr		

Prefilled Sterile Water Syringe for Catheter Balloon Inflation

Prefilled, ready-to-use sterile water syringe designed to ensure accurate, consistent balloon inflation while minimizing contamination risk. Packaged for convenience and sterility, it supports efficient clinical workflows and reliable catheter performance.



REF	DESCRIPTION	QTY/UNIT
KL20100	10 ml Prefilled Sterile Water Syringe for Foley Catheter Balloon Inflation, Sterile	200/PE bag, 400/Case
KL20101	30 ml Prefilled Sterile Water Syringe for Foley Catheter Balloon Inflation, Sterile	20/PE bag, 240/Case

AMSure® Foley Catheter Trays & Collection



AS89800

AMSure Add-a-Foley Urine Meter Trays

The *AMSure* Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. (HCPCS Code A4354)

- Sterile, not manufactured with natural rubber latex

Trays include: 400 ml Urine Meter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag

REF	DESCRIPTION	PREFILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30cc	10/Case

AMSure Add-a-Foley Urine Meter Trays with *FlexGRIP*®



AS89900

Trays include: 400 ml Urine Meter, waterproof and fenestrated drapes, (2) Nitrile gloves, 3 PVP swab sticks, 1 BZK wipe, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, 1 *FlexGRIP* securement device, SteriGuide IFU.

REF	DESCRIPTION	PREFILLED SYRINGE	QTY/UNIT
AS89900	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag, <i>FlexGRIP</i>	10cc	10/Case
AS89930	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag, <i>FlexGRIP</i>	30cc	10/Case

AMSure® Foley Catheter Trays & Collection



AMSure Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or (6) Povidone-iodine swab sticks, 1 prefilled saline syringe, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

REF	PREFILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS880 ©	10cc	•		20/Case
AS890 ©	30cc	•		
AS886	10cc		•	
AS896	30cc		•	

AMSure Foley Insertion Trays with FlexGRIP®

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), 1 prefilled saline syringe, lubrication jelly, fenestrated drape, 1 FlexGRIP Securement Device. SteriGuide IFU.

REF	PREFILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS882	10cc	•		20/Case
AS892	30cc	•		

AMSure® Foley Catheter Trays & Collection

AMSure Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube. (HCPCS code A4358)

- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex



AS322 Pear



AS326 Oval



AS32600 Low-Profile

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312 ©	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322 ©	Pear	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	
AS332 ©	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	
AS302 ©	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	
AS32200 ©	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600* ©	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	
AS326* ©	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	

*100% sterile packed in peel pouch

AMSure® Foley Catheter Trays & Collection

AMSure Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. (HCPCS code A4358)

- **Anti-Reflux Chamber helps improve patient safety by reducing urine backflow**
- **Twist valve drains at the bottom of the urine meter for easy drainage**
- **Leak-proof seams ensure fluid containment**
- **Needle-free sampling port reduces the risk of needle stick injuries**
- **Air vents help to maintain consistent urine flow delivering safer patient outcomes**
- **Sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS327200 ©	200cc Urine Meter with 2000 mL Urine Drainage Bag	10/Case
AS328 ©	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case



ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)

- **Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag**
- **Needle-free sampling port with luer lock access**
- **Leak-proof seams ensure fluid containment**
- **Meter bottle capacity 150 mL**
- **Adjustable bed sheet clamp for securing drainage bag tubing**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS327100 ©	ACCU-MEASURE Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case




AMSure® Foley Catheters

AMSure Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
5cc	AS41014H©	14 Fr	16 in. (40.6 cm)	10/Box
	AS41016H©	16 Fr		
	AS41018H©	18 Fr		
30cc	AS42014H	14 Fr	16 in. (40.6 cm)	10/Box
	AS42016H©	16 Fr		
	AS42018H©	18 Fr		

Coude Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5cc Capacity: 5cc	AS41012HC	12 Fr. × 16 in. (40.6 cm)	10/Box
	AS41014HC	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016HC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018HC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS41020HC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS41022HC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS41024HC	24 Fr. × 16 in. (40.6 cm)	10/Box
30cc Capacity: 30cc	AS42016HC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS42018HC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS42020HC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS42022HC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS42024HC	24 Fr. × 16 in. (40.6 cm)	10/Box

AMSure® Foley Catheters

AMSure Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
3cc	AS40008 ©	8Fr	10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr		
5cc	AS41012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS41014 ©	14 Fr		
	AS41016 ©	16 Fr		
	AS41018 ©	18 Fr		
	AS41020 ©	20 Fr		
	AS41022 ©	22 Fr		
	AS41024 ©	24 Fr		
	AS41026 ©	26 Fr		
	AS41028 ©	28 Fr		
	AS41030 ©	30 Fr		
30cc	AS42012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS42014 ©	14 Fr		
	AS42016 ©	16 Fr		
	AS42018 ©	18 Fr		
	AS42020 ©	20 Fr		
	AS42022 ©	22 Fr		
	AS42024 ©	24 Fr		
	AS42026 ©	26 Fr		
	AS42028 ©	28 Fr		
	AS42030 ©	30 Fr		

3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	AS43016 ©	16 Fr	16 in. (40.6 cm)	10/Box
	AS43018 ©	18 Fr		
	AS43020 ©	20 Fr		
	AS43022 ©	22 Fr		
	AS43024 ©	24 Fr		
	AS43026 ©	26 Fr		

AMSure® Foley Catheters

AMSure 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
3cc	AS40008S ©	8Fr	10.6 in. (27cm)	10/Box
	AS40010S ©	10 Fr		
5cc	AS41012S ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS41014S ©	14 Fr		
	AS41016S ©	16 Fr		
	AS41018S ©	18 Fr		
	AS41020S ©	20 Fr		
	AS41022S ©	22 Fr		
	AS41024S ©	24 Fr		
	AS41026S ©	26 Fr		
30cc	AS42016S ©	16 Fr	16 in. (40.6 cm)	10/Box
	AS42018S ©	18 Fr		
	AS42020S ©	20 Fr		
	AS42022S ©	22 Fr		
	AS42024S ©	24 Fr		


3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	AS43018S ©	18 Fr	16 in. (40.6 cm)	10/Box
	AS43020S ©	20 Fr		
	AS43022S ©	22 Fr		
	AS43024S ©	24 Fr		
	AS43026S ©	26 Fr		

FlexGRIP® Foley Securement Device

Next-Generation Securement Technology

FlexGRIP combines innovative stabilization technology with a patient-centered design to create a securement solution clinicians can trust. The large adhesive footprint and soft cushioning materials provide exceptional hold while remaining gentle on the skin.



REF	DESCRIPTION	QTY/UNIT
FOL1118	Foley Catheter Fixation System	25/Box

AMSure® Intermittent Urinary Catheters

AMSure Urethral Insertion Trays and Kits

AMSure Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection**

AMSure PVC Catheter Tray

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS870

AMSure Red Rubber Catheter Trays

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
AS870 ©	14 Fr PVC*	20/Case
AS87114 ©	14 Fr Red Rubber Latex	

*Not manufactured with natural rubber latex.

AMSure® Intermittent Urinary Catheters



AS89214R

AMSure Closed System Urethral Drainage Kits

AMSure Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection**
- **Sterile**

AMSure Red Rubber Catheter Kits

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

AMSure PVC Catheter Kits

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	LENGTH	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	AS89214R ©	14 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	20/Case
	AS89214	14 Fr	16 in. (40.6 cm)	PVC Latex- Free	
BZK Swab Sticks	AS89214RK ©	14 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	20/Case
	AS89216RK ©	16 Fr	16 in. (40.6 cm)	Red Rubber (Contains Latex)	

AMSure® Intermittent Urinary Catheters



AS85010

AMSure Self-Catheterization Kits with R-POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with *R-Polished* eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. *R-Polished* technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Introducer tip reduces the potential for urinary tract infection (UTI)**
- **Collection bag with built-in thumb holes allows for ease of use**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

SWAB STICK	REF	SIZE	LENGTH	QTY/UNIT
PVP Swab Sticks	AS85008 ©	8 Fr	16 in. (40.6 cm)	50/Case
	AS85010 ©	10 Fr		
	AS85012 ©	12 Fr		
	AS85014 ©	14 Fr		
	AS85016 ©	16 Fr		
BZK Swab Sticks	AS85008K ©	8 Fr	16 in. (40.6 cm)	50/Case
	AS85010K ©	10 Fr		
	AS85012K ©	12 Fr		
	AS85014K ©	14 Fr		
	AS85016K ©	16 Fr		

AMSure® Intermittent Urinary Catheters

AMSure PVC Intermittent Urethral Catheter with R-POLISHED™ Eyelets


AMSure PVC Intermittent Catheters with *R-Polished* technology use radio frequency to create ultra-smooth eyelets for optimal patient comfort. They are designed to provide a smooth surface in a catheter that is firm enough for insertion, yet soft enough to be slightly pliable. These urethral catheters are designed for intermittent or short-term use and are available in various French sizes, as well as female, pediatric and male lengths. (HCPCS code A4351)

- Advanced *R-Polished* eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- PVC material eliminates latex allergy concerns
- Clear shaft provides easy viewing as you measure
- Short-term catheter use helps decrease infections
- Individually packaged, sterile
- BPA-free, not manufactured with natural rubber latex



FEATURES	REF	SIZE	LENGTH	QTY/UNIT
Female	AS860608 ©	8 Fr	6 in. (15.24 cm)	50/Case
	AS860610 ©	10 Fr		
	AS860612 ©	12 Fr		
	AS860614 ©	14 Fr		
	AS860616 ©	16 Fr		
Pediatric	AS861006 ©	6 Fr	10 in. (25.40 cm)	50/Case
	AS861008 ©	8 Fr		
	AS861010 ©	10 Fr		
Male	AS861608 ©	8 Fr	16 in. (40.64 cm)	50/Case
	AS861610 ©	10 Fr		
	AS861612 ©	12 Fr		
	AS861614 ©	14 Fr		
	AS861616 ©	16 Fr		
	AS861618 ©	18 Fr		
Coudé Tip	AS861610C	10 Fr	16 in.(40.65 cm)	30/Case
	AS861612C	12 Fr		
	AS861614C	14 Fr		
	AS861616C	16 Fr		
	AS861618C	18 Fr		

AMSure® Intermittent Urinary Catheters



AMSure Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free


REF	SIZE	LENGTH	QTY/UNIT
AS44008 ©	8 Fr	11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr		
AS44012 ©	12 Fr	16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr		
AS44016 ©	16 Fr		
AS44018 ©	18 Fr		
AS44020 ©	20 Fr		
AS44022 ©	22 Fr		
AS44024 ©	24 Fr		
AS44028 ©	28 Fr		

AMSure® Irrigation Sets

AMSure Cystoscopy & TUR Irrigation Sets

The *AMSure* Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
ABL400 ©	Single Lead	20/Case



REF	TYPE	QTY/UNIT
ABL402 ©	Two Leads	20/Case



REF	TYPE	QTY/UNIT
ABL404 ©	Four Leads	20/Case

AMSure® Irrigation Trays and Kits

AMSure Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need. (HCPCS code A4320)

AMSure Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

AMSure Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector



AS121



AS126

AMSure Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121 ©		•		30/Case
AS126 ©	•			
AS127 ©	•		•	

AMSure Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		



AS136

AMSure® Ancillary & Patient Care Products

AMSure Syringes (HCPCS code A4322)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap

AS010

AS016NS

AS115

TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS010 ©	Bulk	Non-Sterile			50/Case
	AS011 ©	Poly Pouch	Sterile (EO)	•		50/Case
	AS011P ©	Form Fill Seal	Sterile (EO)	•		50/Case
Piston Thumb Control Ring	AS015 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS015CA*	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016NS	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
	AS016NSCA*	Resealable Hanger Bag	Non-Sterile	•	•	50/Case
Piston Flat Top	AS115 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116 ©	Resealable Pole Bag	Non-Sterile	•	•	30/Case

*Manufactured in Canada

AMSure Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples. (HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
AS340 ©		•		•		•	•	100/Case
AS341 ©			•	•	•	•	•	100/Case
AS343 ©			•				•	500/Case
AS345 ©	•			•			•	100/Case

AS340

AMSure® Ancillary & Patient Care Products

AMSure Urinary Leg Bags (HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex

SIZE	REF	DRAIN PORT STYLE	QTY/UNIT
600 mL	AS306N ©	Push-Pull	48/Case
	AS308N ©	Twist-Turn	
900 mL	AS307N ©	Push-Pull	48/Case
	AS309N ©	Twist-Turn	



AS307N

AMSure Extension Tube for Leg Bags (HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS311 ©	18 in. (45.72 cm) Extension Tube	50/Case



AMSure Pediatric Urine Collector (HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS409 ©	200 mL Pediatric Urine Bag	50/Box



AMSure®

STERILE WATER AND SALINE PRODUCTS

Inhalation, Irrigation, and Injection Solutions



AMSure Respiratory Therapy Systems are designed to deliver optimal FiO₂ levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Inhalation Solutions.....	40
Irrigation Solutions	41
Nebulizer Kits	73
Device and Wound Irrigation	42
Saline Products	43
MedXL Saline Flush Syringes.....	45
Prefilled Syringes	47
PraxiMed Syringe Filling and Labeling System	49

**Manufactured
in USA**

AMSure® Inhalation Solutions

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.

AMSure Sterile Water for Inhalation (Bottles)

- 500 mL and 1000 mL bottles Sterile Water, USP
- For use with nebulizers and humidifiers
- Hanging Spikeable bottle helps solve the problem of look-alike solution “bags” and prevents mishaps from occurring – Use with vented spikes only



AS1065

AS1000

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	—	500 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	
AS1000	AS1000E ©	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	

AMSure Sterile Water Unit Dose Vials

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case (10 Boxes – 100 per Box)
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	

AMSure Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

AMSure® Inhalation Solutions

AMSure 0.9% Sodium Chloride Unit Dose Vials

- 3 mL, 5 mL and 15 mL Sterile Sodium Chloride Inhalation Solution (0.9%), USP
- Single-use vials with graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case (10 Boxes – 100 per Box)
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case (10 Boxes – 100 per Box)
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case (3 Boxes – 48 per Box)



ASR0039



ASR0059



ASR0159

AMSafe® Irrigation Solutions



AS399



AS398

AMSafe Sterile Water and Saline

- 100 mL sterile water or saline bottle with screw cap and inner seal
- Easy open twist screw top container
- Used for device and wound irrigation

REF	DESCRIPTION	QTY/UNIT
AS399 ©	100 mL Sterile Water for Device and Wound Irrigation	48/Case
AS398 ©	100 mL Sterile Saline 0.9% for Device and Wound Irrigation	48/Case

AMSure® Device and Wound Irrigation



AMSure Sterile Water and 0.9% Sodium Chloride for Irrigation, USP

- Used for moistening of wound dressings, wound debridement, and device irrigation
- Bottle shape allows for easy handling and pouring
- Variety of sizes and packaging configurations offer convenience and flexibility
- Manufactured and sterilized in the USA

REF	DESCRIPTION	QTY/UNIT
AS-6240	100 mL 0.9% Sodium Chloride for Irrigation, USP	48/Case
AS-6250	100 mL Sterile Water for Irrigation, USP	48/Case
AS-6270	250 mL 0.9% Sodium Chloride for Irrigation, USP	24/Case
AS-6260	250 mL Sterile Water for Irrigation, USP	24/Case
AS-6380	500 mL 0.9% Sodium Chloride for Irrigation, USP	18/Case
AS-6390	500 mL Sterile Water for Irrigation, USP	18/Case
AS-6481	1000 mL 0.9% Sodium Chloride for Irrigation, USP	9/Case
AS-6491	1000 mL Sterile Water for Irrigation, USP	9/Case

AMSure Sal-Sense™ Single Use Saline Topical Solution 0.9% Sodium Chloride

- Single patient use avoids cross contamination
- Save time from standard practice of drawing solution
- Twist design to open effortlessly
- Easy to irrigate area in any position



REF	DESCRIPTION	QTY/UNIT
AS0159	15 mL Single Use Saline Topical Solution	144/Case (3 Boxes – 48 per Box)
AS0359	35 mL Single Use Saline Topical Solution	120/Case (6 Boxes – 20 per Box)
AS0559	50 mL Single Use Saline Topical Solution	120/Case (6 Boxes – 20 per Box)

AMSafe® Saline Products

AMSafe Prefilled Saline Flush Syringes

- 3 mL 5 mL, 10 mL and 20 mL syringes, 0.9% Sodium Chloride Injection, USP
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared
- Sterile, non-pyrogenic, DEHP-free, PVC-free, not manufactured with natural rubber latex



12 mL and 20 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1203 ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205 ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210 ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020 ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)
IVF0503M	3 mL USP Sodium Chloride Injection Fill in 5 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF0505M	5 mL USP Sodium Chloride Injection Fill in 5 mL Syringe	180/Case (6 Boxes – 30 per Box)

Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
IVF1203SS ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205SS ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210SS ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020SS ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

AMSafe® Saline Products

AMSafe Prefilled Saline Flush Syringes



10 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1003 ©	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1005 ©	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1010 ©	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1010TM ©	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

Prefilled Sterile Water Syringe for Catheter Balloon Inflation

Prefilled, ready-to-use sterile water syringe designed to ensure accurate, consistent balloon inflation while minimizing contamination risk. Packaged for convenience and sterility, it supports efficient clinical workflows and reliable catheter performance.



REF	DESCRIPTION	QTY/UNIT
KL20100	10 ml Prefilled Sterile Water Syringe for Foley Catheter Balloon Inflation, Sterile	200/PE bag, 400/Case
KL20101	30 ml Prefilled Sterile Water Syringe for Foley Catheter Balloon Inflation, Sterile	20/PE bag, 240/Case

Prefilled Saline Flush Syringes

Praxiject™ Prefilled Saline Flush Syringes

Designed to use for flushing only, with sterile fluid path or sterile field ready configuration, *Praxiject* Low Temperature Sterilization Syringes have the advantage of preventing incidents of transitory taste/odor and reducing risks of nausea during flushing. *Praxiject* Syringes are safe and cost-effective for flushing vascular access.



FEATURES AND BENEFITS	Praxiject SP	Praxiject SF
Easy to open	•	•
Ready to use	•	•
Clear labeling	•	•
Meets USP specifications	•	•
Latex and preservative-free	•	•
Zero reflux design	•	•
Pharmaceutical grade plunger	•	•
Pressure of 14 psi	•	•
Easy glide piston	•	•
Low temperature Gamma sterilization	•	•
Minimal Leachable & Extractables (VOC)	•	•
Minimal Transitory taste/odor	•	•
Precise graduation		•
Peel pouch packaging		•
Sterile field compatible		•
Price	++	+++

Sterile Fluid Pathway Flush Syringes

REF	DESCRIPTION	QTY/UNIT
3705-3US1	10 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe	180/Case (6 Boxes – 30 per Box)
IVF1010PX	10 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe	180/Case (6 Boxes – 30 per Box)

Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
37043US	3 ml of 0.9% Sodium Chloride Injection USP in a 5 ml syringe	720/Case (6 Boxes – 120 per Box)
37053US	3 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe	600/Case (6 Boxes – 100 per Box)
3704US	5 ml of 0.9% Sodium Chloride Injection USP in a 5 ml syringe	720/Case (6 Boxes – 120 per Box)
37055US	5 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe	600/Case (6 Boxes – 100 per Box)
3705CUS	10 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe	600/Case (6 Boxes – 100 per Box)
3706US	20 ml of 0.9% Sodium Chloride Injection USP in a 20 ml syringe	360/Case (6 Boxes – 100 per Box)

Prefilled Saline Syringes with Enhanced Packaging

Praxiject™ EP is an hermetically sealed packaging designed for use in kits and set assembly. This packaging offers stability, and excellent gas barrier properties. The material is non breathable and therefore resistant to ethylene oxide (ETO) gas transmission.



REF	DESCRIPTION	QTY/UNIT
3705-1US	10 ml of 0.9% Sodium Chloride Injection USP in a 10 ml syringe with Enhanced Packaging	460/Case



Sodium Citrate Prefilled Syringes

CitraFlow™ 46.7% SF for Vascular Access DeviceLocking

CitraFlow 46.7% SF prefilled syringes are filled with an anticoagulant and antimicrobial 46.7% sodium citrate solution. CitraFlow 46.7% SF is a proven alternative to heparin to lock Vascular Access Devices. It avoids the risks associated with systemic heparinization and is safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia.



REF (Outside US only)	DESCRIPTION	QTY/UNIT
38143 ☺	CitraFlow 46.7% SF, 3ml 46.7% Sodium Citrate solution in 5ml syringe (SF)	900/Case (6 Boxes – 150 per Box)
38143-1 ☺	CitraFlow 46.7% SF, twinpack of two 3ml 46.7% Sodium Citrate solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

CitraFlow 30% SF for Vascular Access DeviceLocking

CitraFlow 30% SF products are anticoagulant and antimicrobial 30% sodium citrate prefilled syringes. They are a proven alternative to heparin for Vascular Access Devices locking as they avoid the risks associated with systemic heparinization and are safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
38243 ☺	CitraFlow 30% SF, 3ml 30% Sodium Citrate solution in 5ml syringe (SF)	900/Case (6 Boxes – 150 per Box)
38243-1 ☺	CitraFlow 30% SF, twinpack of two 3ml 30% Sodium Citrate solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

CitraFlow 4% SF for Vascular Access Device Locking

CitraFlow 4% SF prefilled syringes are a proven alternative to heparin for locking indwelling catheters on patients. Sodium citrate avoids the risks associated with systemic heparinization and is safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia

REF (Outside US only)	DESCRIPTION	QTY/UNIT
38543 ☺	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in a 5ml syringe (SF)	900/Case (6 Boxes – 150/Box)
38543-1 ☺	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in 5ml syringe, twinpack (SF)	600/Case (6 Boxes – 100/Box)
3854E1 ☺	CitraFlow 4% SF, 5ml 4% Sodium Citrate solution USP in 5ml syringe (SF)	720/Case (6 Boxes – 120/Box)
38553 ☺	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in 10ml syringe (SF)	600/Case (6 Boxes – 100/Box)
38555 ☺	CitraFlow 4% SF, 5ml 4% Sodium Citrate solution USP in 10ml syringe (SF)	600/Case (6 Boxes – 100/ Box)

Sodium Citrate Prefilled Syringes



CitraFlow PLUS™ Prefilled Syringes

CitraFlow PLUS prefilled syringes combines in one syringe the anticoagulant properties of sodium citrate PLUS the antimicrobial properties of ethanol. This solution combination prevents biofilm formation, reduces incidences of CRBSI (catheter related blood stream infections) as well as being a safe alternative to heparin.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
3844	CitraFlow PLUS, 5ml 4% Sodium Citrate/30% Ethanol solution in a 5ml syringe (SF)	720/Case (6 Boxes – 120 per Box)
38443	CitraFlow PLUS, 3mL 4% Sodium Citrate/30% Ethanol solution in a 5mL syringe (SF)	900/Case (6 Boxes – 150 per Box)
38443-1	CitraFlow PLUS, twinpack of two 3ml 4% Sodium Citrate/30% Ethanol solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)
38453	CitraFlow PLUS, 3ml 4% Sodium Citrate/30% Ethanol solution in a 10ml syringe (SF)	600/Case (6 Boxes – 100 per Box)
38455	CitraFlow PLUS, 5ml 4% Sodium Citrate/30% Ethanol solution in a 10ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

PraxiFlow™ Prefilled Syringe Kits for Vascular Access Device Locking

Each PraxiFlow kit consists of one Praxiject™ 0.9% NaCl syringe for Flush and one CitraFlow™ Sodium Citrate syringe for Lock of Vascular Access Devices. All the proven benefits of Praxiject and CitraFlow prefilled syringes in one convenient easy to use package.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
3215 ©	1 Praxiject 0.9% NaCl syringe (10ml in 10cc syringe), 1 CitraFlow 46.7% Sodium Citrate syringe (3ml in 5cc syringe) (SF)	450/Case (6 Boxes – 75/Box)
3225 ©	1 Praxiject 0.9% NaCl syringe (10ml in 10cc syringe), 1 CitraFlow 30% Sodium Citrate syringe (3ml in 5cc syringe)	450/Case (6 Boxes – 75/Box)
3255 ©	1 Praxiject 0.9% NaCl syringe (10ml in a 10cc syringe), 1 CitraFlow 4% Sodium Citrate syringe (3ml in 5cc syringe)	450/Case (6 Boxes – 75/Box)

CitraFlow products are not currently cleared by the US FDA.



Fast, Accurate and Repeatable Syringe Filling and Labeling Made Easy

PraxiMed is an automated system for the preparation of medications in syringes. Syringe based systems have a proven track record. Ready to use prefilled syringes represent the safest and most cost-effective method for administering medications. When commercial prefilled syringes are not available the next best option is unit dose syringes prepared by the pharmacy.



PraxiFill™ Automated Syringe Filler

PraxiFill is an automated syringe filler intended to increase productivity of the syringe filling process by eliminating the tedious manual task of aspirating solution into the syringe and capping the filled syringe.

PraxiFill is designed to be used with MedXL pre-capped sterile syringes for increased efficiency and productivity.

3cc, 5cc, 10cc, 20cc and 50cc syringes are available.

Filling speed is up to 8 syringes/minute depending on syringe size.

- **Fast, accurate, repeatable filling and capping**
- **No need for calibration, auto zero on all functions**
- **Compact, small footprint**
- **Can fill 3cc, 5cc, 10cc, 20cc and 50cc syringes**
- **Simple and fast change over between different sizes**
- **Minimum fill level from 10% of syringe nominal capacity**
- **Filling accuracy and repeatability to +/-1%**
- **No drip filling feature**
- **Caps are designed for freezing and thawing**



PraxiLabel™ Semi-Automatic Syringe Labeler

PraxiLabel is a semi-automatic tabletop syringe labeler intended to increase labeling productivity as it eliminates the work related injuries due to frequent hand and wrist movements involved in the syringe labeling operation.

The labels are pregraduated and may be printed on any heat transfer label printer.

PraxiLabel is designed for use with the 3cc, 5cc, 10cc, 20cc and 50cc syringe sizes and labels from medXL.

Only one set of 2 way snap-in syringe support rollers covers all syringe sizes (no tools needed for the change).

The PraxiLabel should be used in conjunction with the PraxiFill for labeling the prefilled syringes.

- **Can label 3cc, 5cc, 10cc, 20cc and 50cc syringes**
- **Simple and fast change over between different sizes**
- **Adjustable label position for accurate labeling**

PraxiMed™ System

PraxiFill™ Order Information

REF	DESCRIPTION	QTY/UNIT
PF01	PraxiFill Automated Syringe Filler	Each
PF103	Triple Transfer Spike with Stopcock	40/Box
PF104	Transfer Spike	100/Box
PF201	Syringe Filler Set with Spike	40/Box
PF201-1	Syringe Filler Set	40/Box
PF3022	PraxiFill Vial Adapter, Large Vial Adapter	1/Pk, 45Pks/Box, 4 Boxes/Case(180)
PF3023	PraxiFill Vial Adapter, Small Vial Adapter	1/Pk, 45Pks/Box, 4 Boxes/Case(180)
PF3025	PraxiFill Vial Adapter, Large Vial Adapter with Swab Valve	1/Pk, 45Pks/Box, 4 Boxes/Case(180)
PR10CC02	PraxiFill 10cc Syringe with Cap (Pharmacy Pack)	50 Syringes/Bag, 4 Bags/Case (200)
PR20CC01	PraxiFill 20cc Syringe with Cap (Pharmacy Pack)	25 Syringes/Bag, 5 Bags/Case (125)
PR30CC01	PraxiFill 30cc Syringe with Cap (Pharmacy Pack)	20 Syringes/Bag, 4 Bags/Case (80)
PR3CC01	PraxiFill 3cc Syringe with Cap (Pharmacy Pack)	40 Syringes/Bag, 10 Bags/Case (400)
PR50CC01	PraxiFill 50cc Syringe with Cap (Pharmacy Pack)	15 Syringes/Bag, 3 Bags/Case (45)
PR50CC02	PraxiFill 50cc Syringe with Cap (Pharmacy Pack)	15 Syringes/Bag, 3 Bags/Case (45)
PR5CC01	5cc Syringe with Cap, Sterile (Pharmacy Pack)	25 Syringes/Bag, 12 Bags/Case (300)
PRCAP01	PraxiFill Replacement Cap, Female Luer Lock Cap with Male Cap	6/Pk, 50Pks/Box, 4 Boxes/Case (200)

PraxiLabel™ Order Information

REF	DESCRIPTION	QTY/UNIT
PL01	PraxiLabel Semi-Automatic Syringe Labeler	Each
PL10CC01	PraxiLabel 10cc Syringe White Band Label	2500/Roll, 9 Rolls/Case (22500)
PL20CC01	PraxiLabel 20cc Syringe White Band Label	2500/Roll, 9 Rolls/Case (22500)
PL30CC01	PraxiLabel 30cc Syringe White Band Label	2500/Roll, 9 Rolls/Case (22500)
PL50CC01	PraxiLabel 50cc Syringe White Band Label	2500/Roll, 9 Rolls/Case (22500)

RECEPTAL®

FLUID MANAGEMENT SYSTEMS

Complete Solutions From Collection to Disposal

RECEPTAL® Fluid Management Solutions


AMSINO is a leading global provider of suction fluid management systems. The breadth of AMSINO's RECEPTAL® suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

EZE-VAC® Suction Canister System	52
EZE-VAC 2.0 Suction Canister System.....	53
RECEPTAL Liners	54
RECEPTAL Twist Liners	55
RECEPTAL Outer Shells	56
ACCU-MEASURE® Critical Measurement Devices	57
RECEPTAL Solidifiers	57
AMSure® Suction Tubing with Connectors.....	58
Suction Yankauers and Kits	59
Pre-Connected Yankauer and Tubing.....	59
SteriCadi™ Suction Device Holster System.....	60
Accessories.....	61
Reusable Hardware	64
iRECEPTAL® Surgical Waste Fluid Management System	67

EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD®, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01 ©	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01 ©	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01 ©	3300 mL Canister, Lid and 2 Elbows	12/Case

EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
43210-01 ©	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID x 6 ft. (5 mm x 1.8 m) + 3/16 in ID x 18 in. (5 mm x 46 cm)	30/Case

RECEPTAL® Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
1000	BULK PACKAGED		
	43010-05 ©	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	
1500	BULK PACKAGED		
	43025-01 ©	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	
	43015S ©	Liner with VAC-GARD and 45g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	
	0L654-01 ©	Liner with VAC-GARD and Pour Spout	
2000	BULK PACKAGED		
	43027-01 ©	Liner with Pour Spout,	50/Case
	43041-01 ©	Liner with Ortho Port	
	43042-01 ©	Liner with VAC-GARD	
	43044-05 ©	Liner with VAC-GARD and Pour Spout	
	43020S ©	Liner with VAC-GARD and 60g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	
	0G905-01 ©	Liner with Ortho Port	
	0P507-01 ©	Liner with Pour Spout	
3000	BULK PACKAGED		
	0L656-01* ©	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	0L655-01* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	
	43030S* ©	Liner with VAC-GARDs, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	

* Dual port allows for dual suctioning during procedures.
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID

RECEPTAL® *Twist* Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.


- **Premeasured, preloaded and contact-free Solidifier Dispenser Cap**
- **Integrated solidifier eliminates the need for open ports**
- **All-in-one solidifier and liner solution reduces costs associated with extra inventory**
- **100% leak tested seamless liners minimize risk of contamination**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**



43121S



43120S shown above
in canister 43445-01
(canister not included)



REF	DESCRIPTION	QTY/UNIT
43110S	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43115S ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
43116S ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	
43120S ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
43121S ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	

RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex



Color-coded by size for easy selection.

REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	
43445-01 ©	2000 mL Outer Shell	
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

ACCU-MEASURE® Critical Measurement Devices

ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

RECEPTAL
ACCU-MEASURE
Device



REF	DESCRIPTION	FEATURES	QTY/UNIT
43306-01 ©	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10 ©	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case

RECEPTAL® Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- No shaking required
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4320SP

Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
4305SP	500 Solidifier Pouch, solidifies 500 mL	96/Case
4312SP	1200 Solidifier Pouch, solidifies 1200 mL	
4315SP	1500 Solidifier Pouch, solidifies 1500 mL	
4320SP	2000 Solidifier Pouch, solidifies 2000 mL	
4330SP	3000 Solidifier Pouch, solidifies 3000 mL	

Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
4305SB	500 Solidifier Bottle, solidifies 500 mL	96/Case
4312SB	1200 Solidifier Bottle, solidifies 1200 mL	
4315SB	1500 Solidifier Bottle, solidifies 1500 mL	
4320SB	2000 Solidifier Bottle, solidifies 2000 mL	
4330SB	3000 Solidifier Bottle, solidifies 3000 mL	
43160SB	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case



4315SB

4320SB

AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures.

- Flexible, low-coil memory allows for easy draping.
- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 ☺ (50/Case)	AS820B (100/Case)
	6 ft. (1.8 m)	AS821 ☺ (50/Case)	AS821B (50/Case)
	10 ft. (3.0 m)	AS822 ☺ (50/Case)	AS822B (50/Case)
	12 ft. (3.7 m)	AS823 ☺ (20/Case)	AS823B (50/Case)
	20 ft. (6.1 m)	AS824 ☺ (20/Case)	AS824B (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 ☺ (50/Case)	AS825B (50/Case)
	10 ft. (3.0 m)	AS826 ☺ (50/Case)	AS826B (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C ☺ (50/Case)	—
	12 ft. (3.7 m)	AS827 ☺ (20/Case)	AS827B (50/Case)
	20 ft. (6.1 m)	AS828 ☺ (20/Case)	AS828B (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M ☺ (50/Case)	AS825BM (50/Case)
	10 ft. (3.0 m)	AS826M ☺ (50/Case)	AS826BM (50/Case)
	12 ft. (3.7 m)	AS827M ☺ (20/Case)	AS827BM (50/Case)
	20 ft. (6.1 m)	AS828M ☺ (20/Case)	AS828BM (50/Case)

AMSure® Suction Yankauers and Kits



AMSure Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			AS830 ④ 50/Case	AS830B 500/Case
	•		•	AS831 ④ 50/Case	AS831B 500/Case
		•		AS832 ④ 50/Case	AS832B 500/Case
		•	•	AS833 ④ 50/Case	AS833B 500/Case
Flexible	•			AS834 ④ 50/Case	—
	•		•	AS835 ④ 50/Case	—
		•		AS836 ④ 50/Case	—

AMSure® Pre-Connected Yankauer and Tubing

AMSure Pre-connected yankauer and tubing sets include a rigid non-vented bulbous tip yankauer pre-attached to suction tubing in a variety of sizes. Sterile.


REF	TUBING DIAMETER	TUBING LENGTH	QTY/UNIT
AS823P ④	3/16 in. (5 mm)	12 ft. (3.7m)	20/Case
AS824P ④	3/16 in. (5 mm)	20 ft. (6.1 m)	20/Case
AS825P ④	1/4 in. (6 mm)	6 ft. (1.8 m)	50/Case
AS826P ④	1/4 in. (6 mm)	10 ft. (3.0 m)	50/Case
AS827P ④	1/4 in. (6 mm)	12 ft. (3.7m)	20/Case
AS828P ④	1/4 in. (6 mm)	20 ft. (6.1 m)	20/Case

SteriCadi™ Suction Device Holster System

SteriCadi is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to support improved infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder, Sterile	30/Case



Accessories



RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
43303-01 ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
43303-02 ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	
43303-05 ©	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	



External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
43309-01 ©	External Shut-Off Valve, ¼ in. (0.64 cm) Male × Female Elbow	10/Case



ACCU-MEASURE® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
43305-10 ©	Canister Attachment Bracket for 250 mL Canister	10/Case

Accessories

Suction Tubing

For permanent Outer Shell-to-vacuum source, or Outer Shell-to-Outer Shell connections.

REF	DESCRIPTION	QTY/UNIT
43455-01 ©	100 ft. (30.48 m) PVC Non-Sterile SuctionTubing 9/32 in. (7 mm) I.D.	1/Case
43457-25 ©	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case



EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.

- Non-sterile, BPA-free, not manufactured with natural rubber latex


REF	DESCRIPTION	QTY/UNIT
43317-25 ©	18 in. (0.45 m) Tandem Tubing	25/Case



RECEPTAL® Outer Shell Connection Tubing

RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum “T” of the RECEPTAL canisters between them for easy tandem connection.


REF	DESCRIPTION	QTY/UNIT
CTM50B	1/4 in. x 20 in. (6.35 mm x 500 mm), Female Connector	100/Case
CTM20B	1/4 in. x 8 in. (6.35 mm x 200 mm), Female Connector	300/Case



Tubing Manifold

For connection of 2, 3 or 4 canisters.

REF	DESCRIPTION	QTY/UNIT
TMF001	Tubing Manifold	10/Case



Accessories

Elbow Adapters and Taper

Use with Patient Port.



REF	DESCRIPTION	QTY/UNIT
43318-01 ©	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
43316-01 ©	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
43494-25 ©	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
43473-25 ©	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case

RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.



REF	DESCRIPTION	QTY/UNIT
44300-05 ©	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	100/Case
43448-25 ©	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case

AMSure® 5-in-1 Straight Connector

AMSure 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS807 ©	5-in-1 Connector, Sterile	100/Case
AS807B	5-in-1 Connector, Non-sterile	500/Case

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case

Reusable Hardware

Roll Stands

REF	DESCRIPTION	QTY/UNIT
43474-01	Large head, stainless steel Adjustable Roll Stand, conductive casters, holds up to 12 <i>RECEPTAL</i> ® or 8 <i>EZE-VAC</i> ® Canisters, use with 1500 mL, 2000 mL and 3000 mL canisters	1/Case
43475-01	Small head, stainless steel Adjustable Roll Stand, conductive casters, holds multiple <i>RECEPTAL</i> and <i>EZE-VAC</i> Canisters	1/Case
43476-P2	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 1500 mL and 2000 mL canisters	1/Case
43476-P3	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 3000 mL canisters	1/Case



Single and In-Line Floor Stands

REF	DESCRIPTION	QTY/UNIT
43431-01	Single Canister Floor Stand with tip-resistant base, use with 1500 mL and 2000 mL <i>RECEPTAL</i> Canisters	1/Case
43432-01	Dual Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL <i>RECEPTAL</i> Canisters	1/Case
43447-01	In-Line Quad Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL <i>RECEPTAL</i> Canisters	1/Case



Reusable Hardware

RECEPTAL® Liner System Tandem Canister Supports

REF	DESCRIPTION	QTY/UNIT
43429-01*	For 1500 mL and 2000 mL canisters	1/Case
43439-01	For 3000 mL canisters	

RECEPTAL Liner System Single Canister Supports

REF	DESCRIPTION	QTY/UNIT
43425-05*	For 1000 mL canister, metal	5/Case
43424-10	For 1000 mL canister support, with mounting plate, metal	5/Case
43426-05	For 1500 mL and 2000 mL canisters, metal	5/Case
43438-10	For 3000 mL canisters, with mounting plate for vacuum regulator, metal	10/Case
43428-05*	For 1500 mL and 2000 mL canisters, with mounting plate for vacuum regulator, metal	5/Case
43436-05	For 3000 mL canisters, metal	5/Case
43426P	For 1000 mL, 1500 mL, and 2000 mL canisters, plastic	20/Case
43436P	For 3000 mL canisters, plastic	10/Case
43424-V	For 1000/1500 mL canister support ring with valve	10/Case
43428-V	For 2000 mL canister support ring with valve	10/Case
43438-V	For 3000 mL canister support ring with valve	10/Case

REF	DESCRIPTION	QTY/UNIT
43467-01	For EZE-VAC 2.0 800 mL and 1200 mL canisters, metal, small: ID=4.4 in. (11.7 cm), wall to outer edge = 8.5 in. (21.6 cm)	12/Case
43469-12*	For EZE-VAC 1200 mL, 2000 mL and 3300 mL canisters, metal, large: ID=5.6 in. (14.2 cm), wall to outer edge = 9.4 in. (23.9 cm)	12/Case
43469-V	For EZE-VAC 1200 mL, 2000 mL, and 3300 mL with valve	10/Case

*Shown on Wall Bracket (REF 43444-05)



Reusable Hardware

Canister Brackets

REF	DESCRIPTION	QTY/UNIT
43444-05	Wall channel bracket for mounting canister, stainless steel	5/Case
43444P	Wall channel bracket for mounting canister, polycarbonate	50/Case
44303-10	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
44304-05	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
RCB001	Metal Rail Bracket with Universal Holder Clamp	30/Case



Choose the right
iRECEPTAL system
for your needs.

	iRECEPTAL 2.5	iRECEPTAL 5.0	iRECEPTAL Mini
			
FEATURE/CAPACITY	High-volume capacity for demanding procedures	Versatile performance for everyday surgical use	Compact and mobile for space-sensitive settings

Capacity	20L/32L	22L	11L
Rover	✓	✓	✓
Docking Station	✓	✓	✓
Manifold Compatible	✓	✓	✓
Smoke Filter Compatible	—	✓	✓
IV Pole Option	✓	✓	✓
CE Version Available	✓	✓	✓
China Version Available	✓	✓	✓
Ideal Use Case	High-volume fluid collection and longer procedures	General surgical procedures and everyday use	Mobile/bedside procedures and limited space

iRECEPTAL 2.5

REF	DESCRIPTION	QTY/UNIT
IRR025-I	<i>iReceptal 2.5 Rover -32L (China/CE/US)</i>	1/Each
IRR025-II	<i>iReceptal 2.5 Rover -20L (China/CE/US)</i>	1/Each
IRD025	<i>iReceptal 2.5 Docking Station (China)</i>	1/Each
IRD025-02*	<i>iReceptal 2.5 Docking Station (CE)</i>	1/Each
IRD025-01*	<i>iReceptal 2.5 Docking Station(US)</i>	1/Each
IRM502-I	<i>iReceptal 5.0 & Mini & 2.5 Manifold (China)</i>	20 Pieces/Case
IRM500	<i>iReceptal 5.0 & Mini & 2.5 Manifold (CE/US)</i>	20 Pieces/Case

iRECEPTAL 5.0

REF	DESCRIPTION	QTY/UNIT
IRR502-30	<i>iReceptal 5.0 Rover (China)</i>	1/Each
IRR502-02	<i>iReceptal 5.0 Rover (CE)</i>	1/Each
IRR502-01	<i>iReceptal 5.0 Rover (US)</i>	1/Each
IRD502-I	<i>iReceptal 5.0 and Mini Docking Station (China)</i>	1/Each
IRD502-02	<i>iReceptal 5.0 and Mini Docking Station (CE)</i>	1/Each
IRD502-01	<i>iReceptal 5.0 and Mini Docking Station (US)</i>	1/Each
IRM502-I	<i>iReceptal 5.0, Mini and 2.5 Manifold (China)</i>	20 Pieces/Case
IRM500	<i>iReceptal 5.0, Mini and 2.5 Manifold (CE/US)</i>	20 Pieces/Case
IRF500	<i>iReceptal 5.0 and Mini Smoke Filter (China/CE/US)</i>	1/Each

iRECEPTAL Mini

REF	DESCRIPTION	QTY/UNIT
IRR552-10	<i>iReceptal Mini Rover (China)</i>	1/Each
IRR552-02 ☺	<i>iReceptal Mini Rover (CE)</i>	1/Each
IRR552-01	<i>iReceptal Mini Rover (US)</i>	1/Each
IRR551-10	<i>iReceptal Mini Rover without IV Pole (China)</i>	1/Each
IRR551-02 ☺	<i>iReceptal Mini Rover without IV Pole (CE)</i>	1/Each
IRR551-01	<i>iReceptal Mini Rover without IV Pole (US)</i>	1/Each
IRD502-I	<i>iReceptal 5.0 and Mini Docking Station (China)</i>	1/Each
IRD502-02 ☺	<i>iReceptal 5.0 and Mini Docking Station (CE)</i>	1/Each
IRD502-01	<i>iReceptal 5.0 and Mini Docking Station (US)</i>	1/Each
IRM502-I	<i>iReceptal 5.0, Mini and 2.5 Manifold (China)</i>	20 Pieces/Case
IRM500 ☺	<i>iReceptal 5.0, Mini and 2.5 Manifold (CE/US)</i>	20 Pieces/Case
IRF500 ☺	<i>iReceptal 5.0 and Mini Smoke Filter (China/CE/US)</i>	1/Each

AMSure®

RESPIRATORY THERAPY SYSTEMS

– Enhancing Respiratory Care



AMSure Respiratory Therapy Systems are designed to deliver optimal FiO₂ levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Respiratory Masks and Nasal Cannulas.....	70
Nebulizer Kits	73
Accessories.....	74
Aerosol Drainage Bag	74
Heat and Moisture Exchangers	75
Tracheostomy Clean and Care Trays	75
Sterile Water and Saline	76
Suction Catheters, Trays and Kits	77
Suction Tubing with Connectors	80
Suction Yankauers and Kits	81
Pre-Connected Yankauer and Tubing	81
SteriCadi™ Suction Device Holster System	82

AMSure® Respiratory Masks and Nasal Cannulas



AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

AMSure Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS74010	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case



AMSure® Respiratory Masks and Nasal Cannulas



AS75280

High Flow Nasal Oxygen Cannula

Amsino's *AMSure* high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most high flow nasal cannula devices
- Single patient use, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75280 ©	High Flow Nasal Oxygen Cannula	Small, Adult	10/Case
AS75285	High Flow Nasal Oxygen Cannula	Medium, Adult	
AS75282 ©	High Flow Nasal Oxygen Cannula	Large, Adult	

Nasal Oxygen Cannula

AMSure Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




AS75080

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75080 ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
AS75085 ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	

Cushion Nasal Oxygen Cannula

AMSure Cushion Nasal Oxygen Cannulas with soft foam ear covers provide low-to-medium flow oxygen delivery. Soft, pliable prongs are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75181 ©	Soft Nasal Oxygen Cannula	Curved Prongs with Cushion, Adult	25/Case

AMSure® Respiratory Masks and Nasal Cannulas

Soft Nasal Oxygen Cannula

AMSure Soft Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75180 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Adult	50/Case
AS75190 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Pediatric	
AS75199 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Infant	



AS75180

Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75010	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
AS75020	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	



AS75010

Pediatric Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75050	Pediatric Aerosol Mask	50/Case



AS75050

AMSure® Respiratory Masks and Nasal Cannulas




Tracheostomy Mask

AMSure Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060	Adult Tracheostomy Mask	50/Case



AS78010



AS78015

AMSure® Nebulizer Kits

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78015 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece	
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	
AS78025 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece, Reusable	

AMSure® Accessories

Respiratory Adapters

- Available in 15 mm I.D. and 22 mm I.D.
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS730100	Universal Cuff Adapter	22 mm ID Connector	50/Case
AS730200	Multi Adapter	15 mm ID / 22 mm OD Connector	



AS730100



AS730200

AMSure® Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- Pre-connected self-locking plastic chain secures the system to a stable object
- Convenient drain port allows for easy one-handed drainage
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75035	750 mL Aerosol Drainage Bag	50/Case



Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
HME001 ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
HME002 ©	Heat and Moisture Exchangers, Sterile	

AMSure® Tracheostomy Clean and Care Trays



AS861

AMSure Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS860 ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
AS861 ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	

**Manufactured
in USA**

AMSure® Sterile Water and Saline

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.

Sterile Water for Inhalation (Bottles)

- 500 mL and 1000 mL bottles Sterile Water, USP
- For use with nebulizers and humidifiers
- Hanging Spikeable bottle helps solve the problem of look-alike solution “bags” and prevents mishaps from occurring – Use with vented spikes only



AS1065

AS1000

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	—	500 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP with non-sterile screw top cap	
AS1000	AS1000E ©	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	

Unit Dose Vials for Inhalation

- Available in three sizes: 3 mL, 5 mL and 15 mL
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation
- Available in both sterile water and sterile sodium chloride (0.9%), USP



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case (10 Boxes – 100 per Box)
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case (3 Boxes – 48 per Box)
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	

Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

AMSure® Suction Catheters, Trays and Kits



AS361

Suction Catheters

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- Opposing staggered eyelets reduce the potential for blocked suction
- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	STRAIGHT PACK	COILED PACK	QTY/UNIT
AS361 ©	6 Fr	17½ in.	•		50/Case
AS361C ©	6 Fr			•	
AS362 ©	8 Fr		•		
AS362C ©	8 Fr			•	
AS363 ©	10 Fr	23 in.	•		
AS363C ©	10 Fr			•	
AS364 ©	12 Fr		•		
AS364C ©	12 Fr			•	
AS365 ©	14 Fr		•		
AS365C ©	14 Fr			•	
AS366 ©	16 Fr		•		
AS366C ©	16 Fr			•	
AS367 ©	18 Fr		•		
AS367C ©	18 Fr			•	

AMSure® Suction Catheters, Trays and Kits



Suction Catheter Kits

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. Kits include coiled packed suction catheters in various French sizes, a pop-up solution cup with sterile inner wrap, and one or two vinyl gloves.

- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	VINYL GLOVES		QTY/UNIT
			ONE GLOVE	TWO GLOVES	
AS370	5 Fr	17 ½ in.	•		50/Case
AS380	5 Fr			•	
AS371 ☺	6 Fr		•		
AS381 ☺	6 Fr			•	
AS372 ☺	8 Fr		•		
AS382 ☺	8 Fr			•	
AS373 ☺	10 Fr	23 in.	•		
AS383 ☺	10 Fr			•	
AS374 ☺	12 Fr		•		
AS384 ☺	12 Fr			•	
AS375 ☺	14 Fr		•		
AS385 ☺	14 Fr			•	
AS376 ☺	16 Fr		•		
AS386 ☺	16 Fr			•	
AS387 ☺	18 Fr			•	

AMSure® Suction Catheters, Trays and Kits



SCT10

Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- Opposing staggered eyelets reduce the potential for blocked suction
- DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients
- Whistle tip enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	DESCRIPTION	DELEE TIP	WHISTLE TIP	QTY/UNIT
SCT08	8 Fr	17 ½ in.	Suction Catheter, Tray and Vinyl Gloves	•		100/Case
SCT10	10 Fr	23 in.		•		
SCT12	12 Fr	23 in.			•	
SCT14	14 Fr	23 in.			•	
SCT16	16 Fr	23 in.			•	
SCT18	18 Fr	23 in.			•	

AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures.

- Flexible, low-coil memory allows for easy draping
- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 © (50/Case)	AS820B (100/Case)
	6 ft. (1.8 m)	AS821 © (50/Case)	AS821B (50/Case)
	10 ft. (3.0 m)	AS822 © (50/Case)	AS822B (50/Case)
	12 ft. (3.7 m)	AS823 © (20/Case)	AS823B (50/Case)
	20 ft. (6.1 m)	AS824 © (20/Case)	AS824B (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 © (50/Case)	AS825B (50/Case)
	10 ft. (3.0 m)	AS826 © (50/Case)	AS826B (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C © (50/Case)	—
	12 ft. (3.7 m)	AS827 © (20/Case)	AS827B (50/Case)
	20 ft. (6.1 m)	AS828 © (20/Case)	AS828B (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M © (50/Case)	AS825BM (50/Case)
	10 ft. (3.0 m)	AS826M © (50/Case)	AS826BM (50/Case)
	12 ft. (3.7 m)	AS827M © (20/Case)	AS827BM (50/Case)
	20 ft. (6.1 m)	AS828M © (20/Case)	AS828BM (50/Case)

AMSure® Suction Yankauers and Kits



AMSure Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			AS830 ④ 50/Case	AS830B 500/Case
	•		•	AS831 ④ 50/Case	AS831B 500/Case
		•		AS832 ④ 50/Case	AS832B 500/Case
		•	•	AS833 ④ 50/Case	AS833B 500/Case
Flexible	•			AS834 ④ 50/Case	—
	•		•	AS835 ④ 50/Case	—
		•		AS836 ④ 50/Case	—

AMSure® Pre-Connected Yankauer and Tubing

AMSure Pre-connected yankauer and tubing sets include a rigid non-vented bulbous tip yankauer pre-attached to suction tubing in a variety of sizes. Sterile.


REF	TUBING DIAMETER	TUBING LENGTH	QTY/UNIT
AS823P ④	3/16 in. (5 mm)	12 ft. (3.7m)	20/Case
AS824P ④	3/16 in. (5 mm)	20 ft. (6.1 m)	20/Case
AS825P ④	1/4 in. (6 mm)	6 ft. (1.8 m)	50/Case
AS826P ④	1/4 in. (6 mm)	10 ft. (3.0 m)	50/Case
AS827P ④	1/4 in. (6 mm)	12 ft. (3.7m)	20/Case
AS828P ④	1/4 in. (6 mm)	20 ft. (6.1 m)	20/Case

SteriCadi™ Suction Device Holster System

SteriCadi is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to support improved infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder, Sterile	30/Case





SURGICAL AND NURSING SOLUTIONS

Everyday Patient Care Solutions



AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

Specimen Containers.....	84
Vaginal Speculums	84
Disposable Umbilical Cord Clamp.....	84
Ear/Ulcer Syringes	85
Cleansing Enema Sets.....	85
Connectors.....	86

AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
AS340 ©		•		•		•	•	100/Case
AS341 ©			•	•	•	•	•	100/Case
AS343 ©			•				•	500/Case
AS345 ©	•			•			•	100/Case



AS340

AMSure® Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- Round-edge design helps to increase patient comfort
- Single-use disposable
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS032S	Small Vaginal Speculum	100/Case (10 Boxes – 10 per Box)
AS032M	Medium Vaginal Speculum	100/Case (10 Boxes – 10 per Box)
AS032L	Large Vaginal Speculum	100/Case (10 Boxes – 10 per Box)



AS032M

AMSure® Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- Safety lock is designed to decrease risk of accidental disconnection
- Grooved teeth help provide an effective grip
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UCC100	Umbilical Cord Clamp, sterile	50/Case





AMSure® Ear/Ulcer Syringes

AMSure Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	SYRINGE COLOR	FEATURES	QTY/UNIT
AS00501 ©	1 oz. syringe	Green	Non-sterile, Bulk	300/Case
AS00502 ©	2 oz. syringe	Blue	Non-sterile, Bulk	200/Case
AS00503 ©	3 oz. syringe	Green	Non-sterile, Bulk	150/Case
AS00501S ©	1 oz. syringe	Green	Sterile, Single Pack	50/Case
AS00502S ©	2 oz. syringe	Blue	Sterile, Single Pack	50/Case
AS00503S ©	3 oz. syringe	Green	Sterile, Single Pack	50/Case



AS330

AMSure® Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS333

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS330 ©	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case

AMSure® Connectors

5-in-1 Straight Connector

AMSure 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3 /16 in., ¼ in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS807 	5-in-1 Connector	Sterile	100/Case
AS807B	5-in-1 Connector	Non-sterile, Bulk	500/Case



AS807

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



AMSINO Kit Components

Designed for integration into medical device kits and procedure trays, Amsino's Kit Components offering includes a variety of sterile and non-sterile products to support efficient clinical preparation and delivery. With flexible packaging formats and dependable supply capabilities, these products help meet the evolving needs of kit assemblers and healthcare providers.



BZK Antiseptic

- Benzalkonium chloride (BZK) antiseptic applicators for gentle, effective skin cleansing
- Available in both wipe and swabstick configurations for procedural prep and general antiseptic

REF	DESCRIPTION	QTY/UNIT
KL19398	BZK Antiseptic Swabsticks, 4", 3-pack, Sterile	250/Case
KL19399	BZK Antiseptic Wipe, Non-Sterile	10,000/Case



Alcohol Prep Pads

- Saturated with 70% isopropyl alcohol for effective skin preparation
- Individually wrapped pads designed for procedural prep and general antiseptic use

REF	DESCRIPTION	QTY/UNIT
APP001	Alcohol Prep Pad, 70% Isopropyl Alcohol, Sterile	10,000/Case

Providone-Iodine

- Povidone-iodine antiseptic applicators for broad-spectrum skin preparation
- Offered in sterile and non-sterile swabsticks, as well as convenient liquid pouches for flexible clinical use

REF	DESCRIPTION	QTY/UNIT
KL19397	Povidone-Iodine Swabsticks, 3-pack, Non-Sterile	250/Case
KL19418	Povidone-Iodine Swabsticks, 3-pack, Sterile	250/Case
KL19425	3/4 oz. Povidone-Iodine Solution, Sterile	250/Case



Wound and Device Irrigation

- Easy open twist screw top container
- Used for device and wound irrigation

REF	DESCRIPTION	QTY/UNIT
AS398	100mL Sterile Saline (0.9%) for Device/Wound Irrigation	48/Case
AS399	100mL Sterile Water for Device/Wound Irrigation	48/Case



Prefilled Syringes – Enhanced Packaging



- Designed specifically for further processing and custom procedure kit assembly applications, with packaging engineered to withstand ethylene oxide (ETO) sterilization processes
- Enhanced hermetically sealed, non-breathable packaging provides excellent gas barrier properties to help maintain product integrity and stability throughout handling, sterilization, and storage
- Available in multiple configurations, including pre-assembled and space-saving component formats, to support flexible kit packing requirements and optimize packaging efficiency

REF	DESCRIPTION	QTY/UNIT
3705-1US	10mL 0.9% Sodium Chloride USP in a 10mL Syringe with Enhanced Packaging	460/Case
IVF1010SSK2HY	Pre-Filled Normal Saline Flush Syringe, 0.9% Sodium Chloride, For Further Processing	150/Case

Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606.....	13	43120S.....	55	ABL404.....	35	AS322.....	25
003201.....	13	43121S.....	55	AC3110.....	7	AS32200.....	25
003202.....	13	4312SB.....	57	AE0106.....	14	AS326.....	25
0G901-01.....	54	4312SP.....	57	AE0107.....	14	AS32600.....	25
0G902-01.....	54	4315SB.....	57	AE0108-HP.....	14	AS327100.....	26
0G903-01.....	54	4315SP.....	57	AE1106.....	14	AS327200.....	26
0G904-01.....	54	43160SB.....	57	AE1112.....	14	AS328.....	26
0G905-01.....	54	43202-01.....	53	AE1138.....	14	AS330.....	85
0L654-01.....	54	43202-820821.....	53	AE3107.....	13	AS332.....	25
0L655-01.....	54	43203-01.....	52	AE3108.....	13	AS333.....	85
0L656-01.....	54	43206-01.....	52	AE3108B.....	13	AS340.....	37, 84
0P507-01.....	54	4320SB.....	57	AFS100.....	8	AS341.....	37, 84
0P508-01.....	54	4320SP.....	57	AFS102.....	8	AS343.....	37, 84
1010001.....	7	43210-01.....	53	AFS107.....	14	AS345.....	37, 84
1010001B.....	9	43213-01.....	52	AFS109.....	14	AS3500.....	13
107201.....	7	43303-01.....	61	AFS208.....	8	AS3515.....	13
108301.....	7	43303-02.....	61	APP001.....	87	AS361.....	77
108302.....	9	43303-05.....	61	AS-6240.....	42	AS361C.....	77
108303.....	7	43305-10.....	57, 61	AS-6250.....	42	AS362.....	77
108304.....	7	43306-01.....	57	AS-6260.....	42	AS362C.....	77
108305.....	8	43309-01.....	61	AS-6270.....	42	AS363.....	77
108306.....	7	4330SB.....	57	AS-6380.....	42	AS363C.....	77
108308.....	8	4330SP.....	57	AS-6390.....	42	AS364.....	77
108309.....	8	43316-01.....	63	AS-6481.....	42	AS364C.....	77
108310.....	7	43317-25.....	62	AS-6491.....	42	AS365.....	77
109602.....	7	43318-01.....	63	AS00501.....	85	AS365C.....	77
1510501.....	9	43423-01.....	56	AS00501S.....	85	AS366.....	77
158305.....	9	43424-10.....	65	AS00502.....	85	AS366C.....	77
158306.....	9	43424-V.....	65	AS00502S.....	85	AS367.....	77
158309.....	9	43425-05.....	65	AS00503.....	85	AS367C.....	77
159602.....	9	43426-05.....	65	AS00503S.....	85	AS370.....	78
3215.....	48	43426P.....	65	AS010.....	4, 37	AS371.....	78
3225.....	48	43428-05.....	65	AS011.....	4, 37	AS372.....	78
3255.....	48	43428-V.....	65	AS011P.....	4, 37	AS373.....	78
37043US.....	45	43429-01.....	65	AS015.....	37	AS374.....	78
3704US.....	45	43431-01.....	64	AS015CA.....	37	AS375.....	78
3705-1US.....	46, 88	43432-01.....	64	AS0159.....	42	AS376.....	78
3705-3US1.....	45	43436-05.....	65	AS016NS.....	37	AS380.....	78
37053US.....	45	43436P.....	65	AS016NSCA.....	37	AS381.....	78
37055US.....	45	43438-10.....	65	AS032L.....	84	AS382.....	78
3705CUS.....	45	43438-V.....	65	AS032M.....	84	AS383.....	78
3706US.....	45	43439-01.....	65	AS032S.....	84	AS384.....	78
38143.....	47	43444-05.....	66	AS0352.....	40, 76	AS385.....	78
38143-1.....	47	43444P.....	66	AS0352E.....	40, 76	AS386.....	78
38243.....	47	43445-01.....	56	AS0359.....	42	AS387.....	78
38243-1.....	47	43447-01.....	64	AS0552.....	40, 76	AS398.....	41
3844.....	48	43448-25.....	63	AS0552E.....	40, 76	AS398.....	88
38443.....	48	43449-05.....	56	AS0559.....	42	AS399.....	41
38443-1.....	48	43455-01.....	62	AS1000.....	40, 76	AS399.....	88
38453.....	48	43457-25.....	62	AS1065.....	40, 76	AS40008.....	28
38455.....	48	43467-01.....	65	AS1065E.....	40, 76	AS40008S.....	29
38543.....	47	43469-12.....	65	AS115.....	4, 37	AS40010.....	28
38543-1.....	47	43469-V.....	65	AS115CA.....	4	AS40010S.....	29
3854E1.....	47	43473-25.....	63	AS116.....	4, 37	AS409.....	38
38553.....	47	43474-01.....	64	AS116CA.....	4	AS41012.....	28
38555.....	47	43475-01.....	64	AS121.....	5, 36	AS41012HC.....	27
43010-05.....	54	43476-P2.....	64	AS126.....	5, 36	AS41012S.....	29
43010S.....	54	43476-P3.....	64	AS127.....	5, 36	AS41014.....	28
43015S.....	54	43494-25.....	63	AS130.....	5, 36	AS41014H.....	27
43020S.....	54	43503-01.....	56	AS136.....	5, 36	AS41014HC.....	27
43025-01.....	54	44300-05.....	63	AS1565.....	40, 76	AS41014S.....	29
43027-01.....	54	44303-10.....	66	AS302.....	25	AS41016.....	28
43030S.....	54	44304-05.....	66	AS306N.....	38	AS41016H.....	27
43041-01.....	54	607201.....	11	AS307N.....	38	AS41016HC.....	27
43042-01.....	54	608301.....	11	AS308N.....	38	AS41016S.....	29
43044-01.....	54	608304.....	11	AS309N.....	38	AS41018.....	28
43044-05.....	54	608305.....	11	AS3100.....	10	AS41018H.....	27
4305SB.....	57	608306.....	11	AS3101.....	10	AS41018HC.....	27
4305SP.....	57	608310.....	11	AS3102.....	10	AS41018S.....	29
43110S.....	55	AA3101.....	9	AS311.....	38	AS41020.....	28
43115S.....	55	ABL400.....	35	AS3110A.....	10	AS41020HC.....	27
43116S.....	55	ABL402.....	35	AS312.....	25	AS41020S.....	29

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
AS41022	28	AS78015	73	AS861610C	33	ASR0059E	41, 76
AS41022HC	27	AS78020	73	AS861612	33	ASR0159	41, 76
AS41022S	29	AS78025	73	AS861612C	33	ASR0159E	41, 76
AS41024	28	AS807	63, 86	AS861614	33	AY0100	15
AS41024HC	27	AS807B	63, 86	AS861614C	33	AY0200	15
AS41024S	29	AS820	58, 80	AS861616	33	AY0206	15
AS41026	28	AS820B	58, 80	AS861616C	33	CTM20B	62
AS41026S	29	AS821	58, 80	AS861618	33	CTM50B	62
AS41028	28	AS821B	58, 80	AS861618C	33	D38301	12
AS41030	28	AS822	58, 80	AS870	30	D38302	12
AS42012	28	AS822B	58, 80	AS87114	30	D38303	12
AS42014	28	AS823	58, 80	AS880	24	D38304	12
AS42014H	27	AS823B	58, 80	AS882	24	D38307	12
AS42016	28	AS823P	59, 81	AS886	24	D38310	12
AS42016H	27	AS824	58, 80	AS890	24	E-1004	3
AS42016HC	27	AS824B	58, 80	AS892	24	E-1200	3
AS42016S	29	AS824P	59, 81	AS89214	31	E-1203	3
AS42018	28	AS825	58, 80	AS89214R	31	ENF1004K	3
AS42018H	27	AS825B	58, 80	AS89214RK	31	ENF1200K	3
AS42018HC	27	AS825BM	58, 80	AS89216RK	31	ENF1203K	3
AS42018S	29	AS825M	58, 80	AS89300	20	ENF1203KD	3
AS42020	28	AS825P	59, 81	AS89314	18	ENF1203KLD	3
AS42020HC	27	AS826	58, 80	AS89314H	18	ENF500K	3
AS42020S	29	AS826B	58, 80	AS89314S	18	ENK127	5
AS42022	28	AS826BM	58, 80	AS89316	18	ENK136	5
AS42022HC	27	AS826C	58, 80	AS89316H	18	ENS015	4
AS42022S	29	AS826M	58, 80	AS89316S	18	ENS015NS	4
AS42024	28	AS826P	59, 81	AS89318	18	ENS115	4
AS42024HC	27	AS827	58, 80	AS89318H	18	ENS115CA	4
AS42024S	29	AS827B	58, 80	AS89318S	18	ENS115NS	4
AS42026	28	AS827BM	58, 80	AS89330	20	ENS115NSCA	4
AS42028	28	AS827M	58, 80	AS89400	20	ENS115T	4
AS42030	28	AS827P	59, 81	AS89414	19	FOL1118	30
AS43016	28	AS828	58, 80	AS89414H	19	FP092	16
AS43018	28	AS828B	58, 80	AS89414S	19	FP092Y1	16
AS43018S	29	AS828BM	58, 80	AS89416	19	FP100BY2	16
AS43020	28	AS828M	58, 80	AS89416H	19	FP100FY12Y1	16
AS43020S	29	AS828P	59, 81	AS89416S	19	FP100FY1	16
AS43022	28	AS829B-01	58, 80	AS89418	19	FP106FY2	16
AS43022S	29	AS830	59, 81	AS89418H	19	FP109BY3	16
AS43024	28	AS830B	59, 81	AS89418S	19	HME001	75
AS43024S	29	AS831	59, 81	AS89430	20	HME002	75
AS43026	28	AS831B	59, 81	AS896	24	iRD025	68
AS43026S	29	AS832	59, 81	AS89800	23	iRD025-01	68
AS44008	34	AS832B	59, 81	AS89814	21	iRD025-02	68
AS44010	34	AS833	59, 81	AS89814H	21	iRD502-01	68
AS44012	34	AS833B	59, 81	AS89814S	21	iRD502-02	68
AS44014	34	AS834	59, 81	AS89816	21	iRD502-I	68
AS44016	34	AS835	59, 81	AS89816H	21	iRF500	68
AS44018	34	AS836	59, 81	AS89816S	21	iRM500	68
AS44020	34	AS85008	32	AS89818	21	iRM502-I	68
AS44022	34	AS85008K	32	AS89818H	21	iRR025-I	68
AS44024	34	AS85010	32	AS89818S	21	iRR025-II	68
AS44028	34	AS85010K	32	AS89830	23	iRR502-01	68
AS730100	74	AS85012	32	AS89900	23	iRR502-02	68
AS730200	74	AS85012K	32	AS89914	22	iRR502-30	68
AS74010	70	AS85014	32	AS89914H	22	iRR551-01	68
AS75010	72	AS85014K	32	AS89914S	22	iRR551-02	68
AS75020	72	AS85016	32	AS89916	22	iRR551-10	68
AS75035	74	AS85016K	32	AS89916H	22	iRR552-01	68
AS75050	72	AS860	75	AS89916S	22	iRR552-02	68
AS75060	73	AS860608	33	AS89918	22	iRR552-10	68
AS75080	71	AS860610	33	AS89918H	22	IVE2500	13
AS75085	71	AS860612	33	AS89918S	22	IVE2501	14
AS75180	72	AS860614	33	AS89930	23	IVF1003	44
AS75181	71	AS860616	33	ASR0035	40, 76	IVF1005	44
AS75190	72	AS861	75	ASR0035E	40, 76	IVF1010	44
AS75199	72	AS861006	33	ASR0039	41, 76	IVF1010PX	45
AS75280	71	AS861008	33	ASR0039E	41, 76	IVF1010SSK2HY	88
AS75282	71	AS861010	33	ASR0055	40, 76	IVF1010TM	44
AS75285	71	AS861608	33	ASR0055E	40, 76	IVF1203	43
AS78010	73	AS861610	33	ASR0059	41, 76	IVF1203SS	43

Index

REF PAGE

IVF1205	43
IVF1205SS.....	43
IVF1210	43
IVF1210SS.....	43
IVF2020	43
IVF2020SS.....	43
IVSK005	15
KL19397.....	88
KL19398.....	87
KL19399.....	87
KL19418.....	88
KL19425.....	88
KL20100.....	22, 44
KL20101.....	22, 44
MS360Y172.....	79
MS3615Y172.....	79
MS361Y172.....	79
MS392Y172.....	79
OM001	70
PE90B12T.....	2
PE90B5T.....	2
PE90ST.....	2
PF01.....	50
PF103.....	50
PF104.....	50
PF201.....	50
PF201-1	50
PF3022.....	50
PF3023.....	50
PF3025.....	50
PG001	2
PL01.....	50
PL10CC01	50
PL20CC01	50
PL30CC01	50
PL50CC01	50
PR10CC02	50
PR20CC01	50
PR30CC01	50
PR3CC01	50
PR50CC01	50
PR50CC02	50
PR5CC01	50
PRCAP01	50
RCB001	66
S-1000-SI.....	3
SCT08	79
SCT10	79
SCT12	79
SCT14	79
SCT16	79
SCT18	79
SIH0001	60, 82
TMF001	62
UCC100.....	84
UYC001	63, 86
UYC001S	63, 86

Notes



Amsino International, Ltd.

708 Corporate Center Drive
Pomona, CA 91768
USA

+1 800.632.6746

info@amsino.com

www.amsino.com

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required. CE mark in the catalog is only based on latest registration status. The CE status may be changed after the latest revision of the catalog.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO, ACCU-MEASURE, AMSafe, AMSafe-3, AMSHIELD, AMSure, AMSure HG, CitraFlow, CitraFlow PLUS, EZE-VAC, iRECEPTAL, MedXL, One Mask, PraxiFill, PraxiFlow, Praxiject, PraxiLabel, PraxiMed, Puggle, R-Polished, Receptal, Safgard, Sal-Sense, SENTINEL, SteriCadi, Sure-Lok, THE-POLE-SYRINGE and VAC-GARD are registered trademarks of Amsino International, Ltd. © Copyright 2026, Amsino International, Ltd. All rights reserved.