

# AMSINO GENERAL CATALOG

## CANADA





# Table of Contents

## AMSure® Enteral Feeding Systems ..... 1

ENFit® and Non-ENFit Sets.....	2
ENFit and Non-ENFit Syringes.....	3
ENFit and Non-ENFit Irrigation Kits and Trays.....	4

## AMSafe® Medication Delivery Systems ..... 5

Gravity Administration Sets.....	6
10 Drops/mL.....	6
15 Drops/mL.....	8
60 Drops/mL.....	9
AMSafe-3® Selectable Drip Chamber	
10, 15 & 60 drops/mL.....	10
Secondary Piggyback Sets.....	11
IV Extension Sets.....	11
Standard Extension Sets.....	11
Microbore Extension Sets.....	12
Flow Controller Extension Sets.....	12
Sure-Lok® Needle-Free Connector.....	13
Sure-Lok Needle-Free IV Pump Administration Sets.....	14

## AMSure Urological Care ..... 15

Foley Catheter Trays & Collection.....	16
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	16
Add-a-Foley Catheter Trays.....	17
Foley Insertion Trays.....	18
Foley Catheters.....	19
Silicone-Coated Latex Foley Catheters.....	19
100% Silicone Foley Catheters.....	20
Intermittent Urinary Catheters.....	21
Urethral Insertion Trays.....	21
Self-Catheterization Kits with R-Polished™ Eye Catheters.....	22
PVC Intermittent Urethral Catheter with R-Polished Eyes.....	23
Red Rubber Latex Urethral Catheter.....	24
Ancillary & Patient Care Products.....	24
Urinary Leg Bags.....	24
Extension Tube for Leg Bags.....	24
Irrigation Trays and Kits.....	25
Syringes.....	26
Specimen Containers.....	26

## Pre-Filled Syringes for Injection ..... 27

PreFilled Saline Flush Syringes.....	28
Ancillary Products.....	31
Self-Catheterization Kits with R-Polished™ Eye Catheters.....	31
PVC Intermittent Urethral Catheter with R-Polished Eyes.....	31

## Waste Fluid Management Systems ..... 33

EZE-VAC® Suction Canister System.....	34
EZE-VAC 2.0 Suction Canister System.....	35
RECEPTAL® Liners.....	36
RECEPTAL Twist Liners.....	37
RECEPTAL Outer Shells.....	38
ACCU-MEASURE® Critical Measurement Devices.....	39
RECEPTAL Solidifiers.....	39
AMSure Suction Tubing with Connectors.....	40
Suction Yankauers and Kits.....	41
Pre-Connected Yankauer and Tubing.....	41
Accessories.....	43
Reusable Hardware.....	46

## AMSure Respiratory Therapy Systems..... 49

Respiratory Masks and Nasal Cannulas.....	50
Nebulizer Kits.....	53
Accessories.....	54
Aerosol Drainage Bag.....	54
Heat and Moisture Exchangers.....	55
Tracheostomy Clean and Care Trays.....	55
Suction Catheters, Trays and Kits.....	56
Suction Tubing with Connectors.....	59
Suction Yankauers and Kits.....	60

## AMSure Surgical & Nursing Systems..... 61

Specimen Containers.....	62
Vaginal Speculums.....	62
Disposable Umbilical Cord Clamp.....	62
Ear/Ulcer Syringes.....	63
Cleansing Enema Sets.....	63
Connectors.....	64

## Index of REFs ..... 65

## Notes ..... 67





## Optimizing Nutritional Delivery

AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The systems are designed to provide best patient outcomes and includes:

ENFit & Non-ENFit Enteral Feeding Sets.....	2
ENFit & Non-ENFit Enteral Syringes .....	3
ENFit & Non-ENFit Irrigation Kits & Trays .....	4

# ENFit® and Non-ENFit Enteral Feeding Sets

AMSure® Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel™ enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free options, not manufactured with natural rubber latex



ENF1200K

REF	ENFit CONNECTOR	STANDARD BORE (DEHP FREE)	LARGE BORE (DEHP FREE)	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
				Spike Only No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/EntraFlo™	
E-1004				•				•	•	•	30/Case
E-1200						•		•	•	•	
E-1203						•	•				
ENF1004K ©	•			•				•	•	•	
ENF1200K ©	•					•		•	•	•	
ENF1203K ©	•					•	•				
ENF500K ©	•				•			•	•	•	

EntraFlo™ is a registered trademark of Medline Industries, Inc.  
Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.  
ENFit® is a registered trademark of GEDSA.

# ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

AS010



## Bulb Syringe

- Flexible soft bulb promotes one-handed administration

ENS015NS



## Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage

AS116



## Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use

ENS115



REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010 ©			•						•		50/Case
AS011 ©			•			•				•	50/Case
AS011P ©			•				•			•	50/Case
AS115 ©		•			•	•				•	50/Case
AS116 ©		•			•			•			30/Case
ENS015 ©	•			•		•				•	50/Case
ENS015NS	•			•				•			30/Case
ENS115 ©	•				•	•				•	50/Case
ENS115NS	•				•			•			30/Case

# ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



## Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	
AS127		•		•		•	•	
ENK127	•	•		•		•	•	

## Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	
ENK136	•	•	•		•	•	•	•	



## MEDICATION DELIVERY SYSTEMS



## Safety First

The AMSafe IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:


Gravity Administration IV Sets.....	6
10 Drops/mL.....	6
15 Drops/mL.....	8
60 Drops/mL.....	9
AMSafe-3® Selectable Drip Chamber, 10, 15 & 60 drops/mL.....	10
Secondary Piggyback Sets .....	11
Standard Extension Sets .....	11
Microbore Extension Sets .....	12
Flow Controller Extension Sets.....	12
Sure-Lok® Needle-Free Connector .....	13
Sure-Lok® Needle-Free IV Pump Administration Sets .....	14



# Gravity Administration IV Sets








Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex



AC3110

## 10 Drops/mL General Primary Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108303	10	83 in. (211 cm)	Sure-Lok® Needle-free Y-site	16 mL	50/Case
					
108305	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
109602	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
					
108306	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
108308	10	83 in. (211 cm)	(2) Sure-Lok Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
					
AC3110	10	105 in. (266 cm)	(2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					

# Gravity Administration IV Sets


## 10 Drops/mL

### Flow Controller Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AFS102</b>	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok Needle-free Y-site	18.5 mL	50/Case



<b>AFS208 ©</b>	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) AMSafe® Needle-free Y-sites	19.2 mL	50/Case
-----------------	----	------------------	--	---------	---------



## 10 Drops/mL


### General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>1010001B</b>	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case







<b>108302</b>	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
---------------	----	-----------------	--	---------	---------



# Gravity Administration IV Sets

## 15 Drops/mL





### General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AA3101</b>	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
<b>158306</b>	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
<b>159602</b>	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
					
<b>158305</b>	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
					

## 15 Drops/mL




### Piggyback Sets with Back Check Valve


- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AS3100</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, Roller clamp	19.7 mL	50/Case
					
<b>AS3101</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
<b>AS3102</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
<b>AS3110A</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					

# Gravity Administration IV Sets

## 60 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					





# AMSafe-3® Selectable Drip Chamber

## 10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>D38301</b> ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site	27.7 mL	50/Case



<b>D38302</b> ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.2 mL	50/Case
-----------------	--------------	-----------------	---	---------	---------



<b>D38303</b> ©	10, 15 or 60	89 in. (226 cm)	Backcheck Valve, Split Septum Injection Site, Needle-Free Injection Site, Rotating Male Luer Lock Adapter	24 mL	50/Case
-----------------	--------------	-----------------	---	-------	---------



<b>D38304</b> ©	10, 15 or 60	103 in. (262 cm)	Backcheck Valve, Pre-pierced Y-site with Pre-attached EXT Set with Pinch Clamp, Sure-lok Needle-free Y-site, Pre-pierced Y-site	34 mL	50/Case
-----------------	--------------	------------------	---	-------	---------



# Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AS3500</b> ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case



AE3108

## IV Extension Sets

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.

### Standard Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>000606</b>	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok® Needle-free Y-site, Slide Clamp	2.5 mL	50/Case



<b>AE3107</b> ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
-----------------	---------------	--	---------	---------



<b>AE3108</b> ©	8 in. (20 cm)	Removable Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
-----------------	---------------	---	--------	----------



<b>003201</b>	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
---------------	----------------	---	--------	---------







<b>003202</b> ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
-----------------	----------------	--	--------	---------




# IV Extension Sets

## Microbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AE0106</b> ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
<b>AE1106</b> ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
				
<b>AE0107</b> ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok Needle-free Y-site	0.9 mL	50/Case
				
<b>AE1112</b>	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
				

## Flow Controller Extension Sets

Amsino's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AFS107</b> ©	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
				

## Sure-Lok® Needle-Free Connector



Our AMSafe® needle free valves are manufactured to provide both clinicians and patients save procedural outcomes.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AY0200 ©	Sure-Lok Needle-free Connector	0.33 mL	100/Case











# IV Pump Administration Sets

## Sure-Lok® Needle-Free IV Pump Administration Sets

- FDA cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>FP092-CA</b>	10	92 in. (234 cm)	Roller Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
<b>FP092Y1-CA</b>	10	92 in. (234 cm)	1 Needle-free Valve, Roller Clamp, 1 Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
<b>FP100BY2-CA</b>	10	100 in. (254 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
<b>FP109BY3-CA</b>	10	109 in. (277 cm)	Back check Valve, 3 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
<b>FP100FY1-CA</b>	10	100 in. (254 cm)	1 Needle-free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
<b>FP106FY2-CA</b>	10	106 in. (269 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					



## Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

<b>Foley Catheter Trays &amp; Collection .....</b>	<b>16</b>
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	16
Add-a-Foley Catheter Trays .....	17
Foley Insertion Trays.....	18
<b>Foley Catheters .....</b>	<b>19</b>
Silicone-Coated Latex Foley Catheters .....	19
100% Silicone Foley Catheters .....	20
<b>Intermittent Urinary Catheters .....</b>	<b>21</b>
Urethral Insertion Trays .....	21
Self-Catheterization Kits with R-Polished™ Eye Catheters .....	22
PVC Intermittent Urethral Catheter with R-Polished™ Eyes.....	23
Red Rubber Latex Urethral Catheter .....	24
<b>Ancillary &amp; Patient Care Products.....</b>	<b>24</b>
Urinary Leg Bags .....	24
Extension Tube for Leg Bags .....	24
Irrigation Trays and Kits .....	25
Syringes .....	26
Specimen Containers .....	26

# Foley Catheter Trays & Collection



AS89316

## AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

**Trays include:** 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S ©	14 Fr	5cc	10/Case
	AS89316S ©	16 Fr		
	AS89318S ©	18 Fr		
Silicone Coated Latex (HCPCS Code A4314)	AS89314 ©	14 Fr	5cc	10/Case
	AS89316 ©	16 Fr		
	AS89318 ©	18 Fr		

# Foley Catheter Trays & Collection



AS89300

## AMSure® Add-a-Foley Catheter Tray

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10cc	10/Case



# Foley Catheter Trays & Collection



## AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or (6) Povidone-iodine swab sticks, 1 prefilled saline syringe, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

REF	PRE-FILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS880 ©	10cc	•		20/Case
AS890 ©	30cc	•		
AS886 ©	10cc		•	
AS896 ©	30cc		•	

## AMSure Foley Insertion Trays with FlexGRIP®

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), 1 prefilled saline syringe, lubrication jelly, fenestrated drape, 1 FlexGRIP Securement Device. SteriGuide IFU.

REF	PRE-FILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS882 ©	10cc	•		20/Case
AS892 ©	30cc	•		

# Foley Catheters

## AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
3cc	AS40008 ©	8Fr	10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr		
5cc	AS41012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS41014 ©	14 Fr		
	AS41016 ©	16 Fr		
	AS41018 ©	18 Fr		
	AS41020 ©	20 Fr		
	AS41022 ©	22 Fr		
	AS41024 ©	24 Fr		
	AS41026 ©	26 Fr		
	AS41028 ©	28 Fr		
	AS41030 ©	30 Fr		
30cc	AS42012 ©	12 Fr	16 in. (40.6 cm)	10/Box
	AS42014 ©	14 Fr		
	AS42016 ©	16 Fr		
	AS42018 ©	18 Fr		
	AS42020 ©	20 Fr		
	AS42022 ©	22 Fr		
	AS42024 ©	24 Fr		
	AS42026 ©	26 Fr		
	AS42028 ©	28 Fr		
	AS42030 ©	30 Fr		

## 3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	AS43016 ©	16 Fr	16 in. (40.6 cm)	10/Box
	AS43018 ©	18 Fr		
	AS43020 ©	20 Fr		
	AS43022 ©	22 Fr		
	AS43024 ©	24 Fr		
	AS43026 ©	26 Fr		

# Foley Catheters

## AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



### Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
3cc	<b>AS40008S</b> ©	8Fr	10.6 in. (27cm)	10/Box
	<b>AS40010S</b> ©	10 Fr		
5cc	<b>AS41012S</b> ©	12 Fr	16 in. (40.6 cm)	10/Box
	<b>AS41014S</b> ©	14 Fr		
	<b>AS41016S</b> ©	18 Fr		
	<b>AS41018S</b> ©	18 Fr		
	<b>AS41020S</b> ©	20 Fr		
	<b>AS41022S</b> ©	22 Fr		
	<b>AS41024S</b> ©	24 Fr		
	<b>AS41026S</b> ©	26 Fr		
30cc	<b>AS42016S</b> ©	16 Fr	16 in. (40.6 cm)	10/Box
	<b>AS42018S</b> ©	18 Fr		
	<b>AS42020S</b> ©	20 Fr		
	<b>AS42022S</b> ©	22 Fr		
	<b>AS42024S</b> ©	24 Fr		

### 3-way Catheters

BALLOON SIZE	REF	FOLEY SIZE	LENGTH	QTY/UNIT
30cc	<b>AS43018S</b> ©	18 Fr	16 in. (40.6 cm)	10/Box
	<b>AS43020S</b> ©	20 Fr		
	<b>AS43022S</b> ©	22 Fr		
	<b>AS43024S</b> ©	24 Fr		
	<b>AS43026S</b> ©	26 Fr		

# Intermittent Urinary Catheters



AS870

## AMSure® Urethral Insertion Trays

AMSure Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection

### PVC Catheter Tray

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex

### Red Rubber Catheter Tray

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
<b>AS870</b> ©	14 Fr PVC*	20/Case
<b>AS87114</b> ©	14 Fr Red Rubber Latex	

\*Not manufactured with natural rubber latex.



# Intermittent Urinary Catheters



## AMSure® Self-Catheterization Kits with R-POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, BPA-free, not manufactured with natural rubber latex

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

REF	SIZE	LENGTH	QTY/UNIT
AS85012 ©	12 Fr	16 in. (40.6 cm)	50/Case
AS85014 ©	14 Fr		
AS85016	16 Fr		

# Intermittent Urinary Catheters

## AMSure® PVC Intermittent Urethral Catheter with R>POLISHED™ Eyelets

AMSure PVC Intermittent Catheters with R-Polish technology use radio frequency to create ultra-smooth eyelets for optimal patient comfort. They are designed to provide a smooth surface in a catheter that is firm enough for insertion, yet soft enough to be slightly pliable. These urethral catheters are designed for intermittent or short-term use and are available in various French sizes, as well as female, pediatric and male lengths. A Coudé tip is also available to assist in more difficult catheterizations where strictures or other anomalies may be present in the urethra. (HCPCS code A4351)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- PVC material eliminates latex allergy concerns
- Clear shaft provides easy viewing as you measure
- Short-term catheter use helps decrease infections
- Individually packaged, sterile
- BPA-free, not manufactured with natural rubber latex



FEATURES	REF	SIZE	LENGTH	QTY/UNIT
Female	AS860608 ©	8 Fr	6 in. (15.24 cm)	50/Case
	AS860610 ©	10 Fr		
	AS860612 ©	12 Fr		
	AS860614 ©	14 Fr		
	AS860616 ©	16 Fr		
Pediatric	AS861008 ©	8 Fr	10 in. (25.40 cm)	50/Case
	AS861010 ©	10 Fr		
Male	AS861608 ©	8 Fr	16 in. (40.64 cm)	50/Case
	AS861610 ©	10 Fr		
	AS861612 ©	12 Fr		
	AS861614 ©	14 Fr		
	AS861616 ©	16 Fr		
	AS861618 ©	18 Fr		

# Intermittent Urinary Catheters

## AMSure® Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free

REF	SIZE	LENGTH	QTY/UNIT
AS44008 ©	8 Fr	11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr		
AS44012 ©	12 Fr	16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr		
AS44016 ©	16 Fr		
AS44018 ©	18 Fr		
AS44020 ©	20 Fr		
AS44022 ©	22 Fr		



## Ancillary & Patient Care Products

### Urinary Leg Bags (HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex

SIZE	REF	DRAIN PORT STYLE	QTY/UNIT
600 mL	AS306N ©	Push-Pull	48/Case
	AS308N ©	Twist-Turn	
900 mL	AS307N ©	Push-Pull	48/Case
	AS309N ©	Twist-Turn	



### Extension Tube for Leg Bags (HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS311 ©	18 in. (45.72 cm) Extension Tube	50/Case



# Ancillary & Patient Care Products



AS121



AS126



AS136

## Irrigation Trays and Kits (HCPCS code A4320)

AMSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need.

### Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

### Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector

## Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			
AS127	•		•	

## Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		

## Syringes (HCPCS code A4322)

### Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

### Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

### Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	<b>AS010</b> ©	Bulk	Non-Sterile			50/Case
	<b>AS011</b> ©	Poly Pouch	Sterile (EO)	•		50/Case
	<b>AS011P</b> ©	Form Fill Seal	Sterile (EO)	•		50/Case
Piston Thumb Control Ring	<b>AS015</b> ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	<b>AS016NS</b>	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	<b>AS115</b> ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	<b>AS116</b> ©	Resealable Pole Bag	Non-Sterile	•	•	30/Case

## AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.  
(HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
<b>AS340</b> ©		•		•		•	•	100/Case
<b>AS341</b> ©			•	•	•	•	•	100/Case
<b>AS343</b> ©			•				•	500/Case
<b>AS345</b> ©	•			•			•	100/Case



AS340





An **AMSINO** Company

## STERILE INJECTION SOLUTIONS



## Pre-Filled Syringes for Injection

Amsino Pre-Filled Syringes are designed for safe, efficient, and sterile line flushing vascular access. These syringes help streamline care and reduce contamination risk. The product portfolio includes:

Pre-Filled Saline Flush Syringe.....	28
Ancillary Products .....	31
Probe Covers .....	31
Pharmacy Products.....	31

# Pre-Filled Saline Flush Syringes



## Praxiject™ Saline Prefilled Syringes

Manufactured in Canada according to ISO Standards, current Good Manufacturing Practices and MedXL Quality Control requirements, MedXL saline prefilled syringes set the standard for quality and safety. These syringes are terminally sterilized and can be used in a sterile field. Saline prefilled syringes are latex and preservative free. MedXL saline prefilled syringes are secure, reliable and accurate.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
3704	Praxiject Saline, 5ml of 0.9% Sodium Chloride in a 5ml syringe (SF)	720/Case (6 Boxes – 120 per Box)
37043	Praxiject Saline, 3ml of 0.9% Sodium Chloride in a 5ml syringe (SF)	720/Case (6 Boxes – 120 per Box)
37055	Praxiject Saline, 5ml of 0.9% Sodium Chloride in a 10ml syringe (SF)	600/Case (6 Boxes – 100 per Box)
3706	Praxiject Saline, 20ml of 0.9% Sodium Chloride in a 20ml syringe (SF)	360/Case (6 Boxes – 60 per Box)



## Praxiject SP Prefilled Syringes

With their sterile path, the ready-to-use Praxiject SP syringes offer you a safe and economical flushing product. With a quality saline solution contained in a reliable and easy-to-use syringe, discover the many advantages that make Praxiject SP a product chosen by many of the largest hospital institutions.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
3705-3	Praxiject SP, 10ml of 0.9% Sodium Chloride Injection USP in a 10ml syringe. For flushing only.	180/Case (6 Boxes – 30 per Box)



## Praxiject SF Prefilled Syringes

Manufactured in Canada according to ISO Standards, current Good Manufacturing Practices and MedXL Quality Control requirements, MedXL saline prefilled syringes set the standard for quality and safety. These syringes are terminally sterilized and can be used in a sterile field. Saline prefilled syringes are latex and preservative free. MedXL saline prefilled syringes are secure, reliable and accurate.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
3705-4	Praxiject SF, 10ml of 0.9% Sodium Chloride in a 10ml syringe (SF)	500/Case (10 Boxes – 50 per Box)

# Pre-Filled Saline Flush Syringes

## CitraFlow™ 46.7% SF for Vascular Access Device Locking

CitraFlow 46.7% SF prefilled syringes are filled with an anticoagulant and antimicrobial 46.7% sodium citrate solution. CitraFlow 46.7% SF is a proven alternative to heparin to lock Vascular Access Devices. It avoids the risks associated with systemic heparinization and is safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia.



REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>38143</b>	CitraFlow 46.7% SF, 3ml 46.7% Sodium Citrate solution in 5ml syringe (SF)	900/Case (6 Boxes – 150 per Box)
<b>38143-1</b>	CitraFlow 46.7% SF, twinpack of two 3ml 46.7% Sodium Citrate solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

## CitraFlow 30% SF for Vascular Access Device Locking

CitraFlow 30% SF products are anticoagulant and antimicrobial 30% sodium citrate prefilled syringes. They are a proven alternative to heparin for Vascular Access Devices locking as they avoid the risks associated with systemic heparinization and are safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>38243</b>	CitraFlow 30% SF, 3ml 30% Sodium Citrate solution in 5ml syringe (SF)	900/Case (6 Boxes – 150 per Box)
<b>38243-1</b>	CitraFlow 30% SF, twinpack of two 3ml 30% Sodium Citrate solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

## CitraFlow 4% SF for Vascular Access Device Locking

CitraFlow™ 4% SF prefilled syringes are a proven alternative to heparin for locking indwelling catheters on patients. Sodium citrate avoids the risks associated with systemic heparinization and is safe for use in patients with suspected or confirmed Heparin Induced Thrombocytopenia

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>38543</b>	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in a 5ml syringe (SF)	900/Case (6 Boxes – 150/Box)
<b>38543-1</b>	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in 5ml syringe, twinpack (SF)	600/Case (6 Boxes – 100/Box)
<b>3854E1</b>	CitraFlow 4% SF, 5ml 4% Sodium Citrate solution USP in 5ml syringe (SF)	720/Case (6 Boxes – 120/Box)
<b>38553</b>	CitraFlow 4% SF, 3ml 4% Sodium Citrate solution USP in 10ml syringe (SF)	600/Case (6 Boxes – 100/Box)
<b>38555</b>	CitraFlow 4% SF, 5ml 4% Sodium Citrate solution USP in 10ml syringe (SF)	600/Case (6 Boxes – 100/Box)

# Pre-Filled Saline Flush Syringes



## CitraFlow PLUS™ Prefilled Syringes

CitraFlow PLUS prefilled syringes combine in one syringe the anticoagulant properties of sodium citrate PLUS the antimicrobial properties of ethanol. This solution combination prevents biofilm formation, reduces incidences of CRBSI (catheter related blood stream infections) as well as being a safe alternative to heparin.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>3844</b>	CitraFlow PLUS, 5ml 4% Sodium Citrate/30% Ethanol solution in a 5ml syringe (SF)	720/Case (6 Boxes – 120 per Box)
<b>38443</b>	CitraFlow PLUS, 3mL 4% Sodium Citrate/30% Ethanol solution in a 5mL syringe (SF)	900/Case (6 Boxes – 150 per Box)
<b>38443-1</b>	CitraFlow PLUS, twinpack of two 3ml 4% Sodium Citrate/30% Ethanol solution in 5ml syringe (SF)	600/Case (6 Boxes – 100 per Box)
<b>38453</b>	CitraFlow PLUS, 3ml 4% Sodium Citrate/30% Ethanol solution in a 10ml syringe (SF)	600/Case (6 Boxes – 100 per Box)
<b>38455</b>	CitraFlow PLUS, 5ml 4% Sodium Citrate/30% Ethanol solution in a 10ml syringe (SF)	600/Case (6 Boxes – 100 per Box)

## PraxiFlow™ Prefilled Syringe Kits for Vascular Access Device Locking

Each PraxiFlow kit consists of one Praxiject™ 0.9% NaCl syringe for Flush and one CitraFlow™ Sodium Citrate syringe for Lock of Vascular Access Devices. All the proven benefits of Praxiject and CitraFlow prefilled syringes in one convenient easy to use package.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>3215</b>	1 Praxiject 0.9% NaCl syringe (10ml in 10cc syringe), 1 CitraFlow 46.7% Sodium Citrate syringe (3ml in 5cc syringe) (SF)	450/Case (6 Boxes – 75/Box)
<b>3225</b>	1 Praxiject 0.9% NaCl syringe (10ml in 10cc syringe), 1 CitraFlow 30% Sodium Citrate syringe (3ml in 5cc syringe)	450/Case (6 Boxes – 75/Box)
<b>3255</b>	1 Praxiject 0.9% NaCl syringe (10ml in a 10cc syringe), 1 CitraFlow 4% Sodium Citrate syringe (3ml in 5cc syringe)	450/Case (6 Boxes – 75/Box)

## Intravenous Therapy Accessories

In order to meet healthcare professionals' needs, MedXL offers accessories that complete its line of intravenous therapy medical devices.



REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>6039BV</b>	Bidirectional valve with Male/Female connectors	200/Case (4 Boxes – 50 per Box)
<b>614035</b>	IV Administration Set. IV med rate control set for BD brand syringes	100/Box

# Ancillary Products



## Probe Covers

MedXL probe covers include a wide selection of products to fill almost any need. There are many different dimensions available whether you require a probe cover with Velcro, ultrasound gel or simply elastic bands, there is a good chance we have what you need. All probe cover products are sterile and sterile field compatible.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>1030</b>	Ultrasound probe covers 7in.x50in. retractable folding, 2 elastic bands – Sterile	25/Box
<b>1032</b>	Ultrasound probe covers 7in.x36in. retractable folding, 2 elastic bands – Sterile	50/Box
<b>1041</b>	Ultrasound probe covers 7in.x24in. retractable folding, 2 elastic bands, gel – Sterile	50/Box



## Pharmacy Products

MedXL pharmacy products is comprised of transfer sets, medication preparation sets and accessories that are practical, secure and reliable for medication preparation and reconstitution.

REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>3002</b>	Transfer Set with check valves, vented spike, clamp and T-connector	50/Box
<b>3004</b>	Transfer Set with 3-way high flow stopcock and vented spike	50/Box
<b>3006</b>	Transfer Spike, non vented spike with bidirectional valve	100/Box
<b>3021-1</b>	Transfer Spike, vented spike with swabbable valve	100/Box
<b>3040</b>	Medication administration kit for reconstitution and administration of BCG – 15 cm	50/Box
<b>3041</b>	Cystoset for intravesical medication administration – 9 cm	100/Box
<b>3046</b>	Transfer set with high flow tubing – 46 cm	100/Box
<b>3047</b>	Batching pump extension set with 0.2 micron capsule filter	5/Box
<b>3048</b>	Batching pump extension set for Fluidose oral liquid pump set – 167 cm	25/Box





# RECEPTAL®

## WASTE FLUID MANAGEMENT SYSTEMS



## Waste Fluid Management Solutions


Amsino is a leading global provider of suction fluid management systems. The breadth of Amsino's RECEPTAL suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

EZE-VAC® Suction Canister System .....	34
EZE-VAC 2.0 Suction Canister System.....	35
RECEPTAL Liners .....	36
RECEPTAL Twist Liners .....	37
RECEPTAL Outer Shells .....	38
ACCU-MEASURE® Critical Measurement Devices .....	39
RECEPTAL Solidifiers .....	39
AMSure® Suction Tubing with Connectors.....	40
Suction Yankauers and Kits .....	41
Pre-Connected Yankauer and Tubing.....	41
Accessories.....	43
Reusable Hardware .....	46

# EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01 ©	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01 ©	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01 ©	3300 mL Canister, Lid and 2 Elbows	12/Case

# EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
43210-01 ©	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID × 6 ft. (5 mm × 1.8 m) + 3/16 in ID × 18 in. (5 mm × 46 cm)	30/Case



# RECEPTAL® Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
1000	BULK PACKAGED		
	43010-05 ©	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	
1500	BULK PACKAGED		
	43025-01 ©	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	
	43015S ©	Liner with VAC-GARD and 45g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	
	0L654-01 ©	Liner with VAC-GARD and Pour Spout	
2000	BULK PACKAGED		
	43027-01 ©	Liner with Pour Spout,	50/Case
	43041-01 ©	Liner with Ortho Port	
	43042-01 ©	Liner with VAC-GARD	
	43044-05 ©	Liner with VAC-GARD and Pour Spout	
	43020S ©	Liner with VAC-GARD and 60g Solidifier	
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	
	0G905-01 ©	Liner with Ortho Port	
	0P507-01 ©	Liner with Pour Spout	
3000	BULK PACKAGED		
	0L656-01* ©	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	0L655-01* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	
	43030S* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	

\* Dual port allows for dual suctioning during procedures.  
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID



# RECEPTAL® Twist Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.


- Premeasured, preloaded and contact-free Solidifier Dispenser Cap
- Integrated solidifier eliminates the need for open ports
- All-in-one solidifier and liner solution reduces costs associated with extra inventory
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



43121S



43120S shown above  
in canister 43445-01  
(canister not included)



REF	DESCRIPTION	QTY/UNIT
<b>43110S</b>	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
<b>43115S</b> ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
<b>43116S</b> ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	
<b>43120S</b> ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	
<b>43121S</b> ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	

# RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex



Color-coded by size for easy selection.


REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	
43445-01 ©	2000 mL Outer Shell	
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

## ACCU-MEASURE® Critical Measurement Devices

ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL® ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

RECEPTAL  
ACCU-MEASURE  
Device



REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>43306-01</b> ©	250 mL Critical Measurement Canister	Tandem Tube	12/Case
<b>43305-10</b> ©	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case

## RECEPTAL Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- No shaking required
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4320SP

### RECEPTAL Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
<b>4305SP</b>	500 Solidifier Pouch, solidifies 500 mL	96/Case
<b>4312SP</b>	1200 Solidifier Pouch, solidifies 1200 mL	
<b>4315SP</b>	1500 Solidifier Pouch, solidifies 1500 mL	
<b>4320SP</b>	2000 Solidifier Pouch, solidifies 2000 mL	
<b>4330SP</b>	3000 Solidifier Pouch, solidifies 3000 mL	

### RECEPTAL Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
<b>4305SB</b>	500 Solidifier Bottle, solidifies 500 mL	96/Case
<b>4312SB</b>	1200 Solidifier Bottle, solidifies 1200 mL	
<b>4315SB</b>	1500 Solidifier Bottle, solidifies 1500 mL	
<b>4320SB</b>	2000 Solidifier Bottle, solidifies 2000 mL	
<b>4330SB</b>	3000 Solidifier Bottle, solidifies 3000 mL	
<b>43160SB</b>	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case



4315SB

4320SB

# AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm)  Green Connectors	18 in. (46 cm)	<b>AS820</b> ☺ (50/Case)	<b>AS820B</b> (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> ☺ (50/Case)	<b>AS821B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> ☺ (50/Case)	<b>AS822B</b> (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> ☺ (20/Case)	<b>AS823B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> ☺ (20/Case)	<b>AS824B</b> (50/Case)
1/4 in. (6 mm)  Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> ☺ (50/Case)	<b>AS825B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> ☺ (50/Case)	<b>AS826B</b> (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> ☺ (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> ☺ (20/Case)	<b>AS827B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> ☺ (20/Case)	<b>AS828B</b> (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> (1 Roll/Case)
9/32 in. (7 mm)  Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> ☺ (50/Case)	<b>AS825BM</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> ☺ (50/Case)	<b>AS826BM</b> (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> ☺ (20/Case)	<b>AS827BM</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> ☺ (20/Case)	<b>AS828BM</b> (50/Case)

## Suction Yankauers and Kits



AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			<b>AS830</b> ④ 50/Case	<b>AS830B</b> 500/Case
	•		•	<b>AS831</b> ④ 50/Case	<b>AS831B</b> 500/Case
		•		<b>AS832</b> ④ 50/Case	<b>AS832B</b> 500/Case
		•	•	<b>AS833</b> ④ 50/Case	<b>AS833B</b> 500/Case
Flexible	•			<b>AS834</b> ④ 50/Case	—
	•		•	<b>AS835</b> ④ 50/Case	—
		•		<b>AS836</b> ④ 50/Case	—

## Pre-Connected Yankauer and Tubing

AMSure Pre-connected yankauer and tubing sets include a rigid non-vented bulbous tip yankauer pre-attached to suction tubing in a variety of sizes. Sterile.

REF	TUBING DIAMETER	TUBING LENGTH	QTY/UNIT
<b>AS823P</b> ④	3/16 in. (5 mm)	12 ft. (3.7m)	20/Case
<b>AS824P</b> ④	3/16 in. (5 mm)	20 ft. (6.1 m)	20/Case
<b>AS825P</b> ④	1/4 in. (6 mm)	6 ft. (1.8 m)	50/Case
<b>AS826P</b> ④	1/4 in. (6 mm)	10 ft. (3.0 m)	50/Case
<b>AS827P</b> ④	1/4 in. (6 mm)	12 ft. (3.7m)	20/Case
<b>AS828P</b> ④	1/4 in. (6 mm)	20 ft. (6.1 m)	20/Case







## Accessories

### RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

43303-01

REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
<b>43303-01</b> ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
<b>43303-02</b> ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	
<b>43303-05</b> ©	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	

### External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>43309-01</b> ©	External Shut-Off Valve, ¼ in. (0.64 cm) Male x Female Elbow	10/Case

### RECEPTAL Accu-Measure® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
<b>43305-10</b> ©	Canister Attachment Bracket for 250 mL Canister	10/Case

# Accessories

## Suction Tubing

For permanent Outer Shell-to-vacuum source, or Outer Shell-to-Outer Shell connections.

REF	DESCRIPTION	QTY/UNIT
<b>43455-01</b> ©	100 ft. (30.48 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	1/Case
<b>43457-25</b>	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case



## EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.

- **Non-sterile, BPA-free, not manufactured with natural rubber latex**


REF	DESCRIPTION	QTY/UNIT
<b>43317-25</b>	18 in. (0.45 m) Tandem Tubing	25/Case



## RECEPTAL Outer Shell Connection Tubing

RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum "T" of the RECEPTAL canisters between them for easy tandem connection.


REF	DESCRIPTION	QTY/UNIT
<b>CTM50B</b> ©	1/4 in. × 20 in. (6.35 mm × 500 mm), Female Connector	100/Case
<b>CTM20B</b> ©	1/4 in. × 8 in. (6.35 mm × 200 mm), Female Connector	300/Case



## Tubing Manifold

For connection of 2, 3 or 4 canisters.

REF	DESCRIPTION	QTY/UNIT
<b>TMF001</b>	Tubing Manifold	10/Case



## Accessories

### RECEPTAL® Elbow Adapters and Taper

Use with Patient Port.



REF	DESCRIPTION	QTY/UNIT
<b>43318-01</b> ©	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
<b>43316-01</b> ©	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
<b>43494-25</b>	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
<b>43473-25</b>	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case

### RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.



REF	DESCRIPTION	QTY/UNIT
<b>44300-05</b> ©	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	100/Case
<b>43448-25</b>	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case

### 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
<b>AS807</b> ©	5-in-1 Connector, Sterile	100/Case

### 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
<b>UYC001S</b>	6-in-1 Universal Y Connector, Sterile	100/Case
<b>UYC001</b>	6-in-1 Universal Y Connector, Non-sterile	500/Case

# Reusable Hardware

## Roll Stands

REF	DESCRIPTION	QTY/UNIT
43474-01	Large head, stainless steel Adjustable Roll Stand, conductive casters, holds up to 12 RECEPTAL® or 8 EZE-VAC® Canisters, use with 1500 mL, 2000 mL and 3000 mL canisters	1/Case
43475-01	Small head, stainless steel Adjustable Roll Stand, conductive casters, holds multiple RECEPTAL and EZE-VAC Canisters	1/Case
43476-P2	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 1500 mL and 2000 mL canisters	1/Case
43476-P3	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 3000 mL canisters	1/Case



## Single and In-Line Floor Stands

REF	DESCRIPTION	QTY/UNIT
43431-01	Single Canister Floor Stand with tip-resistant base, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43432-01	Dual Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43447-01	In-Line Quad Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case



# Reusable Hardware

## RECEPTAL® Liner System Tandem Canister Supports



REF	DESCRIPTION	QTY/UNIT
<b>43429-01*</b>	For 1500 mL and 2000 mL canisters	1/Case
<b>43439-01</b>	For 3000 mL canisters	

## RECEPTAL Liner System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
<b>43425-05*</b>	For 1000 mL canister, metal	5/Case
<b>43424-10</b>	For 1000 mL canister support, with mounting plate, metal	5/Case
<b>43426-05</b>	For 1500 mL and 2000 mL canisters, metal	5/Case
<b>43438-10</b>	For 3000 mL canisters, with mounting plate for vacuum regulator, metal	10/Case
<b>43428-05*</b>	For 1500 mL and 2000 mL canisters, with mounting plate for vacuum regulator, metal	5/Case
<b>43436-05</b>	For 3000 mL canisters, metal	5/Case
<b>43426P</b>	For 1000 mL, 1500 mL, and 2000 mL canisters, plastic	20/Case
<b>43436P</b>	For 3000 mL canisters, plastic	10/Case
<b>43424-V</b>	For 1000/1500 mL canister support ring with valve	10/Case
<b>43428-V</b>	For 2000 mL canister support ring with valve	10/Case
<b>43438-V</b>	For 3000 mL canister support ring with valve	10/Case

## EZE-VAC® Rigid System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
<b>43467-01</b>	For EZE-VAC 2.0 800 mL and 1200 mL canisters, metal, small: ID=4.4 in. (11.7 cm), wall to outer edge = 8.5 in. (21.6 cm)	12/Case
<b>43469-12*</b>	For EZE-VAC 1200 mL, 2000 mL and 3300 mL canisters, metal, large: ID=5.6 in. (14.2 cm), wall to outer edge = 9.4 in. (23.9 cm)	12/Case
<b>43469-V</b>	For EZE-VAC 1200 mL, 2000 mL, and 3300 mL with valve	10/Case

\*Shown on Wall Bracket (REF 43444-05)

# Reusable Hardware

## Canister Brackets

REF	DESCRIPTION	QTY/UNIT
43444-05	Wall channel bracket for mounting canister, stainless steel	5/Case
43444P	Wall channel bracket for mounting canister, polycarbonate	50/Case
44303-10	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
44304-05	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
RCB001	Metal Rail Bracket with Universal Holder Clamp	30/Case







## Enhancing Oxygen Delivery

AMSure® Respiratory Therapy Systems are designed to deliver optimal FiO<sub>2</sub> levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Respiratory Masks and Nasal Cannulas.....	50
Nebulizer Kits .....	53
Accessories.....	54
Aerosol Drainage Bag .....	54
Heat and Moisture Exchangers .....	55
Tracheostomy Clean and Care Trays .....	55
Suction Catheters, Trays and Kits .....	56
Suction Tubing with Connectors .....	59
Suction Yankauers and Kits .....	60

# Respiratory Masks and Nasal Cannulas



## AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

### Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




REF	DESCRIPTION	QTY/UNIT
AS74010	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

# Respiratory Masks and Nasal Cannulas

## Nasal Oxygen Cannula

AMSure® Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75080

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75080</b> ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
<b>AS75085</b> ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	

## Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75010

REF	DESCRIPTION	QTY/UNIT
<b>AS75010</b>	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
<b>AS75020</b>	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	

## Pediatric Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75050

REF	DESCRIPTION	QTY/UNIT
<b>AS75050</b>	Pediatric Aerosol Mask	50/Case



# Respiratory Masks and Nasal Cannulas




## Tracheostomy Mask

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060	Adult Tracheostomy Mask	50/Case

## Nebulizer Kits



AS78010

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

### T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

### Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	

# Accessories

## Respiratory Adapters

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- **Non-sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS730100</b>	Universal Cuff Adapter	22 mm ID Connector	50/Case
<b>AS730200</b>	Multi Adapter	15 mm ID / 22 mm OD Connector	



AS730100

AS730200

## Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- **Pre-connected self-locking plastic chain secures the system to a stable object**
- **Convenient drain port allows for easy one-handed drainage**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS75035</b>	750 mL Aerosol Drainage Bag	50/Case





## Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
<b>HME001</b> ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
<b>HME002</b> ©	Heat and Moisture Exchangers, Sterile	

## Tracheostomy Clean and Care Trays



AS861

AMSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS860</b> ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
<b>AS861</b> ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	

# Suction Catheters, Trays and Kits



## Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- Opposing staggered eyelets reduce the potential for blocked suction
- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	SIZE	LENGTH	STRAIGHT PACK	COILED PACK	QTY/UNIT
AS361 ©	6 Fr	17½ in.	•		50/Case
AS361C ©	6 Fr			•	
AS362 ©	8 Fr		•		
AS362C ©	8 Fr			•	
AS363 ©	10 Fr	23 in.	•		
AS363C ©	10 Fr			•	
AS364 ©	12 Fr		•		
AS364C ©	12 Fr			•	
AS365 ©	14 Fr		•		
AS365C ©	14 Fr			•	
AS366 ©	16 Fr		•		
AS366C ©	16 Fr			•	
AS367 ©	18 Fr		•		
AS367C ©	18 Fr			•	

# Suction Catheters, Trays and Kits

## Suction Catheter Kits

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. Kits include coiled packed suction catheters in various French sizes, a pop-up solution cup with sterile inner wrap, and one or two vinyl gloves.

- Whistle tip port enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex



REF	SIZE	LENGTH	VINYL GLOVES		QTY/UNIT
			ONE GLOVE	TWO GLOVES	
AS370	5 Fr	17 ½ in.	•		50/Case
AS380	5 Fr			•	
AS371 ©	6 Fr		•		
AS381 ©	6 Fr			•	
AS372 ©	8 Fr		•		
AS382 ©	8 Fr			•	
AS373 ©	10 Fr	23 in.	•		
AS383 ©	10 Fr			•	
AS374 ©	12 Fr		•		
AS384 ©	12 Fr			•	
AS375 ©	14 Fr		•		
AS385 ©	14 Fr			•	
AS376 ©	16 Fr		•		
AS386 ©	16 Fr			•	
AS387 ©	18 Fr			•	

## Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- Opposing staggered eyelets reduce the potential for blocked suction
- DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients
- Whistle tip enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex



REF	SIZE	LENGTH	DESCRIPTION	DELEE TIP	WHISTLE TIP	QTY/UNIT
SCT08	8 Fr	17 ½ in.	Suction Catheter, Tray and Vinyl Gloves	•		100/Case
SCT10	10 Fr	23 in.		•		
SCT12	12 Fr	23 in.			•	
SCT14	14 Fr	23 in.			•	
SCT16	16 Fr	23 in.			•	
SCT18	18 Fr	23 in.			•	



# Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.


- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm)  Green Connectors	18 in. (46 cm)	<b>AS820</b> © (50/Case)	<b>AS820B</b> (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> © (50/Case)	<b>AS821B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> © (50/Case)	<b>AS822B</b> (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> © (20/Case)	<b>AS823B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> © (20/Case)	<b>AS824B</b> (50/Case)
1/4 in. (6 mm)  Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> © (50/Case)	<b>AS825B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> © (50/Case)	<b>AS826B</b> (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> © (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> © (20/Case)	<b>AS827B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> © (20/Case)	<b>AS828B</b> (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> (1 Roll/Case)
9/32 in. (7 mm)  Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> © (50/Case)	<b>AS825BM</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> © (50/Case)	<b>AS826BM</b> (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> © (20/Case)	<b>AS827BM</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> © (20/Case)	<b>AS828BM</b> (50/Case)

## Suction Yankauers and Kits



AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			<b>AS830</b> ④ 50/Case	<b>AS830B</b> 500/Case
	•		•	<b>AS831</b> ④ 50/Case	<b>AS831B</b> 500/Case
		•		<b>AS832</b> ④ 50/Case	<b>AS832B</b> 500/Case
		•	•	<b>AS833</b> ④ 50/Case	<b>AS833B</b> 500/Case
Flexible	•			<b>AS834</b> ④ 50/Case	—
	•		•	<b>AS835</b> ④ 50/Case	—
		•		<b>AS836</b> ④ 50/Case	—

## Pre-Connected Yankauer and Tubing

AMSure® Pre-connected yankauer and tubing sets include a rigid non-vented bulbous tip yankauer pre-attached to suction tubing in a variety of sizes. Sterile.

REF	TUBING DIAMETER	TUBING LENGTH	QTY/UNIT
<b>AS823P</b> ④	3/16 in. (5 mm)	12 ft. (3.7m)	20/Case
<b>AS824P</b> ④	3/16 in. (5 mm)	20 ft. (6.1 m)	20/Case
<b>AS825P</b> ④	1/4 in. (6 mm)	6 ft. (1.8 m)	50/Case
<b>AS826P</b> ④	1/4 in. (6 mm)	10 ft. (3.0 m)	50/Case
<b>AS827P</b> ④	1/4 in. (6 mm)	12 ft. (3.7m)	20/Case
<b>AS828P</b> ④	1/4 in. (6 mm)	20 ft. (6.1 m)	20/Case





## Everyday Patient Care Solutions

AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

Specimen Containers.....	62
Vaginal Speculums.....	62
Disposable Umbilical Cord Clamp.....	62
Ear/Ulcer Syringes .....	63
Cleansing Enema Sets.....	63
Connectors.....	64

## Specimen Containers

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	MOLDED PATIENT INFORMATION	QTY/UNIT
<b>AS340</b> 		•		•		•	•	100/Case
<b>AS341</b> 			•	•	•	•	•	100/Case
<b>AS343</b> 			•				•	500/Case
<b>AS345</b> 	•			•			•	100/Case



AS340

## Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- **Round-edge design helps to increase patient comfort**
- **Single-use disposable**
- **Non-sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS032S</b>	Small Vaginal Speculum	100/Case
<b>AS032M</b>	Medium Vaginal Speculum	100/Case
<b>AS032L</b>	Large Vaginal Speculum	100/Case



AS032M

## Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- **Safety lock is designed to decrease risk of accidental disconnection**
- **Grooved teeth help provide an effective grip**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>UCC100</b>	Umbilical Cord Clamp, sterile	50/Case





## Ear/Ulcer Syringes

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS00501</b> ©	1 oz. syringe	Non-sterile, Bulk	300/Case
<b>AS00502</b> ©	2 oz. syringe	Non-sterile, Bulk	200/Case
<b>AS00503</b> ©	3 oz. syringe	Non-sterile, Bulk	150/Case
<b>AS00501S</b> ©	1 oz. syringe	Sterile, Single Pack	50/Case
<b>AS00502S</b> ©	2 oz. syringe	Sterile, Single Pack	50/Case
<b>AS00503S</b> ©	3 oz. syringe	Sterile, Single Pack	50/Case



AS330

## Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

### Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

### Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS333

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS330</b> ©	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
<b>AS333</b>	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case

# Connectors

## 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3 /16 in., ¼ in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS807 	5-in-1 Connector	Sterile	100/Case



AS807

## 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



# Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606.....	11	4305SB.....	39	608305.....	9	AS380.....	57
003201.....	11	4305SP.....	39	608306.....	9	AS381.....	57
003202.....	11	43110S.....	37	614035.....	30	AS382.....	57
0G901-01.....	36	43115S.....	37	AA3101.....	8	AS383.....	57
0G902-01.....	36	43116S.....	37	AC3110.....	6	AS384.....	57
0G903-01.....	36	43120S.....	37	AE0106.....	12	AS385.....	57
0G904-01.....	36	43121S.....	37	AE0107.....	12	AS386.....	57
0G905-01.....	36	4312SB.....	39	AE1106.....	12	AS387.....	57
0L654-01.....	36	4312SP.....	39	AE1112.....	12	AS40008.....	19
0L655-01.....	36	4315SB.....	39	AE3107.....	11	AS40008S.....	20
0L656-01.....	36	4315SP.....	39	AE3108.....	11	AS40010.....	19
0P507-01.....	36	43160SB.....	39	AFS102.....	7	AS40010S.....	20
0P508-01.....	36	43202-01.....	35	AFS107.....	12	AS41012.....	19
1010001B.....	7	43202-820821.....	35	AFS208.....	7	AS41012S.....	20
1030.....	31	43203-01.....	34	AS00501.....	63	AS41014.....	19
1032.....	31	43206-01.....	34	AS00501S.....	63	AS41014S.....	20
1041.....	31	4320SB.....	39	AS00502.....	63	AS41016.....	19
107201.....	6	4320SP.....	39	AS00502S.....	63	AS41016S.....	20
108301.....	6	43210-01.....	35	AS00503.....	63	AS41018.....	19
108302.....	7	43213-01.....	34	AS00503S.....	63	AS41018S.....	20
108303.....	6	43303-01.....	43	AS010.....	3, 26	AS41020.....	19
108305.....	6	43303-02.....	43	AS011.....	3, 26	AS41020S.....	20
108306.....	6	43303-05.....	43	AS011P.....	3, 26	AS41022.....	19
108308.....	6	43305-10.....	39	AS015.....	26	AS41022S.....	20
109602.....	6	43305-10.....	43	AS032L.....	62	AS41024.....	19
158305.....	8	43306-01.....	39	AS032M.....	62	AS41024S.....	20
158306.....	8	43309-01.....	43	AS032S.....	62	AS41026.....	19
159602.....	8	4330SB.....	39	AS115.....	3, 26	AS41026S.....	20
3002.....	31	4330SP.....	39	AS116.....	3, 26	AS41028.....	19
3004.....	31	43316-01.....	45	AS121.....	4, 25	AS41030.....	19
3006.....	31	43317-25.....	44	AS126.....	4, 25	AS42012.....	19
3021-1.....	31	43318-01.....	45	AS127.....	4, 25	AS42014.....	19
3040.....	31	43423-01.....	38	AS130.....	4, 25	AS42016.....	19
3041.....	31	43424-10.....	47	AS136.....	4, 25	AS42016S.....	20
3046.....	31	43424-V.....	47	AS306N.....	24	AS42018.....	19
3047.....	31	43425-05.....	47	AS307N.....	24	AS42018S.....	20
3048.....	31	43426-05.....	47	AS308N.....	24	AS42020.....	19
3215.....	30	43426P.....	47	AS309N.....	24	AS42020S.....	20
3225.....	30	43428-05.....	47	AS3100.....	8	AS42022.....	19
3255.....	30	43428-V.....	47	AS3101.....	8	AS42022S.....	20
3704.....	28	43429-01.....	47	AS3102.....	8	AS42024.....	19
37043.....	28	43431-01.....	46	AS311.....	24	AS42024S.....	20
3705-3.....	28	43432-01.....	46	AS3110A.....	8	AS42026.....	19
3705-4.....	28	43436-05.....	47	AS330.....	63	AS42028.....	19
37055.....	28	43436P.....	47	AS333.....	63	AS42030.....	19
3706.....	28	43438-10.....	47	AS340.....	26, 62	AS43016.....	19
38143.....	29	43438-V.....	47	AS341.....	26, 62	AS43018.....	19
38143-1.....	29	43439-01.....	47	AS343.....	26, 62	AS43018S.....	20
38243.....	29	43444-05.....	48	AS345.....	26, 62	AS43020.....	19
38243-1.....	29	43444P.....	48	AS3500.....	11	AS43020S.....	20
3844.....	30	43445-01.....	38	AS361.....	56	AS43022.....	19
38443.....	30	43447-01.....	46	AS361C.....	56	AS43022S.....	20
38443-1.....	30	43448-25.....	45	AS362.....	56	AS43024.....	19
38453.....	30	43449-05.....	38	AS362C.....	56	AS43024S.....	20
38455.....	30	43455-01.....	44	AS363.....	56	AS43026.....	19
38543.....	29	43457-25.....	44	AS363C.....	56	AS43026S.....	20
38543-1.....	29	43467-01.....	47	AS364.....	56	AS44008.....	24
3854E1.....	29	43469-12.....	47	AS364C.....	56	AS44010.....	24
38553.....	29	43469-V.....	47	AS365.....	56	AS44012.....	24
38555.....	29	43473-25.....	45	AS365C.....	56	AS44014.....	24
43010-05.....	36	43474-01.....	46	AS366.....	56	AS44016.....	24
43010S.....	36	43475-01.....	46	AS366C.....	56	AS44018.....	24
43015S.....	36	43476-P2.....	46	AS367.....	56	AS44020.....	24
43020S.....	36	43476-P3.....	46	AS367C.....	56	AS44022.....	24
43025-01.....	36	43494-25.....	45	AS370.....	57	AS730100.....	54
43027-01.....	36	43503-01.....	38	AS371.....	57	AS730200.....	54
43030S.....	36	44300-05.....	45	AS372.....	57	AS74010.....	50
43041-01.....	36	44303-10.....	48	AS373.....	57	AS75010.....	51
43042-01.....	36	44304-05.....	48	AS374.....	57	AS75020.....	51
43044-01.....	36	6039BV.....	30	AS375.....	57	AS75035.....	54
43044-05.....	36	607201.....	9	AS376.....	57	AS75050.....	51

REF	PAGE	REF	PAGE
AS75060 .....	53	AS880 .....	18
AS75080 .....	51	AS882 .....	18
AS75085 .....	51	AS886 .....	18
AS78010 .....	53	AS890 .....	18
AS78020 .....	53	AS892 .....	18
AS807 .....	45, 64	AS89300 .....	17
AS820 .....	40, 59	AS89314 .....	16
AS820B .....	40, 59	AS89314S .....	16
AS821 .....	40, 59	AS89316 .....	16
AS821B .....	40, 59	AS89316S .....	16
AS822 .....	40, 59	AS89318 .....	16
AS822B .....	40, 59	AS89318S .....	16
AS823 .....	40, 59	AS896 .....	18
AS823B .....	40, 59	AY0200 .....	13
AS823P .....	41, 60	CTM20B .....	44
AS824 .....	40, 59	CTM50B .....	44
AS824B .....	40, 59	D38301 .....	10
AS824P .....	41, 60	D38302 .....	10
AS825 .....	40, 59	D38303 .....	10
AS825B .....	40, 59	D38304 .....	10
AS825BM .....	40, 59	E-1004 .....	2
AS825M .....	40, 59	E-1200 .....	2
AS825P .....	41, 60	E-1203 .....	2
AS826 .....	40, 59	ENF1004K .....	2
AS826B .....	40, 59	ENF1200K .....	2
AS826BM .....	40, 59	ENF1203K .....	2
AS826C .....	40, 59	ENF500K .....	2
AS826M .....	40, 59	ENK127 .....	4
AS826P .....	41, 60	ENK136 .....	4
AS827 .....	40, 59	ENS015 .....	3
AS827B .....	40, 59	ENS015NS .....	3
AS827BM .....	40, 59	ENS115 .....	3
AS827M .....	40, 59	ENS115NS .....	3
AS827P .....	41, 60	FP092-CA .....	14
AS828 .....	40, 59	FP092Y1-CA .....	14
AS828B .....	40, 59	FP100BY2-CA .....	14
AS828BM .....	40, 59	FP100FY1-CA .....	14
AS828M .....	40, 59	FP106FY2-CA .....	14
AS828P .....	41, 60	FP109BY3-CA .....	14
AS829B-01 .....	40, 59	HME001 .....	55
AS830 .....	41, 60	HME002 .....	55
AS830B .....	41, 60	OM001 .....	50
AS831 .....	41, 60	RCB001 .....	48
AS831B .....	41, 60	SCT08 .....	57
AS832 .....	41, 60	SCT10 .....	57
AS832B .....	41, 60	SCT12 .....	57
AS833 .....	41, 60	SCT14 .....	57
AS833B .....	41, 60	SCT16 .....	57
AS834 .....	41, 60	SCT18 .....	57
AS835 .....	41, 60	TMF001 .....	44
AS836 .....	41, 60	UCC100 .....	62
AS85012 .....	22	UYC001 .....	45, 64
AS85014 .....	22	UYC001S .....	45, 64
AS85016 .....	22		
AS860 .....	55		
AS860608 .....	23		
AS860610 .....	23		
AS860612 .....	23		
AS860614 .....	23		
AS860616 .....	23		
AS861 .....	55		
AS861008 .....	23		
AS861010 .....	23		
AS861608 .....	23		
AS861610 .....	23		
AS861612 .....	23		
AS861614 .....	23		
AS861616 .....	23		
AS861618 .....	23		
AS870 .....	21		
AS87114 .....	21		



# Notes





**Amsino International, Ltd.**

708 Corporate Center Drive  
Pomona, CA 91768  
USA

+1 800.632.6746

info@amsino.com

**[www.amsino.com](http://www.amsino.com)**

---

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required. CE mark in the catalog is only based on latest registration status. The CE status may be changed after the latest revision of the catalog.

---

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO, ACCU-MEASURE, AMSafe, AMSafe-3, AMShield, AMSure, AMSure HG, EZE-VAC, iRECEPTAL, MedXL, One Mask, Praxiject, Puggle, R-Polished, Receptal, Safgard, Sal-Sense, SENTINEL, SteriCadi, Sure-Lok, THE-POLE-SYRINGE and VAC-GARD are registered trademarks of Amsino International, Ltd. © Copyright 2025, Amsino International, Ltd. All rights reserved.