

Amsino Passive Heat and Moisture Exchanger (HME)

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.



Product Details

- ISO 9360 Testing – Pass.
- Classification: Class 1.
- ID Port Connector: 15 mm (standard).
- Moisture output: 25 mg H₂O/L @ 1000 mL.
- Dead space: 8 mL.
- Resistance: 2.0 cmH₂O @ 60 L/min.
- Cross Reference: Smith Medical #570016, Teleflex #41112, Vyaire #003014
- Packaging: 138 x 11" x 11", N.W. 3.5 lbs, G.W. 7.7 lbs.
- See packaging for more details.

Item Code	Description	Quantity/Unit
HME001	Heat and Moisture Exchanger, Non-Sterile	50/Case
HME002	Heat and Moisture Exchanger, Sterile	50/Case

Call Amsino sales reps for pricing, lead time and any other order information.
To view our complete product line, visit us at www.amsino.com



Advancing Healthcare Worldwide®

Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768
800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2020, Amsino International, Inc. All rights reserved.