

AMSure® Closed Suction Catheters

Adult Closed Suction Catheter Trays – 72-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

Pediatric Closed Suction Catheter Trays – 72-hour Pediatric Catheter from 5 FR to 8 FR Catheter Kits with MDI Port Adapter.

Infection Control


- Closed system to prevent cross infections in facility
- Clean easily with needle free port

Safety:

- Avoid interruptions of oxygen supply due to suction
- Safety lock to prevent un-expected suction
- Clear depth mark for accurate insertion
- R-polished catheter tip to reduce respiratory tract injury

Convenience:

- Reduce additional step to separate from ventilator before suction
- Flexible connection to ventilator without moving patient



Clear mark on Vacuum Control – Easy for nurses to know if suction is closed or open (safe for patients).



Graduation marks on the Catheter – Easy to identify the depth when catheter is inserted.



Smooth Catheter tip

Integrated closed structure – Effectively protects patients and medical staff.



Adult Closed Suction Catheter Trays

TYPE	REF	SIZE	LENGTH	ADAPTER TYPE		FLEX CONNECTOR		QTY/UNIT
				Double Swivel Elbow	T-Piece	15 mm x 22 mm	22 mm x 22 mm	
Endotracheal Kits	VS363RC124	10 Fr	21.3 in. (54 cm)	•		•		20/Case
	VS364RC124	12 Fr		•		•		
	VS365RC124	14 Fr		•		•		
	VS363TP124	10 Fr			•		•	
	VS364TP124	12 Fr			•		•	
	VS365TP124	14 Fr			•		•	
Tracheostomy Kits	VS363TC124	10 Fr	12 in. (30.5 cm)	•		•		
	VS364TC124	12 Fr		•		•		
	VS365TC124	14 Fr		•		•		
	VS363TTP124	10 Fr			•		•	
	VS364TTP124	12 Fr			•		•	
	VS365TTP124	14 Fr			•		•	

Pediatric Closed Suction Catheter Trays

REF	SIZE	LENGTH	DESCRIPTION	QTY/UNIT
VS360RC72	5 Fr	12 in. (30.5 cm)	Endotracheal Y-Adapters (2.0 and 2.5 mm)	20/Case
VS361RC72	6 Fr		Endotracheal Y-Adapters (2.5, 3.0 and 3.5 mm)	
VS3615RC72	7 Fr		Endotracheal Y-Adapters (3.0, 3.5 and 4.0 mm)	
VS362RC72	8 Fr		Endotracheal Y-Adapters (3.0, 3.5 and 4.0 mm)	



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 USA
+1 800.632.6746 | info@amsino.com

To view our complete product line,
visit us at www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Amsino and AMSure are registered trademarks of Amsino International, Inc. © Copyright 2025, Amsino International, Inc. All rights reserved.