

AMSure® Foley Insertion Trays



AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or (6) Povidone-iodine swab sticks, 1 prefilled saline syringe, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape.

REF	PRE-FILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS880 ©	10cc	•		20/Case
AS890 ©	30cc	•		
AS886 ©	10cc		•	
AS896 ©	30cc		•	

AMSure Foley Insertion Trays with FlexGRIP®

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), 1 prefilled saline syringe, lubrication jelly, fenestrated drape, 1 FlexGRIP Securement Device. SteriGuide IFU.

REF	PRE-FILLED SYRINGE	PVP SWAB STICKS		QTY/UNIT
		1 PACK OF 3	2 PACKS OF 3	
AS882 ©	10cc	•		20/Case
AS892 ©	30cc	•		



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 USA
+1 800.632.6746 | info@amsino.com

To view our complete product line, visit us at www.amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure is a registered trademark of Amsino International, Inc. © Copyright 2026, Amsino International, Inc. All rights reserved.