


INTRODUCING THE UPGRADED *AMSure*[®] Foley Kit with FlexGRIP[®]

Now designed for enhanced comfort, safety, and infection prevention.

- **New FlexGRIP Securement Device**
Improves patient comfort and reduces catheter-related complications.
- **Enhanced Skin Prep**
Now includes both PVP swabs and a BZK wipe for optimal antiseptic coverage.
- **Updated IFU**
Now aligned with the latest Infection Prevention best practices.



AMSure Pre-Connected Drainage Bag Foley Trays

Simplifying delivery of care while maintaining clinical excellence.



- **Pre-Attached Catheter with Tamper-Evident Seal**
Helps prevent disconnection and maintains a sterile fluid pathway.
- **All-in-One Tray**
Includes everything needed for aseptic catheter insertion, streamlining setup, and reducing waste.
- **Infection Control Focus**
Features sterile gloves, PVP swabsticks, a BZK wipe, and our AmSure SteriGuide IFU.
- **Needle-Free Sampling Port**
Minimizes the risk of needle stick injuries and enhances safety.
- **Multiple Catheter Options**
Available in silicone-coated latex, 100% silicone, and hydrophilic-coated latex for patient-specific needs.
- **Cost-Effective and Clinically Efficient**
Designed to optimize inventory and reduce unnecessary spending.

AMSINO
Advancing Healthcare Worldwide[®]

New AMSure® Foley Catheter Trays with FlexGRIP®

REF	DESCRIPTION	QTY/UNIT
AS89414	Foley Tray, 14Fr Foley Silicone-Coated Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89414H	Foley Tray, 14Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89414S	Foley Tray, 14Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89416	Foley Tray, 16Fr Foley Silicone-Coated Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89416H	Foley Tray, 16Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89416S	Foley Tray, 16Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89418	Foley Tray, 18Fr Foley Silicone-Coated Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89418H	Foley Tray, 18Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89418S	Foley Tray, 18Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89914H	Foley Tray, 14Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89914S	Foley Tray, 14Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89916H	Foley Tray, 16Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89916S	Foley Tray, 16Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89918H	Foley Tray, 18Fr Foley Hydrophilic Latex, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89918S	Foley Tray, 18Fr Foley 100% Silicone, 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89400	Add-a-Foley Tray with 10 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89430	Add-a-Foley Tray with 30 cc Pre-Filled Syringe with Sterile Water, 2000 mL Low Profile Urinary Drainage Bag, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89900	Add-a-Foley Tray with 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case
AS89930	Add-a-Foley Tray with 10 cc Pre-Filled Syringe with Sterile Water, Urine Meter, Waterproof and Fenestrated Drapes, Nitrile Powder-Free, 1 BZK pad, 1 pack of 3 PVP Swab Sticks, 1 FlexGRIP	10/Case

To view our complete product line,
visit us at www.amsino.com



Advancing Healthcare Worldwide®

Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA

+1 800.632.6746 | info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure is a registered trademark of Amsino Medical Group Co. Ltd. © Copyright 2025, Amsino Medical Group Co. Ltd. All rights reserved.