

# Urological Care Product Catalog





Amsino Medical Group offers a comprehensive range of urological products, for acute care, long-term care, and home care. Whether short-term or for chronic care, Amsino's urological care products are designed for optimal patient comfort, safety, and infection prevention.

# Table of Contents

<b>Foley Catheter Trays &amp; Collection</b> .....	<b>1</b>
Pre-Connected Foley Catheter Urine Drainage Bag Trays .....	1
Add-a-Foley Catheter Urine Drainage Bag Trays .....	2
Pre-Connected Foley Catheter Urine Meter Kits .....	3
Add-a-Foley Catheter Urine Meter Trays .....	4
Foley Insertion Trays .....	5
Urine Drainage Bags .....	7
Urine Meter Drainage Bags .....	9
<b>Foley Catheters</b> .....	<b>11</b>
Hydrophilic-Coated Latex .....	11
Silicone-Coated Latex .....	12
100% Silicone .....	13
<b>Intermittent Urinary Catheters</b> .....	<b>14</b>
Urethral Insertion Trays and Kits .....	14
Closed System Urethral Drainage Trays .....	15
Self-Catheterization Kits with R-Polished® Eyes .....	16
PVC Intermittent Catheters with R-Polished® Eyes .....	17
Red Rubber Latex Urethral Catheters .....	18
<b>Irrigation Sets</b> .....	<b>19</b>
<b>Ancillary &amp; Patient Care Products</b> .....	<b>20</b>
Leg Bags .....	20
Leg Bag Extension Tubes .....	20
Pediatric Urine Collector .....	21
Irrigation Trays & Kits .....	22
Syringes .....	23
Specimen Containers .....	24
<b>Index of REFs</b> .....	<b>25</b>



# Foley Catheter Trays & Collection



## AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure Pre-Connected Foley Catheter Urine Drainage Bag Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize contamination. **(HCPCS codes: See product chart below)**

- Tray includes all necessary components for aseptic catheter insertion
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

**Trays include:** 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5 cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag.

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S	14 Fr.	5 cc	10/Case
	AS89316S	16 Fr.	5 cc	10/Case
	AS89318S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89314	14 Fr.	5 cc	10/Case
	AS89316	16 Fr.	5 cc	10/Case
	AS89318	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr.	5 cc	10/Case
	AS89316H	16 Fr.	5 cc	10/Case
	AS89318H	18 Fr.	5 cc	10/Case



AS89300



## AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

**Trays include:** 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag.

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10 cc	10/Case
AS89330	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30 cc	10/Case

**BZK  
Swabstick  
Kits**

**COMING SOON!**

## AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs).

**(HCPCS codes: See product chart below)**

**Trays include:** 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, Sterile, not manufactured with rubber latex.



CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89814S	14 Fr.	5 cc	10/Case
	AS89816S	16 Fr.	5 cc	10/Case
	AS89818S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89814	14 Fr.	5 cc	10/Case
	AS89816	16 Fr.	5 cc	10/Case
	AS89818	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89814H	14 Fr.	5 cc	10/Case
	AS89816H	16 Fr.	5 cc	10/Case
	AS89818H	18 Fr.	5 cc	10/Case

## AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. **(HCPCS Code A4354)**

**Trays include:** 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, Sterile, not manufactured with rubber latex.



REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10 cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30 cc	10/Case



## AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. **(HCPCS code A4310)**

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

### Peel-Top

**Trays include:** 1000 mL collection basin with peel-off paper top web, (1-pack of 3) Povidone-iodine swab sticks or benzalkonium chloride prep pad (3 total), vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape.

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	AS880	10 cc	20/Case
	AS890	30 cc	20/Case
BZK Swab Sticks	AS880K	10 cc	20/Case
	AS890K	30 cc	20/Case

### CSR Wrapped

**Trays include:** 1000 mL collection basin, 500mm x 500mm CSR wrapped in header bag, (2-packs of 3) Povidone-iodine swab sticks (6 total), nitrile powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape.

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks (2 packs of 3)	AS886	10 cc	20/Case
	AS896	30 cc	20/Case





## Prevention of Catheter Associated Urinary Tract Infections


Prevention of Catheter Associated Urinary Tract Infections (CAUTIs) is one of the highest concerns whenever a patient is catheterized for long periods of time. In fact, the longer the period of catheterization, the higher the incidence rate of UTIs. According to the US Centers for Disease Control (CDC), UTIs account for over 30% of infections reported by acute care hospitals.<sup>1</sup> Virtually all of these are caused by introduction of catheters or other instruments into the urinary tract.

The use of a closed system with a sterile fluid path, such as Amsino's AMSure® Pre-Connected Foley Trays, aids in reducing CAUTIs by following the CDC Guidelines<sup>2</sup> in helping to eliminate breaks in aseptic technique and disconnections of the catheter from the drainage tube.


The AMSure closed Foley catheter trays are not only preconnected at the catheter-tubing junctions, but also include tamper-evident tape at the junction to aid in revealing if the system has been opened.

Furthermore, evidence points to the use of hydrophilic coated catheters as beneficial to patient in comfort, and helping to reduce trauma and infections.<sup>3,4</sup> Amsino's AMSure Hydrophilic-Coated Foley catheters utilize a proprietary coating method that offer superior lubricity.

See pages 1 and 3 for closed system AMSure Pre-Connected Foley Trays.



Before coating



After coating

Water drop on the surface

# AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube.

**(HCPCS code A4358)**


- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex




Anti-Reflux Drip Device



Anti-Reflux Large Chamber



Anti-Reflux Small Chamber



AMSafe® Needleless Sampling Port, Luer Lock or Luer Slip



Pre-Pierced Blunt Cannula Port



T-Tap Drain Port

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322	Teardrop	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS332	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS302	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	20/Case
AS32200	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600*	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	20/Case
AS326*	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	20/Case

\* 100% sterile packed in peel pouch.



AS312 Diamond



AS302 Pear



AS322 Teardrop



AS32200 Low-Profile




AS32600 Low-Profile



AS332 Teardrop



AS326 Oval



## AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. **(HCPCS code A4358)**

- Anti-Reflex Chamber helps improve patient safety by reducing urine backflow
- Twist valve drains at the bottom of the urine meter for easy drainage
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow delivering safer patient outcomes
- Sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS328	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case




# ACCU-MEASURE® Pediatric Urine Meter

## 2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)

- Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag
- Needle-free sampling port with luer lock access
- Leak-proof seams ensure fluid containment
- Meter bottle capacity 150 mL
- Adjustable bed sheet clamp for securing drainage bag tubing
- Sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS327100	ACCU-MEASURE Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case

# Foley Catheters

## AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



### Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41014H	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016H	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018H	18 Fr. × 16 in. (40.6 cm)	10/Box
30 cc	AS42014H	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS42016H	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS42018H	18 Fr. × 16 in. (40.6 cm)	10/Box

### Coude Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41012HC	12 Fr. × 16 in. (40.6 cm)	10/Box
	AS41014HC	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016HC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018HC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS41020HC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS41022HC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS41024HC	24 Fr. × 16 in. (40.6 cm)	10/Box
30 cc	AS42016HC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS42018HC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS42020HC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS42022HC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS42024HC	24 Fr. × 16 in. (40.6 cm)	10/Box

## AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
3 cc	AS40008	8 Fr. × 10.6 in. (27cm)	10/Box
	AS40010	10 Fr. × 10.6 in. (27cm)	10/Box
5 cc	AS41012	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026	26 Fr. × 16 in. (40.6cm)	10/Box
	AS41028	28 Fr. × 16 in. (40.6cm)	10/Box
	AS41030	30 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS42012	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020	20 Fr. × 16 in. (40.6cm)	10/Box
	AS42022	22 Fr. × 16 in. (40.6cm)	10/Box
	AS42024	24 Fr. × 16 in. (40.6cm)	10/Box
	AS42026	26 Fr. × 16 in. (40.6cm)	10/Box
	AS42028	28 Fr. × 16 in. (40.6cm)	10/Box
	AS42030	30 Fr. × 16 in. (40.6cm)	10/Box

# AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



## Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41012S	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014S	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016S	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018S	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020S	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022S	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024S	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026S	26 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS42016S	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018S	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020S	20 Fr. × 16 in. (40.6cm)	10/Box

## Coude Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41012SC	12 Fr. × 16 in. (40.6 cm)	10/Box
	AS41014SC	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016SC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018SC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS41020SC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS41022SC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS41024SC	24 Fr. × 16 in. (40.6 cm)	10/Box
	30 cc	AS42018SC	18 Fr. × 16 in. (40.6 cm)
AS42020SC		20 Fr. × 16 in. (40.6 cm)	10/Box
AS42022SC		22 Fr. × 16 in. (40.6 cm)	10/Box
AS42024SC		24 Fr. × 16 in. (40.6 cm)	10/Box



# Intermittent Urinary Catheters



AS870

## AMSure® Urethral Insertion Trays and Kits

AMSure Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection

### PVC Catheter Tray

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex

### Red Rubber Catheter Trays

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

**Trays include:** 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
AS870	14 Fr PVC*	20/Case
AS87114	14 Fr Red Rubber Latex	20/Case
AS87118	18 Fr Red Rubber Latex	20/Case

\*Not manufactured with natural rubber latex.

## AMSure® Closed System Urethral Drainage Kits

AMSure Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination.

(HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection
- Sterile

### Red Rubber Catheter Kits

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

### PVC Catheter Kits

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex



AS89214R

**Kit includes:** 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	AS89214R	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89214	14 Fr. × 16 in. (40.6cm)	PVC Latex- Free	20/Case
BZK Swab Sticks	AS89214RK	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89216RK	16 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case

# AMSure® Self-Catheterization Kits with R))POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. **(HCPCS code A4353)**

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, BPA-free, not manufactured with natural rubber latex

**Kit includes:** Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)




SWAB STICK	REF	CATHETER SIZE	QTY/UNIT
PVP Swab Sticks	AS85008	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016	16 Fr. x 16 in. (40.6cm)	50/Case
BZK Swab Sticks	AS85008K	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010K	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012K	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014K	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016K	16 Fr. x 16 in. (40.6cm)	50/Case

# AMSure® PVC Intermittent Urethral Catheter with R))POLISHED™ Eyes

AMSure PVC Intermittent Catheters are designed for short-term, intermittent use and are firm enough for easy insertion, yet soft enough to be slightly pliable. The AMSure Intermittent Catheters utilize Amsino's proprietary radio frequency R-Polish technology to create ultra-smooth eyelets and a smooth surface for optimal patient comfort. They are available in various French sizes as well as female, pediatric and male lengths.. (HCPCS code A4351)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Clear shaft provides easy viewing as you measure
- PVC material eliminates latex allergy concerns
- Designed to minimize trauma to help reduce risk of infections
- Individually packaged, sterile, BPA-free, not manufactured with natural rubber latex



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	AS860608	8 Fr, 6 in. (15.24 cm)	50/Case
	AS860610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS860612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS860614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS860616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS861006	6 Fr, 10 in. (25.40 cm)	50/Case
	AS861008	8 Fr, 10 in. (25.40 cm)	50/Case
	AS861010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS861608	8 Fr, 16 in. (40.64 cm)	50/Case
	AS861610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS861612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS861614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS861616	16 Fr, 16 in. (40.64 cm)	50/Case
	AS861618	18 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	AS861610C	10 Fr, 16" (40.65 cm)	30/Case
	AS861612C	12 Fr, 16" (40.65 cm)	30/Case
	AS861614C	14 Fr, 16" (40.65 cm)	30/Case
	AS861616C	16 Fr, 16" (40.65 cm)	30/Case
	AS861618C	18 Fr, 16" (40.65 cm)	30/Case





## AMSure® Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths.

**(HCPCS code A4351)**

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free


REF	DESCRIPTION	QTY/UNIT
AS44008	8 Fr, 11 in. (27.94 cm)	100/Case
AS44010	10 Fr, 11 in. (27.94 cm)	100/Case
AS44012	12 Fr, 16 in. (40.64 cm)	100/Case
AS44014	14 Fr, 16 in. (40.64 cm)	100/Case
AS44016	16 Fr, 16 in. (40.64 cm)	100/Case
AS44018	18 Fr, 16 in. (40.64 cm)	100/Case
AS44020	20 Fr, 16 in. (40.64 cm)	100/Case
AS44022	22 Fr, 16 in. (40.64 cm)	100/Case
AS44024	24 Fr, 16 in. (40.64 cm)	100/Case
AS44028	28 Fr, 16 in. (40.64 cm)	100/Case

# Irrigation Sets

## AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
ABL400	Single Lead	20/Case



REF	TYPE	QTY/UNIT
ABL402	Two Leads	20/Case



REF	TYPE	QTY/UNIT
ABL404	Four Leads	20/Case

# Ancillary & Patient Care Products

## Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex



AS307N

SIZE	REF	DRAIN PORT STYLE	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER & PROTECTIVE CAP	QTY/UNIT
600 mL	AS306N	Push-Pull		48/Case
	AS308N	Twist-Turn		48/Case
	LBE306N	Push-Pull	•	48/Case
	LBE308N	Twist-Turn	•	48/Case
900 mL	AS307N	Push-Pull		48/Case
	AS309N	Twist-Turn		48/Case
	LBE307N	Push-Pull	•	48/Case
	LBE309N	Twist-Turn	•	48/Case

## Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case



## Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS409	200 mL Pediatric Urine Bag	50/Box





# AMSure® Irrigation Trays and Kits

AMSure Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need.  
(HCPCS code A4320)

## Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector



AS121

## Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap



AS126

**Non-sterile irrigation kits include:** 500 mL graduated container, patient ID label and tip protector.

**Sterile irrigation trays include:** 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector.

## Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

## Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130		•	20/Case
AS136	•		20/Case



AS136

## Syringes

(HCPCS code A4322)

### Bulb Syringe


- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

### Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

### Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS011P	Form Fill & Seal	Sterile (EO)	•	•	50/Case
	AS011	Poly Pouch	Sterile (EO)	•		50/Case
	AS010	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	AS015	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	AS115	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116	Resealable Pole Bag	Non-Sterile	•	•	30/Case

## AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples. **(HCPCS code P9615)**

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	PACKAGING	STERILITY	LABELS	QTY/UNIT
AS340	Poly Pouch	Interior	Patient Information	100/Case
AS341	Bulk	Sterile	Tamper Evident	100/Case
AS343	Bulk	Non-Sterile	None	500/Case
AS345	Form-Fill & Seal	Fully Sterile	None	100/Case



AS340

# Index

REF	HCPCS CODE	PAGE NUMBER
ABL400	A4351	19
ABL402	A4351	19
ABL404	A4351	19
AS010	AS4322	23
AS011	AS4322	23
AS011P	AS4322	23
AS015	AS4322	23
AS016	AS4322	23
AS115	AS4322	23
AS116	AS4322	23
AS121	AS4320	22
AS126	AS4320	22
AS127	AS4320	22
AS130	AS4320	22
AS136	AS4320	22
AS302	A4358	7
AS306N	A4358	20
AS307N	A4358	20
AS308N	A4358	20
AS309N	A4358	20
AS311	AS4331	20
AS312	A4358	7
AS322	A4358	7
AS32200	A4358	7
AS326	A4358	7
AS32600	A4358	7
AS327100	A4358	10
AS328	A4358	8
AS332	A4358	7
AS340	P9615	24
AS341	P9615	24
AS343	P9615	24
AS345	P9615	24
AS40008	A4338	12
AS40010	A4338	12
AS409	AS4328	21
AS41012	A4338	12
AS41012HC	A4338	11
AS41012S	A4344	13
AS41012SC	A4344	13
AS41014	A4338	12

REF	HCPCS CODE	PAGE NUMBER
AS41014HC	A4338	11
AS41014S	A4344	13
AS41014SC	A4344	13
AS41014H	A4338	11
AS41016	A4338	12
AS41016H	A4338	11
AS41016HC	A4338	11
AS41016S	A4344	13
AS41016SC	A4344	13
AS41018	A4338	12
AS41018H	A4338	11
AS41018HC	A4338	11
AS41018S	A4344	13
AS41018SC	A4344	13
AS41020	A4338	12
AS41020HC	A4338	11
AS41020S	A4344	13
AS41020SC	A4344	13
AS41022	A4338	12
AS41022HC	A4338	11
AS41022S	A4344	13
AS41022SC	A4344	13
AS41024	A4338	12
AS41024HC	A4338	11
AS41024S	A4344	13
AS41024SC	A4344	13
AS41026	A4338	12
AS41026S	A4344	13
AS41028	A4338	12
AS41030	A4338	12
AS42012	A4338	12
AS42014	A4338	12
AS42014H	A4338	11
AS42016	A4338	12
AS42016H	A4338	11
AS42016HC	A4338	11
AS42016S	A4344	13
AS42018	A4338	12
AS42018H	A4338	11
AS42018HC	A4338	11
AS42018S	A4344	13

REF	HCPCS CODE	PAGE NUMBER
AS42018SC	A4344	13
AS42020	A4338	12
AS42020HC	A4338	11
AS42020S	A4344	13
AS42020SC	A4344	13
AS42022	A4338	12
AS42022HC	A4338	11
AS42022SC	A4344	13
AS42024	A4338	12
AS42024HC	A4338	11
AS42024SC	A4344	13
AS42026	A4338	12
AS42028	A4338	12
AS42030	A4338	12
AS44008	A4351	18
AS44010	A4351	18
AS44012	A4351	18
AS44014	A4351	18
AS44016	A4351	18
AS44018	A4351	18
AS44020	A4351	18
AS44022	A4351	18
AS44024	A4351	18
AS44028	A4351	18
AS85008	A4353	15
AS85008K	A4353	16
AS85010	A4353	16
AS85010K	A4353	16
AS85012	A4353	16
AS85012K	A4353	16
AS85014	A4353	16
AS85014K	A4353	16
AS85016	A4353	16
AS85016K	A4353	16
AS860608	A4351	17
AS860610	A4351	17
AS860612	A4351	17
AS860614	A4351	17
AS860616	A4351	17
AS861006	A4351	17
AS861008	A4351	17

REF	HCPCS CODE	PAGE NUMBER
AS861010	A4351	17
AS861608	A4351	17
AS861610	A4351	17
AS861610C	A4351	17
AS861612	A4351	17
AS861612C	A4351	17
AS861614	A4351	17
AS861614C	A4351	17
AS861616	A4351	17
AS861616C	A4351	17
AS861618	A4351	17
AS861618C	A4351	17
AS870	A4353	14
AS87114	A4353	14
AS87118	A4353	14
AS880	A4310	5
AS880K	A4310	5
AS886	A4310	5
AS890	A4310	5
AS890K	A4310	5
AS89214	A4353	15
AS89214R	A4353	15
AS89214RK	A4353	15
AS89216RK	A4353	15
AS89300	A4354	2
AS89314	A4314	1
AS89314H	A4314	1
AS89314S	A4315	1
AS89316	A4314	1
AS89316H	A4314	1
AS89316S	A4315	1
AS89318	A4314	1
AS89318H	A4314	1
AS89318S	A4315	1
AS89330	A4354	2
AS896	A4310	5
AS89800	A4354	4
AS89814	A4314	3
AS89814H	A4314	3
AS89814S	A4315	3
AS89816	A4314	3



# Index

REF	HCPCS CODE	PAGE NUMBER
AS89816H	A4314	3
AS89816S	A4315	3
AS89818	A4314	3
AS89818H	A4314	3
AS89818S	A4315	3
AS89830	A4354	4
LBE306N	A4358	20
LBE307N	A4358	20
LBE308N	A4358	20
LBE309N	A4358	20

## Sources

1. Centers for Disease Control and Prevention, <https://www.cdc.gov/infectioncontrol/guidelines/cauti/background.html>. Published November 5, 2015. Accessed September 30, 2022.
2. Gould, MD, MSCR, C., Umscheid, MD, MSCE, C., Agarwal, MD, MPH, R., Kuntz, MSW, MSLIS, G., Pegues, MD, D., & Healthcare Infection Control Practices Advisory Committee (HICPAC). (2019, June 6). GUIDELINE FOR PREVENTION OF CATHETER ASSOCIATED URINARY TRACT INFECTIONS 2009. [Www.Cdc.Gov/Infectioncontrol/Pdf/Guidelines/Cauti-Guidelines-H.Pdf](https://www.cdc.gov/infectioncontrol/Pdf/Guidelines/Cauti-Guidelines-H.Pdf). Retrieved September 30, 2022, from <https://www.cdc.gov/infectioncontrol/guidelines/cauti/>
3. De Ridder, D. J., Everaert K, Fernandez LG, et al. Intermittent catheterisation with hydrophilic-coated catheters (SpeediCath) reduces the risk of clinical urinary tract infection in spinal cord injured patients: A prospective randomised parallel comparative trial. *EurUrol.* 2005;48(6):991-995.
4. Stensballe J, Looms D, Nielsen PN, Tvede M. Hydrophilic-coated catheters for intermittent catheterisation reduce urethral micro trauma: A prospective, randomised, participant-blinded, crossover study of three different types of catheters. *EurUrol.* 2005;48(6):978-983.



Advancing Healthcare Worldwide®

**Amsino Medical Group Co. Ltd.**

708 Corporate Center Drive  
Pomona, CA 91768  
USA

P: +1 800.632.6746

F: +1 909.626.3888

info@amsino.com

**[www.amsino.com](http://www.amsino.com)**