

Swabable, Needle-free Valve for Vial Access

Swabable Septum

Silicone

Valve Body

Copolyester Tritan™




Features

- Microbial barrier swabable silicone stem with a split septum, flat seal to allow enhanced disinfection
- Compatible with ISO 11040-4, ISO 80369-7, Luer lock, Luer slip and prefilled glass syringes
- Transparent housing to allow clear view into the fluid pathway
- Compatible with TPN solutions, DMSO, blood and oncology drugs
- ETO Sterilized
- Shelf Life: 5 Years
- In compliance with MDD 93/42/EEC, EC Certificate No: 1984-MDD-18-524, FDA cleared
- Latex-free, PVC-free, DEHP-free

Order Information

REF	DESCRIPTION	LENGTH	PRIMING VOLUME	OVERFILL VOLUME	QTY/UNIT
AY0100	Needle-free Valve for Vial Access	3.1 cm	0.08 mL	0.111 mL	50/Box



To view our complete product line,
visit us at www.amsino.com



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 USA
+1 800.632.6746 | info@amsino.com

Tritan is a trademark of Eastman Chemical Company.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO is a registered trademark of Amsino International, Inc. © Copyright 2025, Amsino International, Inc. All rights reserved.