


NEW!

FlexGRIP® Foley Securement Device



A Smarter Approach to Foley Securement Starts Here

The FlexGRIP® Foley Securement Device is designed to provide superior stabilization for indwelling Foley catheters while enhancing patient comfort and helping reduce catheter-related complications. Built with soft, low-profile materials and a high-performance hypoallergenic adhesive, FlexGRIP delivers strong, reliable fixation without compromising catheter flow rates.

Why FlexGRIP Foley?

- **Secure stabilization for 11–18 Fr Foley catheters**
- **Helps reduce catheter migration and accidental dislodgement**
- **Skin-safe hypoallergenic adhesive designed to help prevent MARS!**
- **Low-profile, soft materials improve patient comfort**
- **Easy-to-use, pre-mounted fixation leash**
- **Additional locking mechanism for added security**
- **Positionable anywhere along the catheter shaft — not limited to the bifurcation**
- **Designed to maintain flow rates, even with soft silicone Foley catheters**
- **Sterile, single-use design**

Next-Generation Securement Technology

FlexGRIP combines innovative stabilization technology with a patient-centered design to create a securement solution clinicians can trust. The large adhesive footprint and soft cushioning materials provide exceptional hold while remaining gentle on the skin.

REF	DESCRIPTION	QTY/UNIT
FOL1118	Foley Catheter Fixation System	25/Box



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 USA
+1 800.632.6746 | info@amsino.com

To view our complete product line,
visit us at [Amsino.com](https://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO and FlexGRIP are registered trademarks of Amsino International, Inc. © Copyright 2026, Amsino International, Inc. All rights reserved.