

ENFit® Enteral Connection System

Purpose-built for safer enteral nutrition delivery

The AMSure® ENFit enteral syringe is intended for enteral nutrition and enteral medication administration through compatible enteral access devices in patients requiring enteral feeding support due to chronic physical impairment, swallowing dysfunction, gastrointestinal disorders, or other conditions requiring enteral nutritional support.

The device is designed exclusively for enteral applications and is not intended for IV or other non-enteral routes of administration.

The ENFit Connector:

- ISO 80369-3 compliant enteral connector design
- Helps reduce the risk of tubing misconnections with IV or other delivery systems
- Secure twist connection helps minimize accidental disconnections
- Supports safer enteral care practices across acute and home care settings



 **GEDSA®**
stay connected®

GEDSA supports global enteral safety standards.

For more information, visit
www.stayconnected.org
or scan the QR code:



Amsino Enteral Feeding Products are Rx Only and intended for patients requiring long-term or chronic enteral nutrition support. They are commonly utilized in long-term care facilities, nursing homes, rehabilitation centers, and home care environments.

Order Information

Item Code	Image	Description	Features	Qty/Unit
ENF1004K		AMSure Enteral Feeding Administration Kit Pump Set with ENFit Connector, compatible with Sentinel® and Compat™/ EntraFlo™ Pumps	Spike-Only No Bag	30/Case
ENF1200K		AMSure Enteral Feeding Administration Kit Pump Set with ENFit Connector, compatible with Sentinel® and Compat™/ EntraFlo™ Pumps	1200 mL Feeding Bag	30/Case
ENF1203K		AMSure Enteral Feeding Administration Kit Gravity Set with ENFit Connector	1200 mL Feeding Bag	30/Case
ENF1203KD		AMSure Enteral Feeding Administration Kit Gravity Set with ENFit Connector, Standard Bore (DEHP-Free)	1200 mL Feeding Bag	30/Case
ENF1203KLD		AMSure Enteral Feeding Administration Kit Gravity Set with ENFit Connector, Large Bore (DEHP-Free)	1200 mL Feeding Bag	30/Case
ENF500K		AMSure Enteral Feeding Administration Kit Pump Set with ENFit Connector, compatible with Sentinel® and Compat™/ EntraFlo™ Pumps	500 mL Feeding Bag	30/Case
ENS015		AMSure® 60 mL Thumb Controlled Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS015NS		AMSure 60 mL Thumb Controlled Piston Syringe with ENFit Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENS115T		AMSure 60 mL Flat Top Piston Syringe with ENFit Tip, transition adapter included	Non-Sterile Resealable Pole Bag*	30/Case
ENS115		AMSure 60 mL Flat Top Piston Syringe with ENFit Tip	Sterile Poly Pouch	50/Case
ENS115NS		AMSure 60 mL Flat Top Piston Syringe with ENFit Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENK127		AMSure Enteral Feeding Irrigation Kit – Includes: 500 mL Graduated Container, 60 mL Thumb Controlled Piston Syringe with ENFit Tip, Patient Label	Non-Sterile Resealable Pole Bag*	30/Case
ENK136		AMSure Enteral Feeding Irrigation Tray – Includes: 1000cc Outer Tray, 500cc Graduated Container, 60cc Thumb Control Ring Piston Syringe with ENFit Tip, tip protector, alcohol prep pad, large moisture-proof underpad, latex-free, sterile	Sterile Poly Pouch	20/Case

*Resealable pole bag designed for safe and convenient storage.

All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 USA

+1 800.632.6746 | info@amsino.com

To view our complete product line, visit us at **Amsino.com**

ENFit® is a registered trademark of GEDSA. *EntraFlo™* is a registered trademark of Medline Industries, Inc. *Compat™* is a registered trademark of Société des Produits Nestlé, S.A Corporation.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol „®“ indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Amsino and AMSure are registered trademarks of Amsino International, Inc. © Copyright 2026, Amsino International, Inc. All rights reserved.