


AMSafe® IV Catheter

IV Catheter is inserted into a patient's vascular system to sample blood, monitor blood pressure or administer fluids. IV Catheter provides access to a vein. IV Catheter consists of several parts: catheter hub assembly, needle hub assembly, protective cap assembly, and sheath.



ITEM CODE	SIZE	QUANTITY/UNIT
C5429732	14GX5 1/2" Central I.V. Catheter	500
C5439732	16GX5 1/2" Central I.V. Catheter	500
G7139732	14GX2" I.V. Catheter	500
G7149732	16GX2" I.V. Catheter	500
G7159732	18GX2" I.V. Catheter	500
G7169732	18GX1 1/4" I.V. Catheter	500
G7179732	20GX1 1/4" I.V. Catheter	500
G7189732	20GX2" I.V. Catheter	500
G7199732	22GX1 1/4" I.V. Catheter	500
G7209732	24GX3/4" I.V. Catheter	500
G9449732	26GX3/4" I.V. Catheter	500
P292	19GA BUTTERFLY WINGED NEEDLE	500
P293	21GA BUTTERFLY WINGED NEEDLE	500
P294	23GA BUTTERFLY WINGED NEEDLE	500
P295	25GA BUTTERFLY WINGED NEEDLE	500
P296	27GA BUTTERFLY WINGED NEEDLE	500
P299	21GA BUTTERFLY WINGED NEEDLE	500
P300	23GA BUTTERFLY WINGED NEEDLE	500
P301	25GA BUTTERFLY WINGED NEEDLE	500