

IMPROVED!

RECEPTAL®


3000 mL Canister and Liner

The RECEPTAL® 3000 mL canister and liners are engineered to manage high volume waste fluid collection up to 36 liters of suction fluid capacity without disconnection.

- **3000 mL Volume Indicator:** Color coded printed volume indicator matches color coding on RECEPTAL canisters for enhanced set-up safety and efficiency
- **Easily Identifiable Marked Ports:** Helps to minimize the potential of tubing mis-connections for improved patient care
- **Dual Port (Ortho/Standard):** Allows for simultaneous dual suction for improved clinical efficacy
- **Large Side Handle:** Designed for easy and safe handling of heavy fluids
- **Pour Spout:** Included on all 3L liners and can be used with specimen trap for convenient sample collection



Item Code	Description	Qty/Unit
43503-01	3000 mL RECEPTAL® Canister	8/Case
0L655-01	3000 mL Liner with VAC-GARD®, Dual Ports (Ortho/Standard), Pour Spout	50/Case
0L656-01	3000 mL Liner with Dual Ports (Ortho/Standard), Pour Spout	50/Case



Safety.
Cost Control.
Clinical Efficacy.

RECEPTAL®

1000 mL Canister and Liner



RECEPTAL® Waste Fluid Management Systems help hospitals manage suction waste by providing simple, safe and complete solutions. The 1000 mL canister and liner have been developed for low-volume waste fluid collection. The small, compact size helps maximize space efficiency and provides a cost-effective solution for portable suction systems.

- **VAC-GARD® Shut-Off Valve:** 99% bacterial eliminating filter prevents damage and contamination of vacuum system, reducing risk of exposure and cost
- **100% Leak-Tested Liner:**¹ Minimizes risk of exposure and contamination
- **Anti-Foaming Emulsion:** Prevents froth formation and aerosolization, reducing the risk of exposure and contamination while maximizing fluid capacity


Item Code	Description	Qty/Unit
43449-05	1000 mL RECEPTAL® Canister	10/Case
43010-05	1000 mL RECEPTAL® Liner with VAC-GARD® and Pour Spout	10/Case



Unique color coding on canister and liner lid designed for safety and improving outcomes.



White canister rim seal maximizes vacuum power to provide the best clinical performance.



Yellow vacuum connector for intuitive, simplified system set-up.

To view our complete product line, visit us at www.amsino.com

AMSINO

Advancing Healthcare Worldwide®

Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768

800.632.6746 • fax: 909.626.3888 • info@amsino.com

1. Data on file. Amsino International, Inc.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Receptal and Vac-Gard are registered trademarks of Amsino International, Inc. © Copyright 2018, Amsino International, Inc. All rights reserved.