

	Manufacturer	 Store between -10°C and +40°C	Temperature limit
	Authorized representative in the European Community		Consult instructions for use
	Manufacture Date		Caution
	Use-by date		Contains or presence of natural rubber latex
	Batch code		Contains or presence of phthalate
	Catalogue number		Non-pyrogenic
	Sterilized using ethylene oxide		Drops per millilitre
	Sterilized using irradiation		Fragile, handle with care
	Do not re-use		Keep away from sunlight
	Do not re-sterilize		Keep dry
	Not made with DEHP		Keep up
	Not made with natural rubber latex		Federal (USA) law restricts the use of this device to sale by or on the order of a physician.
	Non-sterile		STAR TUBE
	Do not use if package is damaged		Number of Units
	Sterile fluid path		Gravity type
	NEEDLE FREE		Not made with PVC
	NOT FOR IV USE		Medical Device
	EU importer		Serial number
	Refer to instruction manual/booklet Signifies that the instruction manual/booklet must be read		Class II equipment (degree of protection against electrical shock), double insulated
	Type BF Part		No protection against ingress of water
	Restriction of Hazardous Substances		WEEE Signifies waste from electrical and electronic equipment
	Liquid filter with pore size		Marks for UL's Classification Service
	Humidity limitations		Atmospheric pressure limitation